EN
3

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Bruksanvisning

T _ K N T C _ R E V 3

KAYE NEPHROSTOMY TAMPONADE BALLOON CATHETER AND STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions prior to using this device.

DEVICE DESCRIPTION

- 17.0 Fr radiopaque catheter 25 cm long with balloon 15 cm long by 12 mm (36.0 Fr) diameter
- Non-radiopaque flexible stylet
- 5.0 Fr radiopaque polyurethane stent 75 cm long
- 14.0 Fr polyvinylchloride Y-connecting tube
- 10 mL syringe

NOTE: Product will accept up to an .038 inch diameter wire guide.

NOTE: Set components may vary.

INTENDED USE

- The Kaye Nephrostomy Tamponade Balloon Catheter and Stent Set is used for nephrostomy drainage, ureteral drainage, and dynamic tamponade to prevent hemorrhage in nephrolithotomy and percutaneous stone removal procedures.
- The Kaye Nephrostomy Tamponade Balloon Catheter Set is used for nephrostomy drainage, and dynamic tamponade to prevent hemorrhage in nephrolithotomy and percutaneous stone removal procedures.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk for percutaneous catheterization and/or dynamic tamponade of a nephrostomy tract.

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.
- Not recommended for dilation of a nephrostomy tract.
- Balloon Catheter may remain in place for a total of 48 hours.
- Stent may remain in place for a total of 96 hours.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized, and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

NOTE: Remove all instrumentation related to previous procedure, with the exception of the wire guide. If the wire guide was removed during the procedure, replace it under endoscopic vision.

1. Withdraw protective sheath from balloon of catheter. Discard sheath.
2. Using syringe, inject 3-4 mL of dilute contrast medium into the balloon to assist with accurate placement.


NOTE: Always inflate balloon with sterile liquid. Never inflate with air, carbon dioxide, or any other gas.

NOTE: Refer to product label for maximum inflation pressure. Do not exceed the maximum inflation pressure for this balloon device.

3. Pass the catheter over previously placed wire guide (not in set) into the renal pelvis.
4. Under fluoroscopic control, position the catheter so that 2-3 cm of the balloon is within the collecting system.

NOTE: It is inadvisable to have the balloon catheter too far into the pelvis or down the ureter since the balloon may occlude infundibulae, or the catheter tip may irritate the pelvis (**Figure A**).

Figure A

5. Once the balloon catheter is in the desired location, the inner stylet can be removed.
6. To effect dynamic tamponade, inflate the balloon with dilute contrast medium to a maximum pressure of 40 psi (275 kPa). Only moderate pressure is required for tamponade.
7. Before placing the stent, remove the Tuohy-Borst adapter and the stent-to-catheter adapter or side-arm adapter. Advance the stent over the previously placed wire guide, through the balloon catheter, and into the desired location within the ureter.
8. Replace the stent-to-catheter adapter or side-arm adapter, and the Tuohy-Borst adapter, and adjust the length of the stent in the ureter as desired. Secure the locking cap on the stent-to-catheter adapter or side-arm adapter finger tight around the stent. **NOTE:** Since the distal 33 cm of the Kaye stent is sideported, the stent may be advanced fully into the bladder, thereby enhancing drainage from the renal pelvis into the bladder.
9. Check positioning of balloon catheter and stent with fluoroscopy. If satisfied, remove the wire guide.
10. Reconfirm position with an intraoperative nephrostogram.
11. Secure balloon catheter with tape and dressings as needed.

NOTE: Do not place tape directly on balloon.

NOTE: Leave the inflation valve in an accessible position.

12. Lock one leg of the Y-connecting tube onto the side-arm fitting for nephrostomy drainage. Lock the other leg onto the Tuohy-Borst adapter secured on the end of the stent for stent lumen drainage (**Figure B**).


Figure B

SUGGESTED REMOVAL INSTRUCTIONS

- A)
 1. Completely deflate the balloon; remove entire balloon catheter and stent as one. Evaluate the wound site for bleeding.

B)

1. Provided ureteral drainage is still desirable and tract tamponade has been effected, the Kaye balloon catheter may be removed (see **Removal Instruction A**) and the Kaye stent left to drain internally. Secure the stent to the skin and clamp the proximal end.
2. Using a nephrostogram, check the positioning of the stent to assure accurate position and the distal end of the stent is in the bladder. **NOTE:** The Kaye stent may also drain externally, thereby serving as a small nephrostomy tube. Drainage can be accomplished by reattaching the Tuohy-Borst adapter to the proximal end of the stent and connecting it to a connecting tube. Should bleeding occur while the stent is in place, reinsert a wire guide through the stent and reintroduce the balloon catheter over the stent.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.