EN
3

Kwart Retro-Inject™ Stent Set

Instructions for Use

CS
4

Souprava stentu Kwart Retro-Inject™

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DA
6

Kwart Retro-Inject™ stentsæt

Brugsanvisning

DE
8

Kwart Retro-Inject™ Stentset

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EL
10

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
KWART RETRO-INJECT™ STENT SET

SO[®]-FLEX®, LSe SO[®]-FLEX®, AND SO[®]-FLEX AQ®

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque stent with tether
- Radiopaque inserter
- Appropriate diameter TFE-coated stainless steel wire guide, 145 cm long with 3 cm flexible tip
- Radiopaque release sleeve
- Female Luer lock to Tuohy-Borst adapter

NOTE: Set components may vary.

INTENDED USE

Used for retrograde injection during E.S.W.L. and leaving an indwelling ureteral stent post-E.S.W.L.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Complications of ureteral stent placement are documented in literature. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- The Sof-Flex, LSe Sof-Flex, and Sof-Flex AQ stents must not remain indwelling more than six months.
- These stents are not intended as permanent indwelling devices.
- Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Quality Assurance department at the address/phone listed at www.cookmedical.com.

INSTRUCTIONS FOR USE

NOTE: Prior to use, immerse AQ® stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

NOTE: Do not remove tether prior to stent placement.

1. Under X-ray or fluoroscopic control, cystoscopically pass the wire guide through the ureter to the renal pelvis.
2. Pass the stent over the wire guide through the cystoscope.
3. Pass the inserter over the wire guide and push it into the end of the stent. **NOTE:** The inserter is pushed into the stent to the 5 mm index mark of the inserter.
4. Under direct vision, advance the stent into the ureter with the inserter.
5. When the renal pelvis has been entered (this can be confirmed fluoroscopically), gently remove the wire guide while holding the inserter. The stent coil will form spontaneously in the renal pelvis. **NOTE:** The graduations on the stent, when viewed cystoscopically, aid in proper stent placement.
6. Attach the female Luer lock to the Tuohy-Borst adapter on the inserter. Secure the tether by taping it to the inner thigh to avoid inadvertent displacement.
7. Dilute contrast medium may be injected intermittently throughout the E.S.W.L. procedure to disclose stones or stone fragments for targeting.
NOTE: The inserter can be removed post-E.S.W.L., or left in place if desired and utilized for external drainage or injection.

Removal of the Inserter Post-E.S.W.L.

1. Remove the female Luer lock to Tuohy-Borst adapter from the inserter and pass the release sleeve over the inserter until resistance to advancement is felt. **NOTE:** A wire guide may be readvanced into the inserter and through the stent before the release sleeve is positioned. The wire guide would serve as a safety wire should the stent become inadvertently dislodged during the inserter removal.
2. Hold the release sleeve in position to prevent the stent from being dislodged. If a wire guide has been placed, it should be removed prior to removing the inserter. Gently remove the inserter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.