EN
3 **Kwart Retro-Inject™ Ureteral Stent Set**

Instructions for Use

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Kullanma Talimatı

KWART RETRO-INJECT™ URETERAL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque stent with tether
- Radiopaque inserter
- Appropriate diameter TFE-coated stainless steel wire guide, 145 cm long with 3 cm flexible tip
- Radiopaque release sleeve
- Adapter: Female Luer lock to Tuohy-Borst adapter

NOTE: Set components may vary.

INTENDED USE

Used for retrograde injection during Extracorporeal Shock Wave Lithotripsy (E.S.W.L.) and leaving an indwelling ureteral stent post-E.S.W.L.

CONTRAINDICATIONS

None known

WARNINGS

Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

PRECAUTIONS

- Complications of ureteral stent placement are documented in literature. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- Stents must not remain indwelling more than six months.
- These stents are not intended as permanent indwelling devices.
- The tether should be removed if the stent is to remain indwelling longer than 14 days.
- Do not force components during removal or replacement. If resistance is encountered, stop. Do not continue without first determining the cause of resistance and taking remedial action.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Quality Assurance department at the address/phone number listed at www.cookmedical.com.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Kwart Retro-Inject Stent is **MR Conditional**. A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or 1.5 Tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Kwart Retro-Inject Stent is not expected to induce RF heating greater than that of biological tissue.

The image artifact extends approximately 1.9 mm from the Kwart Retro-Inject Stent as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

POTENTIAL ADVERSE EVENTS

Complications of ureteral stent placement are documented. These complications include, but are not limited to:

- Extravasations
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract
- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function

INSTRUCTIONS FOR USE

NOTE: Prior to use, immerse AQ® stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Under x-ray or fluoroscopic control, cystoscopically pass the wire guide through the ureter to the renal pelvis.
2. Pass the stent over the wire guide through the cystoscope.
3. Pass the inserter over the wire guide and push it into the end of the stent. **NOTE:** The inserter is pushed into the stent to the 5 mm index mark of the inserter.
4. Under direct vision, advance the stent into the ureter with the inserter.
5. When the renal pelvis has been entered (this can be confirmed fluoroscopically), gently remove the wire guide while holding the inserter. The stent coil will form spontaneously in the renal pelvis.
6. Attach the adapter on the inserter and secure the tether by taping it to the inner thigh to avoid inadvertent displacement.
7. Dilute contrast medium may be injected intermittently throughout the E.S.W.L. procedure to disclose stones or stone fragments for targeting.

Removal of the Inserter Post-E.S.W.L.

1. Remove the adapter from the inserter and pass the release sleeve over the inserter until resistance to advancement is felt. **NOTE:** A wire guide may be readvanced into the inserter and through the stent before the release sleeve is positioned. The wire guide would serve as a safety wire should the stent become inadvertently dislodged during the inserter removal.
2. Hold the release sleeve in position to prevent the stent from being dislodged. If a wire guide has been placed, it should be removed prior to removing the inserter. Gently remove the inserter.

Stent Removal

The stent may be removed by gentle withdrawal traction using endoscopic forceps or tether.

Tether Options:

A. Removal after stent positioning:

After stent placement, hold the knot of the tether and cut one strand. While holding the knot, gently pull on the tether to remove, maintaining stent position with positioner and wire guide.

B. Stent removal:

The tether may be used for subsequent stent removal. To avoid possible tether encrustation, the tether should be removed if the stent is to remain indwelling longer than 14 days.

C. Removed prior to stent placement:

Remove tether by holding the knot, cutting one strand, and while holding the knot, gently pull on the tether.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.