- EN**
3 **Topel Endoscopic Cyst Aspirator Set**
Instructions For Use
- CS**
5 **Souprava Topel pro endoskopickou aspiraci cyst**
Návod k použití
- DA**
7 **Topel endoskopisk cysteaspiratorsæt**
Brugsanvisning
- DE**
10 **Topel endoskopisches Zystenaspirator-Set**
Gebrauchsanweisung
- EL**
13 **Σετ συσκευής ενδοσκοπικής αναρρόφησης κύστης Topel**
Οδηγίες χρήσης
- ES**
16 **Equipo aspirador de quistes endoscópico Topel**
Instrucciones de uso
- FR**
18 **Set d'aspirateur de kyste endoscopique de Topel**
Mode d'emploi
- HU**
21 **Topel endoszkópos cisztaleszívató készlet**
Használati utasítás
- IT**
24 **Set con aspiratore endoscopico di cisti Topel**
Istruzioni per l'uso
- NL**
26 **Topel endoscopische cyste-aspiratorset**
Gebruiksaanwijzing
- NO**
29 **Topel endoskopisk cyste-aspiratorsett**
Bruksanvisning
- PL**
32 **Zestaw ssaka Topel do endoskopowego odsysania torbieli**
Instrukcja użycia
- PT**
34 **Conjunto de aspirador de quistos endoscópico Topel**
Instruções de utilização
- SV**
37 **Topel endoskopiskt cystaspiratorsett**
Bruksanvisning
- ZH**
40 **Topel内窥镜囊肿抽吸套装**
使用说明

TOPEL ENDOSCOPIC CYST ASPIRATOR SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

Used to stabilize a cyst or other anatomic structure under direct endoscopic vision and allow aspiration of its fluid content without spillage.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Should product label contain the PHT symbol this warning statement applies; the potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

PRECAUTIONS

- This product must be cleaned immediately after use. Blood, mucus and other debris, if allowed to dry, can render the device inoperable.

INSTRUCTIONS FOR USING TOPEL ENDOSCOPIC CYST ASPIRATOR SET

1. Connect the cyst aspirator to full wall suction and close the vacuum control valve to prevent aspiration of distention media (CO₂) during introduction. The cyst aspirator is passed through a 10 mm trocar channel into the abdominal cavity. In most instances, the aspirator enters from the side opposite the cyst.
2. An endoloop suture is inserted from the contralateral side and placed over the surface of the cyst.
3. The cyst aspirator is placed through the endoloop and gently held against the surface of the cystic structure.
4. After proper positioning of the cyst aspirator, wall suction is applied by opening the vacuum control valve; the surface of the cyst becomes firmly attached to the aspirating channel.
5. A 14 or 17 gauge aspirating needle is passed through the center channel and advanced to the first depth mark on the proximal end of the needle. At this point the needle tip rests at the distal end of the channel.
6. A large syringe is then attached to the stopcock on the aspirating needle.
7. The aspirating needle is then advanced into the cystic structure and fluid is aspirated. The depth of penetration can be determined by the centimeter marks located at the proximal end of the needle.
8. After the cystic structure is totally evacuated and collapsed, the aspirating needle is withdrawn. The endoloop is drawn tightly below the cyst aspirator to ligate the puncture site. The suture end is then cut and withdrawn.
9. Using a grasping instrument, the cystic tissue inside the channel is gently detached from the cyst aspirator. At times, wall suction can be slightly decreased to facilitate this task.
10. The cyst aspirator is then immediately withdrawn from the abdominal cavity to prevent collapse of the pneumoperitoneum.
 - It is important not to break suction until the aspirator has been completely withdrawn. This will prevent fluid from dripping back down the aspirating channel into the abdominal cavity.
11. The cystic structure is now ready for endoloop ligation and removal from the abdomen.

CLEANING AND STERILIZATION INSTRUCTIONS

- These instructions apply only to the reusable stainless steel cyst aspirator. The needles and the disposable cyst aspirator are intended for one-time use.
1. Disassemble the cyst aspirator.
 - **IMPORTANT:** Do not attempt to remove the screw from the vacuum control valve.
 2. Remove all residual organic matter, blood and irrigation solution to allow the sterilization or disinfection medium to contact the surface of the device.
 - Ethylene oxide (ETO) gas and disinfection solutions do not effectively penetrate residual organic matter.
 3. Clean all exterior surfaces with a solution of mild detergent and water. Use a soft brush to remove all residual organic matter, blood and irrigation solution from the device.
 4. Ultrasonic cleaning of the components may be performed. Ultrasonic cleaning should only be performed with non-abrasive and non-corrosive solutions.
 5. Rinse thoroughly with clean water.
 6. Dry all exterior surfaces with a soft cloth or an air pistol.
 7. Dry all openings with an air pistol.
 - Prior to sterilization and reuse, this product should be inspected to verify that it functions properly. All edges and surfaces should be smooth and free of abrasions.

ETHYLENE OXIDE GAS STERILIZATION

Routine ethylene oxide (ETO) gas sterilization is the sterilization method of choice.

ETO gas sterilization and aeration should be performed in accordance with AAMI Standards and Recommended Practices with a validated sterilization cycle.

AUTOClave STERILIZATION

This device can be sterilized by steam autoclave. Autoclave sterilization should be performed in accordance with AAMI Good Hospital Practices with a validated sterilization cycle. Autoclave temperature should not exceed 275°F.

DISINFECTION

Chemical sterilization should be performed in accordance with AAMI Good Hospital Practices.

- Only EPA registered hard surface disinfectants capable of high-level disinfection should be used.
1. Carefully follow the disinfectant manufacturer instructions.
 2. Devices to be disinfected should be placed in a plastic basin to avoid scratches and electrolytic corrosion.
 3. Rinse disinfected devices thoroughly to remove all traces of residue.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

H. C. Topel, M.D., Department of Obstetrics and Gynecology, Lutheran General Hospital, Park Ridge, Illinois.
U.S. Patent Number: 4,958,621.

AAMI/ISO 11135-1994-Medical Devices-Validation and Routine Control of Ethylene Oxide Sterilization.
ANSI/AAMI ST79:2006-Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

ANSI/AAMI ST41:1999-Ethylene Oxide Sterilization in Health Care Facilities: Safety and Effectiveness.

ANSI/AAMI ST58:2005-Chemical Sterilization and High-Level Disinfection in Health Care Facilities.

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ANSI/AAMI ST79:2006-Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities (卫生保健机构蒸汽灭菌及灭菌保证综合指南)。

ANSI/AAMI ST41:1999-Ethylene Oxide Sterilization in Health Care Facilities: Safety and Effectiveness (卫生保健机构环氧乙烷灭菌: 安全及有效性)。

ANSI/AAMI ST58:2005-Chemical Sterilization and High-Level Disinfection in Health Care Facilities (卫生保健机构化学灭菌及强力消毒说明)。

MANUFACTURER

COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

www.cookmedical.com
© COOK 2019

EC

REP

EC REPRESENTATIVE

Cook Medical Europe Ltd
O'Halloran Road
National Technology Park
Limerick, Ireland

2019-09
T_KTEA_REV1