- EN**  
**3** **LithAssist®**  
**Suction Control for Laser Lithotripsy**  
Instructions for Use
- CS**  
**4** **LithAssist®**  
**Kontrola odsávání pro laserovou litotripsii**  
Návod k použití
- DA**  
**5** **LithAssist®**  
**Sugekontrol til laserlitotripsi**  
Brugsanvisning
- DE**  
**6** **LithAssist®**  
**Saugkraftregler für die Laser-Lithotripsie**  
Gebrauchsanweisung
- EL**  
**7** **Συσκευή ελέγχου αναρρόφησης για λιθοτριψία με λέιζερ**  
**LithAssist®**  
Οδηγίες χρήσης
- ES**  
**8** **LithAssist®**  
**Control de aspiración para litotricia por láser**  
Instrucciones de uso
- FR**  
**9** **LithAssist®**  
**Dispositif de contrôle de l'aspiration pour lithotripsie laser**  
Mode d'emploi
- HU**  
**10** **LithAssist®**  
**Szívásszabályzó lézeres kőzúzáshoz**  
Használati utasítás
- IT**  
**12** **LithAssist®**  
**Dispositivo ausiliario ad aspirazione controllata per litotripsia laser**  
Istruzioni per l'uso
- NL**  
**13** **LithAssist®**  
**Zuigkrachtregeling voor laserlithotripsie**  
Gebruiksaanwijzing
- NO**  
**14** **LithAssist®**  
**Sugekontroll for laserlitotripsi**  
Bruksanvisning
- PL**  
**15** **LithAssist®**  
**Regulator ssania do litotrypsji laserowej**  
Instrukcja użycia
- PT**  
**16** **LithAssist®**  
**Controlo de aspiração para litotripsia a laser**  
Instruções de utilização
- SV**  
**17** **LithAssist®**  
**Sugkontroll för laserlitotripsi**  
Bruksanvisning

## LITHASSIST®

### SUCTION CONTROL FOR LASER LITHOTRIPSY

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

#### DEVICE DESCRIPTION

The LithAssist is a rigid device with a handle for control, used in conjunction with a holmium laser for lithotripsy. The device has a dual lumen cannula. The small lumen provides access for the introduction of a laser fiber for use during lithotripsy of kidney or bladder stones. The large lumen provides suction assistance to extract fluid, soft tissue and stone fragments when used with an independent vacuum source. The device features a trigger mechanism that enables the user to manually control the amount of suction.

#### INTENDED USE

The LithAssist is intended to provide access for laser fiber introduction and suction for fluid, soft tissue, and stone fragment removal during laser lithotripsy.

#### CONTRAINDICATIONS

None known

#### WARNINGS

None known

#### PRECAUTIONS

- This device is conductive. Avoid contact with any electrified instrument.
- Take normal precautions when using a laser (i.e., safety glasses).
- This product is intended for use by physicians trained and experienced in urologic lithotripsy techniques. Standard techniques for use of an endoscope and laser fiber should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

#### INSTRUCTIONS FOR USE

1. Attach suction tubing to the universal adapter at the bottom of the LithAssist. Adjust suction to preference.
2. Using previously established kidney/bladder access, maintain direct vision with a rigid endoscope.
3. Insert the LithAssist into the working channel of the endoscope until the distal end is visible.
4. Maintain direct visualization of the stone, and insert a laser fiber through the Tuohy Borst adapter at the rear of the LithAssist.
5. Once the laser fiber extends beyond the distal end of the LithAssist to the desired distance, tighten the Tuohy Borst adapter around the laser fiber to secure it in place.
6. Control suction via the handle trigger while lasing the stone. **NOTE:** A fully open trigger provides the highest suction; a fully closed trigger provides no suction. The trigger may be locked into the fully closed position via the sliding switch immediately above the handle.  
**NOTE:** The device may become clogged with stone material when suctioning. If this occurs, remove the device from the patient, fully depress/lock the handle to shut off the suction, and flush the obstructed lumen with a saline-filled syringe via the Luer lock on the rear of the device.
7. Repeat steps 3 through 6 as necessary.

#### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a

dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd