Northstar® Lumax® Flex and White Lumax® Guiding Catheters

Instructions for Use

EN
2

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T _ L F G _ R E V . 2

NORTHSTAR® LUMAX® FLEX AND WHITE LUMAX® GUIDING CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Northstar® Lumax® Flex Guiding Catheter

The Northstar Lumax Flex Guiding Catheter includes a braided nylon guiding catheter and a coaxial inner catheter tapered to either a .035 or .038 inch endhole, with a soft radiopaque tip. In products whose order numbers begin with the prefix LMGN, the distal 50 cm of the guiding catheter's outer surface is hydrophilically coated.

White Lumax® Guiding Catheter

The White Lumax Guiding Catheter includes a Lumax Flex guiding catheter with a radiopaque tip and a coaxial Slip-Cath® inner catheter. Both guiding catheter and coaxial catheter are constructed of radiopaque braided nylon. The distal 50 cm of the coaxial Slip-Cath catheter is hydrophilically coated.

INTENDED USE

Northstar Lumax Flex and White Lumax Guiding Catheters are intended for the delivery of angioplasty balloons and other types of interventional devices.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- Do not completely block flow within vessel or inadvertently occlude flow to small vessels.
- Use caution when injecting contrast. Due to the catheter's large internal diameter, excessive flow can be generated by using little force on the syringe.
- Northstar Lumax Flex and White Lumax Guiding Catheters have been designed to provide maximum backup support for complex interventional procedures. Care should be taken during manipulation of the catheter to prevent injury or dissection to the vessel wall.
- **White Lumax Guiding Catheter components must be introduced through a percutaneous introducer sheath.**

PRECAUTIONS

- Manipulate catheter using fluoroscopic control.
- Catheter tip is manufactured from a heat-sensitive material. Do not attempt to reshape catheter curves.
- Due to extreme thinwall construction of catheter, use through synthetic vascular grafts or known highly calcified lesions should be avoided. Damage to catheter may result.
- Extreme care must be exercised during manipulation and withdrawal to prevent pulling catheter apart.
- If an interventional cardiovascular device should become lodged during use, the guiding catheter and interventional device should be withdrawn as a unit through the introducer sheath.
- Activate hydrophilic coating on distal segment of outer catheter (Northstar Lumax Flex Guiding Catheter) or inner catheter (White Lumax Guiding Catheter) by wetting the surface of the catheter with sterile water or saline. For best results, maintain wetted condition during use.

PRODUCT RECOMMENDATIONS

French Size	Guiding Catheter Internal Diameter	Recommended Wire Guide Diameter for Coaxial Introduction
6.0	.060 inch	.035 inch
7.0	.073 inch	.038 inch
8.0	.086 inch	.038 inch
9.0	.099 inch	.038 inch

INSTRUCTIONS FOR USE

1. Insert the inner catheter into the guiding catheter and lock into position.
2. After accessing the vascular system with a percutaneous introducer sheath, introduce the coaxial guiding catheter assembly over an appropriately sized standard wire guide.
3. After introducing the coaxial guiding catheter assembly into the primary vascular system, the inner catheter may be removed.
4. The guiding catheter may then be manipulated in a standard fashion. **NOTE:** If resistance is encountered during manipulation, stop and determine the cause before proceeding any further.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATENTS

US Patent Numbers: 5,221,270; 5,300,048; 4,979,959; 4,973,493; 5,002,582; 5,217,492; 5,258,041; 5,263,992; Australian Patent Numbers: 211631; 652254; European Patent Number: EP0520692. Other patents pending.

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Bloomington, IN 47404 U.S.A.