**COOK**

**MEDICAL**

**EN**  
**3**

## **Guiding Catheter Lumax®**

Instructions for Use

**FR-CA**  
**5**

## **Cathéter guide Lumax®**

Mode d'emploi

T \_ L F G \_ R E V 5

## LUMAX® GUIDING CATHETER

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Lumax Guiding Catheter includes either a 7 French guiding catheter with a radiopaque tip or an 8 French guiding catheter with a radiopaque tip and a coaxial Slip-Cath® inner catheter. Both guiding catheter and coaxial catheter are constructed of radiopaque braided nylon. The distal 40 cm of the coaxial Slip-Cath catheter is hydrophilically coated.

### INTENDED USE

Lumax Guiding Catheters are intended for the intravascular introduction of interventional devices into the peripheral, including pulmonary, vasculature.

### CONTRADICTIONS

None known

### WARNINGS

- Testing conducted does not support the use of this device in patients who cannot be systemically anticoagulated.
- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Do not completely block flow within vessel or inadvertently occlude flow to small vessels.
- Use caution when injecting contrast. Due to the catheter's large internal diameter, excessive flow can be generated by using little force on the syringe.
- Lumax Guiding Catheters have been designed to provide maximum backup support for complex interventional procedures. Care should be taken during manipulation of the catheter to prevent injury or dissection to the vessel wall.
- Lumax Guiding Catheter components must be introduced through a percutaneous introducer sheath.
- Administer appropriate anticoagulation therapy to reduce potential thrombosis. If the patient is not appropriately anti-coagulated, thrombus formation may occur. The safety of this product when used without anticoagulant has not been evaluated.

### PRECAUTIONS

- Recommended use of an .038-inch wire guide or smaller for 8 French coaxial introduction.
- Manipulate catheter using fluoroscopic control.
- Catheter tip is manufactured from a heat-sensitive material. Do not attempt to reshape catheter curves.
- Due to extreme thinwall construction of catheter, use through synthetic vascular grafts or known highly calcified lesions should be avoided. Damage to catheter may result.
- Extreme care must be exercised during manipulation and withdrawal to prevent pulling catheter apart.
- If an interventional cardiovascular device should become lodged during use, the guiding catheter and interventional device should be withdrawn as a unit through the introducer sheath.
- Activate hydrophilic coating on distal segment of inner catheter by wetting the surface of the catheter with sterile water or saline. For best results, maintain wetted condition during use.

### INSTRUCTIONS FOR USE

1. If using the 8 French guiding catheter, begin by inserting the inner catheter into the guiding catheter and lock into position.
2. After accessing the vascular system with a percutaneous introducer sheath, introduce either the 7 French guiding catheter or the 8 French coaxial guiding catheter assembly over an appropriately sized standard wire guide.
3. After introduction of the 8 French coaxial guiding catheter assembly into the primary vascular system, the inner catheter may be removed.
4. The guiding catheter may then be manipulated in a standard fashion.  
**NOTE:** If resistance is encountered during manipulation, stop and determine the cause before proceeding any further.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

**REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

**PATENTS**

US Patent Numbers: 5,221,270; 5,300,048; 4,979,959; 4,973,493; 5,002,582; 5,217,492; 5,258,041; 5,263,992; Australian Patent Numbers: 211631; 652254; European Patent Number: EP0520692. Other patents pending.

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Bloomington, IN 47404 U.S.A.