**COOK**

MEDICAL

**CE**  
0123**EN**  
3**Multipurpose Drainage Catheter  
With Locking Luer Fixation**

Instructions for Use

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Bruksanvisning

T \_ L M C L \_ R E V 4

## MULTIPURPOSE DRAINAGE CATHETER WITH LOCKING LUER FIXATION

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

Multipurpose Drainage Catheters are constructed from Ultrathane® or polyethylene and come in a range of French sizes, lengths and sideport quantities.

When ordered with the suffix –HC, the distal part of the catheter, sideports included, has an AQ® hydrophilic coating.

### INTENDED USE

Multipurpose drainage catheters are intended for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess), either by direct stick or Seldinger access technique.

### CONTRAINDICATIONS

- Bleeding diathesis and uncontrolled hypertension
- Anticoagulant use

### WARNINGS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- This product contains **NATURAL RUBBER LATEX**, which may cause allergic reactions.

### PRECAUTIONS

- These products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous drainage catheters should be employed.
- Manipulation of products requires ultrasound, fluoroscopy, or other imaging guidance.
- When inserting a stiffening cannula into a catheter with retention suture, hold suture during cannula insertion to avoid bunching or tangling of suture.
- A TFE-coated wire guide must be used with Ultrathane® catheters.
- Activate the hydrophilic coating, if present, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.
- Catheters should be irrigated on a routine basis to ensure function.
- Patients with indwelling drainage catheters should be evaluated routinely to ensure continuous function of the catheter.
- Traction on the locking suture, if present, should be sufficient to ensure adequate retention of the tip, but should not be overly tight. Verify catheter tip configuration by fluoroscopy.
- It is recommended to use a wire guide when removing a Locking Loop catheter.
- The Peel-Away® Pigtail Straightener, if present, is not to be used as a vascular introducer sheath.

### POTENTIAL ADVERSE EVENTS

- Arrhythmia (transient), including tachycardia
- Bleeding (including chronic hematuria, post-procedural, intraprocedural, puncture site, hemoperitoneum)
- Bile leak/biloma
- Catheter dislocation/dislodgement
- Catheter malposition/Transcolonic catheter placement
- Catheter occlusion/obstruction
- Damage to adjacent organs
- Death (due to uncontrolled hypotension)
- Hypoxia
- Inflammation/Infection (including cellulitis at puncture site, Fever/chills within 24 hours, rigors, sepsis during catheter exchange, septic shock leading to death, urinary tract infection, peritonitis)
- Mechanical catheter damage
- Pain (including intraprocedural pain, post-procedural pain)
- Perihepatic fluid collection
- Small bowel injury
- Urinary leak

### MRI SAFETY INFORMATION

The Multipurpose Drainage Catheter with Locking Luer Fixation is **MR Safe** according to ASTM F2503. The accessory components have not been tested for MRI compatibility.

### INSTRUCTIONS FOR USE

#### Catheter Placement

1. Under fluoroscopic control, perform standard techniques for placement of percutaneous drainage catheters, either by Seldinger access or trocar access.
2. Once catheter is in desired location, remove any wire guides, trocars, or stiffeners, allowing the catheter to form its configuration.
3. Pull back on the drawstring to form the catheter's loop configuration, then lock the loop in place by securing the secondary hub onto the catheter's female luer fitting.

**NOTE:** To unlock the catheter loop, release the drawstring by removing the secondary hub.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd