EN
3

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Bruksanvisning

T _ L S _ R E V 2

LAWSON RETROGRADE NEPHROSTOMY WIRE PUNCTURE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions prior to using this device.

DEVICE DESCRIPTION

- .038 inch (0.97 mm) diameter stainless steel wire guide 145 cm long with 3 cm flexible tip
- 7.0 French radiopaque polyethylene Torcon® catheter 82 cm long with .045 inch (1.14 mm) diameter tip deflecting wire guide and three-ring handle
- 3.0 French TFE sheath 85 cm long with .017 inch (0.43 mm) diameter stainless steel puncture wire 145 cm long with plastic pin vise actuator
- Disposable coaxial needle: inner needle with stylet 22 gage 22.5 cm long and outer needle 18 gage 14 cm long

INTENDED USE

Used to gain precise percutaneous access to the kidney by means of controlled fine wire puncture from within the collecting system for patients aged 12 years and older.

CONTRAINDICATIONS

- No suitable calyces found on pre-op workup.
- Other situations where the risk benefit ratio for retrograde nephrostomy is not favorable.

WARNINGS

None known

PRECAUTIONS

- Stone burden in the ureter and/or the calyces that might impede device tracking and/or deployment.
- Urinary stones completely filling or tightly adhering to the calyces walls.

POTENTIAL ADVERSE EVENTS:

- Hemorrhage
- Retroperitoneal hematoma
- Unintended organ perforation
- Kidney injury and laceration
- Sepsis
- Pleural effusion
- Extravasation/perforation
- Tissue injury

INSTRUCTIONS FOR USE

1. Under fluoroscopic control, carefully pass the .038 inch (0.97 mm) diameter flexible tip wire guide cystoscopically and retrograde well into the renal pelvis.
2. Advance the 7.0 French catheter over the wire guide until the catheter tip resides in the renal pelvis. Remove the wire guide.
3. Advance the .045 inch (1.14 mm) diameter tip deflecting wire guide into the catheter until the assembly can be twist-locked together. (It will not pass completely through the catheter tip.)
4. Use combined twisting and deflecting of the catheter to maneuver the catheter tip well into the selected calyx. (A middle or lower pole posterior calyx is generally selected, depending on the procedure planned [Figures A and B]. An upper pole calyx [Figure C] is generally avoided due to twelfth-rib obstruction, undesirable intercostal puncture, or puncture of adjacent viscera such as liver, spleen, or lung.)


Figure A


Figure B

After IVP or Retrograde Disclosure, A-P Fluoroscopy will show:

- (1) Cup-like structures: Anterior Calyces**
- (2) End-on, round concentrations: Posterior Calyces**


Figure C
(Generally Avoided)

CAUTION: Viewed A/P — approximately 69% of Right kidneys are Brodel type (Figure A) and end-on calyces are Posterior. Approximately 75% of Left kidneys are Hodson type and Lateral calyces are Posterior calyces.

5. Remove the tip deflecting wire guide from the catheter, taking care to maintain the inner-calyx position of the catheter tip. Inject a bolus of contrast at this point, if desired.

CAUTION: Inspect the puncture wire provided in this set to be certain the wire tip does not protrude beyond the sheath but is about even with the end of the sheath. (The pin vise actuator is locked on the wire when shipped.) It is equally important that the puncture wire point not be retracted more than 1 mm from the tip of the sheath.

6. Aided by an assistant, carefully advance the puncture wire (Figure D) into the catheter until the puncture wire tip is observed exiting the tip of the outer catheter. At this point, careful guidance is necessary to be certain the puncture wire extension is positioned in the desirable calyx. Careful angling of the patient or the film plane may be advisable to confirm proper positioning.

Figure D

NOTE: The outer catheter may back out of the selected calyx or “cobra” away from the calyx as the puncture wire approaches or exits the catheter tip. Once the puncture wire is advanced carefully and as far into the catheter as possible, the pin vise actuator is twist-locked onto the catheter hub for added stability (Figure E).

Figure E

7. Firmly grasp the cap of the pin vise actuator between the thumb and index finger (Figure F), and free the puncture wire by un-twisting the handle counterclockwise.


Figure F

Maintaining the thumb/index finger grasp, slowly advance the fine wire to perform the puncture into and through parenchyma, capsule, fat, muscle, fascia, and skin (Figure G). **NOTE:** The wire point may “tent” just beneath the skin exit site. A scalpel nick is sometimes necessary.


Figure G

CAUTION: Close A-P and lateral fluoroscopic positioning of the catheter tip and guidance of the advancing puncture wire is necessary to avoid undesirable passage into or through the adjacent viscera (especially the bowel, liver, spleen, lung, twelfth rib, or intercostal space).

8. After puncturing the skin, advance the wire from below (within) until approximately 30 cm of wire is externally visible.

CAUTION: The wire must not be pulled out since this could lacerate the parenchyma or act to “saw” the puncture tract. An exit site anterior to the mid-axillary line must not be used to develop a tract. Dilation of such a tract significantly increases the risk of puncturing adjacent viscera. The wire must be withdrawn and the catheter repositioned so the puncture wire will exit in the mid-to-lower-posterior quadrant.

9. Stabilize the wire at the skin surface with a hemostat. Remove the stylet of the assembled needle set. Advance the coaxial needle assembly over the puncture wire until the 22 gage needle tip is observed in the collecting system and is positioned end to end with the tip of the ureteral catheter.
10. Unlock the 18 gage needle with a counterclockwise twist; advance the needle over the puncture wire obturator 22 gage needle until the external etch mark indicates the needle tips are congruent.
11. Carefully maintaining needle-point-to-catheter-tip contact and allowing no gap, retract the puncture wire antegrade well into the catheter. Remove the catheter cystoscopically only after the wire is completely sheathed in the 7.0 French ureteral catheter. It is essential to maintain close approximation of the needle tip and 7.0 French catheter tip under fluoroscopy as the puncture wire is withdrawn back into the catheter.

CAUTION: Failure to observe this procedure could result in serious injury to the kidney by laceration of parenchyma with the puncture wire.

12. Remove the 22 gage needle. Position the .038 inch (0.97 mm) diameter “working” wire guide through the 18 gage needle, well within the renal pelvis and upper ureter. Dilate the punctured tract. Introduce an additional safety wire guide in preparation for the planned procedure.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

- K. W. Kaye, D. B. Reinke: "Detailed Calyceal Anatomy for Endourology," Journal of Urology, 132 (1984), 1085-1088.
- R. K. Lawson, J. B. Murphy, A. J. Taylor, S. C. Jacobs: "Retrograde Method for Percutaneous Access to Kidney," Urology, 22 (1983), 580-582.
- R. K. Lawson, M.D., The Medical College of Wisconsin, Froedtert Memorial Lutheran Hospital, Milwaukee, Wisconsin.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd