**EN**  
**3**

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**Uretärlaserkateter**  
Bruksanvisning

## LASER URETERAL CATHETER

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.**

### DEVICE DESCRIPTION

The Laser Ureteral Catheter is constructed of radiopaque tubing that is available in 7.1 French size, with 40 cm in length. The catheter will accept a laser fiber up to an overall outer diameter of 1000  $\mu\text{m}$ .

### INTENDED USE

The Laser Ureteral Catheter is intended for protection and delivery of a laser fiber in the urinary tract.

### CONTRAINDICATIONS

This device is contraindicated for use in patients unsuitable for catheterization.

### WARNINGS

- Do not use if package is open/damaged.
- Do not extend fiber beyond catheter prior to insertion of catheter.
- Laser fiber should not be fired unless fiber tip extends beyond catheter device.

### PRECAUTIONS

- This device is intended for use by physicians trained and experienced in placement of ureteral catheters. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- Stop catheter manipulation if you encounter resistance when advancing or withdrawing the catheter. Determine the cause of the resistance before proceeding.

### POTENTIAL ADVERSE EVENTS

Adverse events that may result from ureteral catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

### INSTRUCTIONS FOR USE

1. Remove the product protector from the laser ureteral catheter.
2. Preload the laser fiber of choice into the catheter so that the laser fiber tip is just proximal to the distal end of the laser ureteral catheter. Do not advance the laser fiber beyond the distal tip of the laser ureteral catheter. Tighten the connector cap.

3. With a previously placed scope of choice in the desired location, use direct vision to advance the laser ureteral catheter and laser fiber as one unit to the desired location beyond the distal end of the scope.
4. Loosen the connector cap to allow advancement of the laser fiber to the desired location beyond the laser ureteral catheter. NOTE: to manipulate the laser ureteral catheter and laser fiber as one unit you can tighten the connector cap.
5. Advance the laser fiber so that it extends beyond the distal tip of the catheter.
6. Operate the laser fiber per the manufacturer's instructions for use.

## REMOVAL INSTRUCTIONS

To remove the device, withdraw the laser fiber back inside the laser ureteral catheter, and then remove the fiber/catheter assembly from the scope.

**NOTE:** If the connector cap was tightened during the use of the laser fiber you will need to loosen the connector cap to allow the laser fiber to be withdrawn back into the laser ureteral catheter.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.