COOK

MEDICAL

CE
0088

EN
4

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- 1. Orificio lateral
- 1. Orifice latéral
- 1. Foro laterale
- 1. Zijopening

- 1. Orificio lateral
- 1. Sideport

MARX-COPE GASTROJEJUNOSTOMY CATHETERS AND SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Marx-Cope Gastrojejunostomy Set consists of a pigtail gastrostomy catheter, and a jejunostomy catheter with hydrophilic coating on the distal 50 cm. The set also contains appropriately sized dilators, a Peel-Away® introducer sheath, an introducer needle, a wire guide, and accessories for percutaneous placement of a gastrojejunostomy catheter.

Sets also include a Cope Gastrointestinal Suture Anchor Set, consisting of an introducer needle with a preloaded suture anchor, a special .018/.035 inch wire guide with a spring coil tip, and a second suture anchor in a loading cannula.

INTENDED USE

The Marx-Cope Gastrojejunostomy Catheters and Sets are intended for percutaneous placement of a coaxial catheter system for gastroenteric feeding.

The Cope Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in gastrojejunostomy techniques. Standard techniques for percutaneous placement of gastrojejunostomy catheters should be employed.
- Manipulation of the product requires fluoroscopic control.
- The use of TFE-coated wire guides for catheter placement is recommended.
- The Peel-Away introducer in this set has a uniquely fabricated transition between dilator and sheath. Do not turn or pull back dilator from sheath prior to introduction. Manipulation of dilator will result in misalignment of sheath and dilator. Do not attempt to flush sheath prior to insertion.
- Activate hydrophilic coating on jejunostomy catheter by wetting catheter surface with sterile water or saline. For best results, maintain wetted condition of device during placement.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Placement of the Cope Gastrointestinal Suture Anchor Set

1. Insufflate stomach with approximately 500 cc of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Attach a 10 ml syringe, half filled with contrast medium, to the introducer needle that has been preloaded with the suture anchor.
4. Advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. **Air should be freely aspirated. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 ml syringe from the needle hub. While maintaining slight tension on the trailing suture, introduce the distal spring coil .035 inch portion of the wire guide into the needle and use it to push the suture anchor out of the needle into the stomach cavity. **Maintaining slight tension on the suture during anchor introduction will prevent it from bunching and lodging within the needle lumen.**
6. Remove the introducer needle over the wire guide and trailing suture needle. **The .018 inch part of the wire guide must be positioned adjacent to the suture needle to allow clearance and removal of the introducer needle.**
7. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
8. While maintaining traction on the suture anchor, secure the suture to the skin under slight tension.

9. Remove the wire guide.
10. To place the second suture anchor, reintroduce the introducer needle as described in step 4.
11. Advance the loading cannula containing the second suture anchor into the hub of the introducer needle.
12. Advance the wire guide through the cannula into the introducer needle, pushing the second suture anchor into the abdominal cavity.
13. Follow steps 6 through 8 with the second suture anchor.
14. The wire guide may now be used to facilitate passage of fascial dilators. Take care not to advance the entire distal spring coil segment of the wire guide into the abdominal cavity.

NOTE: The suture may be left in place for two weeks while tract formation occurs, or it may be cut after gastrostomy catheter placement. This releases the anchor into the stomach, allowing its passage via the gastrointestinal system.

Placement of Gastrojejunostomy Catheter

1. Using the supplied introducer needle, make a puncture of the stomach adjacent to the suture anchor. **NOTE:** This puncture should be performed with the needle oriented at a 45-degree angle toward the pylorus.
2. Advance the supplied extra-stiff wire guide through the introducer needle and allow it to coil in the antrum of the stomach.
3. Remove the introducer needle.
4. Using the supplied dilators, sequentially dilate up to 12.0 French. **NOTE:** Ensure that each dilator is inserted through the gastric wall and into the gastric lumen. Fluoroscopic observation and maintenance of tension on the suture anchor during insertion can help ensure intraluminal position.
5. Remove the dilator and advance the Peel-Away introducer over the wire guide.
6. Ensure that the entire introducer is within the stomach. Again, fluoroscopic visualization and maintenance of tension on the suture anchor during insertion can help ensure intraluminal position.
7. Turn the dilator one-half turn to break the seal between sheath and dilator.
8. Pull back the dilator to remove from sheath, leaving the sheath and wire guide in place.
9. Advance the pigtail gastrostomy catheter over the wire guide through the sheath.
10. Under fluoroscopic guidance, withdraw the wire guide until its tip is immediately proximal to the sideport on the gastrostomy catheter's outer curvature. Allow the catheter to form a pigtail, and secure the loop in appropriate fashion with the retention suture.
11. The Peel-Away sheath may be removed at this point or, alternatively, left in place to add stiffness to the transperitoneal tract during catheterization of the duodenum.
12. Under fluoroscopic guidance, position the pigtail of the gastrostomy catheter so the sideport on its outer curvature is pointed toward the pylorus. Advance the wire guide through the sideport. **NOTE:** Alternately, the seeking catheter can be advanced over the wire guide and the catheter tip can be used to exit the sideport. (**Fig. 1**)
13. If the seeking catheter has not already been introduced, advance it over the wire guide. Using the seeking catheter and wire guide as a unit to traverse the pylorus and duodenum, position the tip of the wire guide beyond the Ligament of Treitz. **NOTE:** Non-ionic contrast and air may be injected through the gastric sideport of the gastrostomy catheter to aid in visualizing the pylorus.
14. Remove the seeking catheter, leaving the wire guide in place, and flush the gastrostomy catheter with saline through the sideport.
15. Activate the hydrophilic coating on the jejunostomy catheter by wetting the distal half of the catheter with sterile water or saline. The catheter surface will become very lubricious.
16. Advance the jejunostomy catheter over the wire guide through the gastrostomy catheter until its tip is beyond the Ligament of Treitz.
17. Luer lock the jejunostomy catheter onto the gastrostomy catheter. Flush and cap both ports. **NOTE:** Alternately, the gastric port may be left open to gravity drainage for gastric decompression.
18. If the Peel-Away sheath has been left in place, remove it at this point.
19. If suture anchors have been placed, the suture may be cut to release the anchors into the stomach. **NOTE:** Alternately, sutures may be left in place for two weeks to ensure tract formation.
20. Secure the gastrostomy catheter with the supplied fixation disc. After applying the fixation disc, flush both ports to ensure that the disc is not pulled too tightly. A too-tight pull tie may occlude the jejunal catheter lumen. **NOTE:** If using a Luer tip syringe for flushing purposes, remove feeding adapter and connect Luer converter to the hub of the jejunal catheter. Remove Luer converter and connect the feeding adapter after flushing the catheter.

Gastrojejunostomy Catheter Removal or Exchange

1. Remove the jejunostomy tube, if present.
2. Under fluoroscopic guidance, advance the wire guide into the gastrostomy catheter, ensuring the wire guide exits the distal catheter tip.
3. Uncover the suture by folding back the latex sleeve, and release the suture.
4. Withdraw the catheter over the wire guide.

NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open pouches. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.