- EN**
2 Multi-Sideport Catheter Infusion Set with Dual Check Valve
Instructions for Use
- DA**
6 Kateterinfusionssæt med multi-sideporte og dobbeltcheckventil
Brugsanvisning
- DE**
10 Infusionsset mit Multi-Seitenloch-Katheter und Doppelrückschlagventil
Anleitung zum Gebrauch
- EL**
13 Σετ έγχυσης καθετήρα πολλαπλών πλευρικών θυρών με ανεπίστροφη βαλβίδα
Οδηγίες χρήσης
- ES**
17 Equipo de catéter con múltiples orificios laterales para infusión con válvula de no retorno
Instrucciones de uso
- FR**
21 Set de cathéter de perfusion multiperforé à double valve anti-retour
Mode d'emploi
- HU**
24 Több oldalnyílással rendelkező katéteres infúziós készlet kettős visszacsapó szeleppel
Használati utasítás
- IT**
28 Set per infusione dotato di catetere con fori laterali multipli e valvola di non ritorno
Istruzioni per l'uso
- NL**
32 Infusiekatheterset met meerdere zijopeningen en dubbele keerklep
Gebruiksaanwijzing
- PT**
35 Conjunto de perfusão com cateter com múltiplos orifícios laterais e válvula de segurança unidirecional dupla
Instruções de utilização
- SV**
39 Kateterinfusionssats med flera sidoportar och dubbelspärventil.
Bruksanvisning

MULTI-SIDEPORT CATHETER INFUSION SET WITH DUAL CHECK VALVE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

MULTI-SIDEPORT CATHETER INFUSION ASSEMBLY


Fig. 1

The Multi-Sideport Catheter Infusion Set consists of a catheter with distal and proximal radiopaque markers indicating the location of the infusion segment, a Tuohy-Borst side-arm adapter, a catheter tip occluder with filiform tip and a dual check valve with syringe. (**Fig. 1**) The dual check valve is designed to facilitate uncomplicated loading of an infusate from a larger syringe source to a small 1 mL delivery syringe. This permits subsequent injection without disconnection from the delivery catheter assembly. The use of the catheter tip occluder assembly enhances controlled localized dispersion of infusate through the sideports of the infusion catheter.

Infusion Segment Length/Sideports

7.0 cm long	14 sideports spiraled over 7.0 cm length
15.0 cm long	30 sideports spiraled over 15.0 cm length
20.0 cm long	40 sideports spiraled over 20.0 cm length

INTENDED USE

Intended to administer infusions of various therapeutic solutions and contrast media into the peripheral vasculature. (Refer to manufacturer's package insert for complete prescribing and dosing information.)

CONTRAINDICATIONS

None known

WARNINGS

- Do not inject contrast with a power injector through the Tuohy-Borst adapter, or when the tip occluder is in place. Damage to the vessel may occur.
- Due to the thinwall construction of the multi-sideport catheter, care must be exercised during manipulation and withdrawal to prevent catheter damage.
- Insertion through the wall of a synthetic vascular graft should be avoided. If introduction through the wall of a synthetic vascular graft is necessary, an introducer sheath must be used to prevent catheter damage.
- Ensure that all connections are secure before use. Do not over tighten as excessive force may damage the product.
- All therapeutic agents to be infused must be according to the manufacturer's instructions of use.
- This procedure should be repeated for the duration of therapy as directed by the physicians.

PRECAUTIONS

- Manipulation of products requires fluoroscopic control.
- Do not insert, manipulate or withdraw the multi-sideport catheter without the use of a wire guide.
- Do not use a power injector with the dual check valve. Damage to the dual check valve may occur.
- If resistance is felt between the vessel and the catheter, or between the catheter and the wire guide, do not attempt to use force, as the vessel, the catheter or the wire guide may become damaged.
- Verify the cause of resistance by fluoroscopy, or remove the catheter and wire guide together. If there is any doubt whether the catheter or wire guide has been damaged, use a new catheter or wire guide.
- The product is intended for use by physicians trained and experienced in selective catheterization techniques and interventional procedures. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

INSTRUCTIONS FOR USE

1. Prepare and drape the access site using standard procedure.
2. Introduce the infusion catheter over the wire guide to the desired catheter position.
3. With catheter at desired position, remove the wire guide and introduce the tip occluder into the catheter. Advance the tip occluder to the tip of the catheter until resistance is felt, indicating occlusion of the catheter.


Fig. 2

4. Connect the dual check valve to the clear Tuohy-Borst adapter. The Tuohy-Borst valve should be closed completely. Connect a 20 mL heparinized saline-filled syringe to the sideport of the dual check valve adapter and slowly purge air from the injection port. (Fig. 2)


Fig. 3

5. Purge the remaining air by placing a thumb over the injection port. Slowly flush the remaining air bubbles out of the dual check valve and Tuohy-Borst adapter. (Fig. 3)


Fig. 4

6. Attach the 1 mL syringe to the injection port and aspirate 0.5 mL of solution. Detach the 1 mL syringe and purge all air bubbles from the syringe. Leave a meniscus of solution on the injection port before the 1 mL syringe is reattached by depressing the 20 mL syringe slightly to prime the injection port. Reattach the 1 mL syringe to the injection port. (Fig. 4)


Fig.5

7. To attach the dual check valve and Tuohy-Borst adapter to the infusion catheter, loosen the Tuohy-Borst valve and slide the assembly over the tip occluder. Before attaching the assembly to the catheter, hold the 1 mL syringe plunger stationary and depress the 20 mL syringe, making sure all air is purged from the assembly. Leave a meniscus of solution on the Luer lock fitting of the Tuohy-Borst adapter before connecting to the infusion catheter. **(Fig. 5)**


Fig.6

8. Purge any remaining air bubbles by holding the plunger on the 1 mL syringe stationary, and slowly depressing the 20 mL syringe to force any remaining air bubbles out. Once all air has been removed from the system, tighten the Tuohy-Borst valve around the tip occluder. **(Fig. 6)**


Fig.7

9. To start the procedure for desired administration of therapeutic solution, remove the 20 mL syringe with heparinized saline solution, leaving a meniscus of solution on the sideport fitting. Before attaching another fluid-filled syringe to the sideport, make sure all air bubbles are purged. Reattach with a therapeutic solution-filled syringe. **(Fig. 7)**


Fig.8

10. To infuse solution, draw desired dose into 1 mL syringe and inject. Reload 1 mL syringe and repeat infusion technique as necessary. (Fig. 8)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.