**EN**  
**3** **Malecot Nephrostomy Catheter/Stent Set**

Instructions for Use

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**SV**  
**25** **Malecot nefrostomikateter/stentset**

Bruksanvisning

## MALECOT NEPHROSTOMY CATHETER/STENT SET

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

### DEVICE DESCRIPTION

The Malecot Nephrostomy Catheter/Stent Set consists of:

- Radiopaque Malecot catheter/stent and two flexible stylets
- Retention disc with pull tie
- Connecting tube

**NOTE:** Set components may vary.

**NOTE:** Reserve the individually packaged, flexible stylet for straightening the Malecot catheter before removal at a later date.

### INTENDED USE

The Malecot Nephrostomy Catheter/Stent Set is intended for use as a nephrostomy drainage catheter and ureteral stent. The catheter is placed percutaneously through an existing nephrostomy tract.

### MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Malecot Nephrostomy Catheter/Stent (including the retention disc with pull tie and connecting tube) is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions. **NOTE:** The flexible stylets have not been evaluated for use in an MR environment.

- Static magnetic field of 3.0 Tesla or 1.5 Tesla only
- Maximum magnetic field spatial gradient of 1900 Gauss/cm (19 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of  $\leq 2.0$  W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Malecot Nephrostomy Catheter/Stent is not expected to experience RF induced heating above background heating.

The image artifact extends approximately 0.20 mm from the Malecot Nephrostomy Catheter/Stent as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

This product is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.

### POTENTIAL ADVERSE EVENTS

- Catheter occlusion and/or dislodgment
- Encrustation
- Hemorrhage or hematoma

- Infection
- Pain

## INSTRUCTIONS FOR USE

1. Perform standard techniques for establishing a nephrostomy tract under fluoroscopic guidance.
2. Leave a wire guide in the collecting system, and ensure fluoroscopically that the wire guide travels the length of the ureter and is well into the bladder. **NOTE:** The catheter/stent will accept a wire guide up to .038 inch in diameter.
3. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract. If you are replacing a catheter that is already in place, skip this step and proceed to the next step.
4. Place the blunt stylet inside the Malecot catheter and secure its position with the Luer lock to straighten the Malecot wings.
5. Pass the stent end of the Malecot catheter over the external end of the wire guide. Gradually advance the catheter end into the collecting system and place the desired length of the stent segment into the ureter. Confirm position fluoroscopically.
6. Release the Luer lock to open the Malecot wings, then carefully remove the blunt stylet. **NOTE:** Take care to avoid withdrawing the wire guide itself.
7. While holding the shaft of the Malecot catheter securely in position with one hand, withdraw the wire guide with the other hand. When you have confirmed appropriate position of the Malecot via fluoroscopic visualization, remove the wire guide.
8. Tape or suture the retention disc to the skin, then secure the catheter to the disc using the plastic pull tie. A dressing may be applied.
9. Use the connecting tube to connect the Malecot catheter to a drainage bag or leg bag.

## Catheter Removal

1. Place the straightening stylet into the catheter.
2. If desired, pass a wire guide into the collecting system to maintain access after the stylet is engaged and before catheter removal.
3. Remove the catheter.

**CAUTION:** The indwelling time for this catheter should not exceed four weeks.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd