COOK

MEDICAL

CE
0088

EN
3

Maglinte Enteroclysis Catheter Set

Instructions for Use

CS
4

Souprava katetru Maglinte pro enteroklýzu

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DA
5

Maglinte enteroclyse-katetersæt

Brugsanvisning

DE
6

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T _ M E C _ R E V 5

MAGLINTE ENTEROCLYSIS CATHETER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Maglinte Enteroclysis Catheter Set is a polyvinylchloride balloon catheter with a TFE-coated stainless steel torque cable. The catheter can be introduced nasally or orally.

INTENDED USE

The Maglinte Enteroclysis Catheter Set is intended to administer contrast medium to the small intestine.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in interventional and gastroenterology techniques. Standard techniques for placement of nasogastric tubes should be employed.
- Manipulation of product requires fluoroscopic control.
- A topical anesthetic spray to the throat is helpful, but over-anesthetization of the pharyngeal region should be avoided so that the patient can swallow during the intubation process.
- It is recommended that the patient be premedicated with a prokinetic agent prior to the procedure.
- At no time should the catheter be forced against resistance.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Nasal Intubation

1. Place the patient in a supine position on the fluoroscopic table.
2. Introduce the torque cable into the catheter, then withdraw the torque cable 5 cm to make the catheter tip flexible.
3. Under fluoroscopy, advance the catheter through the nares into the oropharynx.
4. Have the patient flex his or her neck. Advance the catheter gently while the patient swallows. **NOTE:** If intubation is performed in the recumbent position, and resistance is met, place the patient in the right lateral position following passage of the catheter tip beyond the pharynx.
5. Continue inserting the catheter. As the catheter is advancing, retract the tip of the torque cable so that it does not go beyond the pyloric aperture.
6. Occasionally, the catheter will coil in the fundus of the stomach. When this happens, withdraw the catheter to the gastroesophageal junction and reinsert the angled tip of the torque cable for directional control.
7. If possible, position the balloon of the catheter in the horizontal duodenum. If this is not possible, position the balloon anywhere in the duodenum or proximal jejunum.
8. When the desired position of the balloon is achieved, introduce 15 ml of air through the air lumen. Perform fluoroscopy to check the size of the balloon. **NOTE:** 15 ml of air will distend the balloon to approximately 3 cm in width. To prevent bowel wall necrosis, avoid prolonged use of an over-distended balloon.
9. Tape the proximal end of the catheter in place, as increased peristalsis can carry the catheter tip distally.
10. After evaluation of the small intestine, deflate the balloon and withdraw the catheter.

Oral Intubation

1. With the patient sitting up, apply a light anesthetic spray to the throat.
2. Introduce the torque cable into the catheter, then withdraw the torque cable 5 cm to make the catheter tip flexible.
3. With the patient's neck hyperextended, quickly introduce the catheter while the patient swallows continuously.
4. If resistance is met, stop intubation and perform fluoroscopy. Reposition the patient if necessary.
5. From this point, oral passage of the catheter is similar to nasal passage. Follow steps 5-10 in the Nasal Intubation section.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.