EN
3

Multi-Length Stent Set

Instruction for Use

CS
4

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
Figure B

MULTI-LENGTH STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Set contains (*set components may vary*)
- Radiopaque stent
- Radiopaque stent positioner
- Radiopaque ureteral catheter
- Appropriate diameter stainless steel wire guide

INTENDED USE

Used for temporary internal drainage from the ureteropelvic junction to the bladder. The multi-length stents will accommodate ureters from approximately 8-32 cm long.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

None known

PRECAUTIONS

- Complications of ureteral stent placement are documented in literature. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures
- The Endo-Sof™ stents must not remain indwelling more than twelve (12) months. The Sof-Flex®, LSe Sof-Flex®, C-Flex®, and AQ® stents must not remain indwelling more than six (6) months. If the patient's status permits, the stent may be replaced with a new stent.
- These stents are not intended as permanent indwelling devices.
- Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- Prior to use, immerse AQ® stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

RECOMMENDED WIRE GUIDE

Stent Size	Accepts Wire Guide	Stent Size	Accepts Wire Guide
3.0 French	.018 inch	7.0 French	.038 inch
3.7 French	.025 inch	8.0 French	.038 inch
4.0 French	.025 inch	8.2 French	.038 inch
4.7 French	.035 inch	10.0 French	.038 inch
5.0 French	.035 inch	12.0 French	.038 inch
6.0 French	.038 inch		

FOR USING MULTI-LENGTH STENT SET

1. Pass the flexible wire guide tip beyond the obstruction to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using the wire guide and an open-end ureteral catheter in combination. Remove catheter before attempting to place stent.
2. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
3. When the renal pelvis has been entered (**Figure A**) (this can be confirmed fluoroscopically), gently remove the wire guide while firmly holding the positioning catheter. The stent coil will form spontaneously in the renal pelvis. Carefully remove the positioning catheter (**Figure B**).
 - For ureters longer than 22 cm, endoscopic forceps can be used to adjust the stent position, leaving approximately equal coils in the renal pelvis and bladder.

NOTE: If problems occur using this device, please call your Cook Urological® sales representative or contact our Customer Quality Assurance department at the address/phone listed within this booklet.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

M. F. Camacho, R. Pereira, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 31 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

C-FLEX® is a registered trademark of Concept Polymer Technologies, Inc.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.