EN 3 Malecot Nephrostomy Catheter Set

Instructions for Use

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T _ M N C S _ R E V 1

MALECOT NEPHROSTOMY CATHETER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Wire guide
- Radiopaque dilators (optional)
- Radiopaque Malecot catheter and flexible stylet
- Retention disc with pull tie
- Connecting tube

NOTE: Set components may vary.

INTENDED USE

Used for nephrostomy drainage in the renal pelvis. The Malecot catheter is placed percutaneously through an existing nephrostomy tract.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The indwelling time of this device should not exceed four weeks. Evaluation and replacement beyond four weeks must be considered.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Perform the standard techniques for establishing a nephrostomy tract under fluoroscopic control, leaving a Safe-T-J® wire guide in the collecting system.
2. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided. If replacing a Malecot catheter for another, skip this step and proceed with step 3.
3. Place the blunt stylet inside the Malecot catheter and secure its position with the Luer lock to straighten the Malecot wings.
4. Pass the Malecot end of the catheter over the external end of the wire guide. Gradually advance the catheter end well into the collecting system. Confirm the position fluoroscopically. After releasing the Luer lock in order to open the Malecot wings, carefully remove the blunt stylet. **NOTE:** Care must be exercised to avoid withdrawal of the wire guide itself.
5. While holding the shaft of the Malecot catheter securely in position with one hand, withdraw the wire guide with the other. When the appropriate position of the Malecot is assured by fluoroscopic visualization, remove the wire guide.
6. Tape or suture the provided retention disc to the skin; secure the catheter to the disc with the plastic pull tie. A dressing may be applied. Use the connecting tube to connect the Malecot catheter to a drainage bag or leg bag.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.