EN 3 Percutaneous Malecot Nephrostomy Set

Instructions for Use

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Bruksanvisning

T _ M N S _ R E V 1

PERCUTANEOUS MALECOT NEPHROSTOMY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Disposable needles, 20 cm long
- .038 inch (0.97 mm) diameter TFE-coated stainless steel wire guide, 80 cm long, with 3 mm Safe-T-J® tip
- Radiopaque dilators, 20 cm long
- Radiopaque Malecot catheter and two flexible stylets
- Silicone retention disc with pull tie
- Clear polyvinylchloride connecting tube, 30 cm long, with one-way stopcock and drainage bag connector

NOTE: Set components may vary.

INTENDED USE

Used for percutaneous placement of a Malecot catheter in the renal pelvis.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The indwelling time of this device should not exceed four weeks. Evaluation and replacement beyond four weeks must be considered.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Using preliminary plain film, I.V.P., retrograde pyelogram, ultrasound or CAT scan, localize the kidney to be drained. A prone or slightly oblique position is preferred. Prepare and drape the flank in the usual fashion. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.
2. Introduce the 22 gage (black hub) needle into the collecting system at a right angle to the skin. Obtain a urine sample, if desired. Visualization of the pyelocalyceal system may be further enhanced by the injection of contrast medium.

NOTE: Maintain the 22 gage needle at this location and position throughout the remainder of the procedure. This fixes the kidney in position and allows further injection of the pyelocalyceal system, if desired.

3. Rotate the patient into a more lateral position with the posterior axillary line of the ipsilateral flank in the superior position. Pass the 18 gage (pink hub) needle vertically through a small incision in the skin (made with a scalpel blade); insert the needle into the appropriate part of the pyelocalyceal system. Free flow of fluid from the needle after removing the obturator confirms a satisfactory location for the tip of the needle cannula.
4. Pass the flexible end of the Safe-T-J® wire guide through the hub of the 18 gage needle cannula; insert the wire guide three more inches into the collecting system. Confirm the position fluoroscopically. Maintain the position of the wire guide as the needle cannula is withdrawn over it. **NOTE:** Care must be exercised to avoid withdrawal of the wire guide itself.
5. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided.
6. Place the blunt stylet inside the Malecot catheter and secure its position with the Luer lock to straighten the Malecot wings.

7. Pass the Malecot end of the catheter over the external end of the wire guide; gradually advance the catheter end well into the collecting system. Confirm the position fluoroscopically. After releasing the Luer lock in order to open the Malecot wings, carefully remove the blunt stylet.
8. While holding the shaft of the Malecot catheter securely in position with one hand, withdraw the wire guide with the other. When the appropriate position of the Malecot is assured by fluoroscopic visualization, remove the 22 gage needle.
9. Tape or suture the provided retention disc to the skin; secure the catheter to the disc with the plastic pull tie. A dressing may be applied. Use the connecting tube to connect the Malecot catheter to a drainage bag or leg bag. The one-way stopcock is useful for testing the drainage capability of the collecting system, if desired.

NOTE: The individually packaged, flexible catheter stylet should be reserved for straightening the Malecot catheter prior to removal at a later date.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.