COOK**MEDICAL****CE**
0088

EN 4	Peripherally Inserted Central Venous Catheter Sets and Trays With Micropuncture® Peel-Away® Introducers Instructions for Use
DA 9	Perifert indførte centralvenekatetersæt og bakker med Micropuncture® Peel-Away® indførere Brugsanvisning
DE 15	Sets mit peripher einführbaren zentralen Venenkathetern mit Micropuncture® Peel-Away® Einführschleusen Gebrauchsanweisung
EL 21	ΣΕΤ και δίσκοι ΚΕΝΤΡΙΚΩΝ ΦΛΕΒΙΚΩΝ ΚΑΘΗΤΗΡΩΝ που εισάγονται ΠΕΡΙΦΕΡΙΚΑ με εισαγωγείς Micropuncture® Peel-Away® Οδηγίες χρήσης
ES 28	Equipos y bandejas de catéteres venosos centrales de inserción periférica con introductores Micropuncture® Peel-Away® Instrucciones de uso
FR 34	Sets et plateaux de cathéter veineux central inséré en périphérie avec introducteurs Micropuncture® Peel-Away® Mode d'emploi
IT 40	Set standard e set completi con catetere venoso centrale ad inserimento periferico con introduttori Micropuncture® Peel-Away® Istruzioni per l'uso
NL 47	Perifeer ingebrachte centraal-veneuze kathetersets en -pakketten met Micropuncture® Peel-Away® introducers Gebruiksaanwijzing
PT 53	Conjuntos e tabuleiros de catetere venoso central de introdução periférica com introdutores Micropuncture® Peel-Away® Instruções de utilização
SV 60	Set och brickor med perifert inlagd central venkateter med Micropuncture® Peel-Away®-införare Bruksanvisning
ZH 65	带 Micropuncture® Peel-Away® 撕脱式导入器的外周插入式中心静脉导管组套和托盘 使用说明

T _ M P I C _ R E V 4


Fig. 1


Fig. 2

PERIPHERALLY INSERTED CENTRAL VENOUS CATHETER SETS AND TRAYS

With Micropuncture® Peel-Away® Introducers

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Silicone Peripherally Inserted Central Venous Catheters (PICCs) are available in single-and double-lumen configurations and four lengths; they can be trimmed to fit the patient's anatomy. The catheters feature a soft silicone suture wing for secure fixation to the skin and reinforced silicone extensions to reduce catheter damage during catheter manipulation. The catheters may be ordered with plastic clamp(s), which help prevent air aspiration due to inadvertent hub dislodgement.

PICC sets and trays also contain a Peel-Away introducer, Micropuncture EchoTip® needle, wire guides, and other accessories for percutaneous vascular placement.

Single-Lumen Catheters

French Size	Lumen Volume	Lumen Diameter	Length
3.0	0.5 ml	.022 inch	50 cm
4.0	0.6 ml	.026 inch	60 cm
5.0	0.8 ml	.036 inch	60 cm

Double-Lumen Catheters

French Size	Small Lumen Diameter	Small Lumen Volume	Large Lumen Diameter	Large Lumen Volume	Length
6.0	.018 inch	0.5 ml	.026 inch	0.6 ml	58 cm
7.0	.019 inch	0.7 ml	.025 inch	0.8 ml	60 cm
7.0 ¹	.019/.045 inch	0.7 ml	.025/.050 inch	0.8 ml	60 cm

¹ D-shaped lumens

INTENDED USE

Silicone Peripherally Inserted Central Venous Catheter (PICC) Sets and Trays are intended for delivery of whole blood or blood products, drug administration, and blood sampling.

CONTRAINDICATIONS

None known

WARNINGS

- Silicone Peripherally Inserted Central Venous Catheters are **not** designed for power injection of contrast medium. Catheter rupture may occur. Use of a 10 ml or larger syringe will reduce the risk of catheter rupture.
- **Peripherally Inserted Central Venous Catheters play an important role in treatment of critically ill patients. However, catheter tips can erode or perforate vascular walls. Extreme caution must be used in placement and monitoring of catheters.**
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system. Catheter tip position should be verified by X-ray and monitored on a routine basis. Periodic lateral view X-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall. **(Reference 1)**

PRECAUTIONS

- This product is intended for use by health care practitioners trained and experienced in proper positioning of catheters in the central venous system using percutaneous entry (Seldinger) technique. Standard technique for placement of vascular access sheaths, catheters and wire guides should be employed.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Patient movement can cause catheter tip displacement. Catheters placed via an antecubital vein have shown tip movement of up to 10 cm with motion of the extremity.
- Catheter size should be as small as the use will allow.

PRODUCT RECOMMENDATIONS

Catheter Size and Puncture Site


Preliminary reports indicate that catheter size can influence clotting. Larger diameter catheters have more tendency to promote clots. As reported by Amplatz, Gianturco and others, clot formation has less relation to type of catheter material than to size of catheter.

The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming X-rays, and found that an incident angle of the catheter to vessel wall greater than 40 degrees was more likely to perforate. **(Reference 2)**

The following variables must also be considered in selecting appropriate catheter and length:

1. Patient history
2. Patient age and size
3. Access site available
4. Unusual anatomical variables
5. Proposed use and duration of treatment plan

Catheter Tip Positioning

Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. **Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.**

Catheter Maintenance

Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization. After catheter placement and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, catheter tip position should be immediately reevaluated by physician.** If catheter is not to be used immediately, its lumen should be maintained by continuous saline or heparinized saline drip or locked with heparinized saline solution. **NOTE:** If CLC-2000, MicroClave or other needleless endcaps approved for saline only lock are used, saline only catheter lock may be used. Catheter heparinization should be determined by institutional protocol and clinical judgement. Heparin concentrations of 10 Units/ml to 100 Units/ml have been reported adequate to maintain lumen patency. **(Reference 3)** Heparin lock should be reestablished after every use or at least every 24 hours if unused. Before using catheter lumen already locked with

heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumen should be flushed with normal saline between administration of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing heparin lock.

INSTRUCTIONS FOR USE

Catheter Obturator Preparation

1. Flush catheter with heparinized saline solution or sterile water. **NOTE:** The catheter may be trimmed if a shorter length is required.
2. Attach a syringe with heparinized saline solution or sterile water to the Luer lock fitting of the catheter obturator holder.
3. Inject enough solution to wet the obturator surface entirely. This will activate the AQ® hydrophilic coating, making the obturator surface very lubricious.
4. Remove the obturator from its holder and insert it into the catheter, locking it into position. If the catheter has been trimmed, only advance the obturator to the distal end of the catheter. **NOTE:** If the surface of the obturator becomes dry after removal from the holder, wetting with additional heparinized saline or sterile saline will renew the hydrophilic effect.
5. The catheter/obturator assembly can now be introduced as described in the following section.

Catheter Placement (Fluoroscopic Method)

6. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to help visualize the tip of the needle during vessel access.
7. Using fluoroscopic guidance, introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
8. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.
9. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
10. Using fluoroscopic control, determine the correct catheter length by advancing the wire guide to the desired catheter tip location. Once the wire guide tip is in proper position, mark the length by clamping forceps onto the wire guide at the skin site.
11. Withdraw the wire guide and measure it from the forceps to the distal tip to determine correct catheter length. Trim catheter if necessary. **NOTE:** Remove obturator prior to trimming the catheter; reinsert for catheter introduction.
12. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the cuffed proximal end of the sheath.
13. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 7 cm distal to the catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter.

14. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
15. After catheter is in final position, remove obturator, secure catheter to skin and dress in standard fashion.
16. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.

Catheter Placement (Non-Fluoroscopic Method)

1. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to help visualize the tip of the needle during vessel access.
2. Introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
3. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.
4. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
5. Using a Tyvek® tape measure, clinical judgment, or other institutional protocol, determine the correct catheter length and trim catheter as needed. **NOTE:** Remove obturator prior to trimming catheter; reinsert for catheter introduction.
6. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the proximal end of the sheath.
7. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 7 cm distal to the catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter.
8. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
9. After catheter is in final position, remove obturator, secure catheter to skin and dress in standard fashion.
10. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Tocino IM, Watanabe A: "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition," *American Journal of Roentgenology*, 1986; 146:487-490.
2. Welch RH, Gravenstein N, Blackshear RH: "Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study," *Journal of Clinical Monitoring*, 1997; 13:75-59.
3. Kelly C, Dumenko L, McGregor SE, McHutchion ME: "Change in Flushing Protocols of Central Venous Catheters," *Oncology Nursing Forum*, 1992; 19:599-605.

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Bloomington, IN 47404 U.S.A.