COOK

MEDICAL

- EN**
3 **Micropuncture®** Introducer Set
Instructions for Use
- DA**
4 **Micropuncture®** indføringsæt
Brugsanvisning
- DE**
5 **Micropuncture®** Einführschleusen-Set
Gebrauchsanweisung
- EL**
6 **ΣΕΤ εισαγωγή Micropuncture®**
Οδηγίες χρήσης
- ES**
7 **Equipo introductor Micropuncture®**
Instrucciones de uso
- FR**
8 **Set d'introducteur Micropuncture®**
Mode d'emploi
- IT**
9 **Set con introduttore Micropuncture®**
Istruzioni per l'uso
- NL**
10 **Micropuncture®** introductieset
Gebruiksaanwijzing
- PT**
11 **Conjunto introdutor Micropuncture®**
Instruções de utilização
- SV**
12 **Micropuncture®** införarsset
Bruksanvisning

T _ M P I S N O C E _ R E V O

MICROPUNCTURE® INTRODUCER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Micropuncture Introducer Set contains a 21 gage introducer needle, a .018 inch wire guide, and an introducer catheter.

INTENDED USE

The Micropuncture Introducer Set is intended for the placement of wire guides up to 0.038 inch diameter into the peripheral vascular system when a small 21 gage needle stick is desired.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- When inserting, manipulating or withdrawing a device through an introducer, always maintain introducer position.
- Do not attempt to insert or withdraw the wire guide and/or introducer if resistance is felt.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through a needle tip may result in breakage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to, the following:

- Bleeding
- Edema
- Extravasation
- Hematoma
- Inflammation, necrosis, or scarring
- Laceration of a vessel, or viscus
- Pain in region
- Vessel perforation
- Infection at access site
- Wire or catheter embolism

INSTRUCTIONS FOR USE

1. Insert the needle into the vessel.
WARNING: Place a finger over the hub of the needle to minimize blood loss and the risk of air aspiration.
2. Advance the .018 inch wire guide through the introducer needle.
CAUTION: Do not withdraw the wire guide through the introducer needle; breakage may result. If the wire guide tip must be withdrawn while the needle is inserted, remove both the needle and the wire as a unit.
3. Withdraw the introducer needle, leaving the wire guide in place.
4. Advance the introducer/dilator pair over the wire guide.
5. Remove the dilator and the wire guide, leaving the introducer catheter in place. **WARNING: Place a finger over the orifice of the introducer catheter to minimize blood loss and the risk of air aspiration.**
6. Advance a wire guide (up to .038 inch diameter) through the introducer catheter.
7. Remove the introducer catheter, leaving the wire guide in place.
8. Proceed with the planned intervention.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, and cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2024
Bloomington, IN 47404 U.S.A.