- EN**
3 **Micropuncture® Pedal Introducer Access Set**
Instructions for Use
- DA**
4 **Micropuncture® indføringsæt til pedaladgang**
Brugsanvisning
- DE**
5 **Micropuncture® Einführschleusen-Set für die A. dorsalis pedis**
Gebrauchsanweisung
- EL**
6 **Σετ πρόσβασης ποδικού εισαγωγέα Micropuncture®**
Οδηγίες χρήσης
- ES**
7 **Equipo introductor para acceso pedio Micropuncture®**
Instrucciones de uso
- FR**
8 **Set d'introducteur pour accès pédieux Micropuncture®**
Mode d'emploi
- IT**
9 **Set con introduttore per accesso pedidio Micropuncture®**
Istruzioni per l'uso
- NL**
10 **Micropuncture® introductieset voor voettoegang**
Gebruiksaanwijzing
- PT**
11 **Conjunto introdutor Micropuncture® para acesso podal**
Instruções de utilização
- SV**
12 **Åtkomstset för Micropuncture® pedalintroducer**
Bruksanvisning

MICROPUNCTURE® PEDAL INTRODUCER ACCESS SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Micropuncture® Pedal Introducer Access Set contains a 21 gage introducer needle, a .018 inch wire guide, an introducer catheter, and a hemostasis valve.

INTENDED USE

This device is intended to introduce up to an 0.038" wire guide into the peripheral vasculature. It is also used to introduce diagnostic and/or interventional devices able to pass through the ID of the introducer.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- When inflating a balloon at, or close to, the introducer tip, ensure the balloon is not inside the distal tip of the introducer.
- When inserting, manipulating or withdrawing a device through an introducer, always maintain introducer position.
- Do not attempt to insert or withdraw the wire guide and/or introducer if resistance is felt.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through a needle tip may result in breakage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to, the following:

- Bleeding
- Edema
- Extravasation
- Hematoma
- Inflammation, necrosis, or scarring
- Laceration of a vessel or viscus
- Pain in region
- Vessel perforation
- Infection at access site
- Wire or catheter embolism

INSTRUCTIONS FOR USE

1. Insert the needle into the vessel.
WARNING: Place a finger over the hub of the needle to minimize blood loss and the risk of air aspiration.
2. Advance the .018 inch wire guide through the introducer needle.
CAUTION: Do not withdraw the wire guide through the introducer needle; breakage may result. If the wire guide tip must be withdrawn while the needle is inserted, remove both the needle and the wire as a unit.
3. Withdraw the introducer needle, leaving the wire guide in place.
4. Advance the pedal introducer/dilator pair over the wire guide.
5. Remove the dilator and the wire guide, leaving the pedal introducer in place. **WARNING: Place a finger over the orifice of the pedal introducer to minimize blood loss and the risk of air aspiration.**
6. Attach the hemostatic valve to the pedal introducer.
7. Proceed with the planned intervention.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, and cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd