**EN**  
**4**

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

**Fig. 1**

**Fig. 2**

# MREYE® BIOPSY AND SPECIAL PURPOSE NEEDLES

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

## DEVICE DESCRIPTION

MReye Needles are manufactured from a non-ferromagnetic material, which prevents occurrence of significant image distortion during Magnetic Resonance Imaging (MRI). Needles are packaged in a protective sleeve, the Centimeasure™ Biopsy Depth Marker, which has graduated centimeter markings beginning at the proximal end.

## INTENDED USE

MReye Needles are intended for initial puncture, positioning or biopsies during MRI. Biopsy needles can be used for cytologic and histologic specimens, depending upon the needle type and diameter.

## CONTRAINDICATIONS

None known

## WARNINGS

MReye Needles provide optimal artifact (distortion) results when a 1.5 Tesla (magnet) rate is used. **Do not exceed a magnet strength of 1.5 Tesla.**

## PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for needle biopsy should be employed.
- Refer to package label for needle specifications.
- Circumferential image distortion may occur around the body of the needle when the needle is oriented perpendicular to the magnetic field.
- RF heating of the needle may occur during MRI procedures. Heating effects are theoretically greatest when the needle is oriented along the length of the bore, near the bore wall.
- If using fast scan protocols, there is a theoretical possibility that peripheral nerve stimulation may occur under worst-case conditions (needle with inserted trocar oriented along the length of the bore, near the bore wall).
- Before using MReye Needles, physicians should have extensive training and experience with interventional Magnetic Resonance Imaging.
- Evaluation of image distortion was performed using a 1.5 Tesla scanner and SE (Spin Echo) and GRASS (Gradient-Recalled Acquisition in the Steady State) MRI sequences, with the device in two different orientations (0° and 90°) and the MRI parameters as shown in the following table. Results of this testing show that under worst-case conditions in terms of artifact size and needle position, circumferential image distortion of up to 2 cm may occur around the body of the needle when the needle is oriented perpendicular to the magnetic field using the GRASS sequence. Under these conditions, error in determining needle tip position may be as much as 1 cm.
- Lung puncture may result in an air embolus, which could lead to ischemia or infarction of major organs including the brain or cardiac system.

### MRI Parameters for Image Distortion Testing

| Parameter             | MRI Pulse Sequence                       |                                                       |
|-----------------------|------------------------------------------|-------------------------------------------------------|
|                       | SE (Spin Echo)                           | GRASS (Gradient-Recalled Acquisition in Steady State) |
| Static Field Strength | 1.5 Tesla                                | 1.5 Tesla                                             |
| Needle Orientation    | 0°, 90°                                  | 0°, 90°                                               |
| TR                    | 300 msec                                 | 100 msec                                              |
| TE                    | 16 msec                                  | 20 msec                                               |
| Flip Angle            | 90°                                      | 30°                                                   |
| Bandwidth             | 16 kHz                                   | 16 kHz                                                |
| Field of View         | Axial, Sagittal: 24x12<br>Coronal: 40x30 | Axial, Sagittal: 24x12<br>Coronal: 40x30              |
| Frequency             | Axial (0°): "x"                          | Axial (0°): "x"                                       |
| Encoding Gradient     | Sagittal: (90°): "z"                     | Sagittal: (90°): "z"                                  |
| Direction             | Coronal: (0°, 90°): "z"                  | Coronal: (0°, 90°): "z"                               |
| Slice Thickness       | 5 mm                                     | 5 mm                                                  |
| Data Acquisition      | 256x256                                  | 256x256                                               |
| Matrix Size           |                                          |                                                       |
| Maximum Gradient      | 1 G/cm                                   | 1 G/cm                                                |
| Field Strength        |                                          |                                                       |
| Phantom Filler T1     | 2 <sup>1</sup> sec                       | 2 <sup>1</sup> sec                                    |
| Phantom Filler T2     | 2 <sup>1</sup> sec                       | 2 <sup>1</sup> sec                                    |

<sup>1</sup> Aqueous solution of 0.45% saline (phantom filler) is assumed to have T1 and T2 value of approximately 2 seconds (similar to water).

## POTENTIAL ADVERSE EVENTS

- Infection
- Internal bleeding and/or bleeding at the puncture site
- Pneumothorax
- Puncture or injury to target organ or to a nearby organ that is traversed by the needle
- Hemoptysis
- Peritonitis
- Damage of tissue adjacent to the needle

## INSTRUCTIONS FOR USE

1. Using MRI guidance, determine the depth of lesion.
2. Adjust needle position using the Centimeasure sleeve and needle depth marker. Final position of needle tip indicates corresponding depth of lesion. **(Fig. 1)**
3. Maintaining needle depth marker position on needle, remove Centimeasure sleeve. **(Fig. 2)**
4. With needle depth marker in position, introduce needle using standard protocol.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.