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McLean-Ring肠饲管套组

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**ZH
25**

McLEAN-RING ENTERAL FEEDING TUBE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The McLean-Ring Enteral Feeding Tube Set consists of:

- 1 TFE-coated, stainless steel Coons interventional wire guide
- 1 Ultrathane® feeding tube
- 1 plastic multipurpose tubing adapter

The feeding tube has a weighted distal end to facilitate passage. The feeding tube endhole allows placement over a wire guide if ineffective peristalsis or obstruction exists.

INTENDED USE

The McLean-Ring Enteral Feeding Tube is intended for enteral feeding.

The product is intended for use by physicians trained and experienced in the placement of feeding tubes. Standard techniques for placement of enteral feeding tubes should be employed.

CONTRAINDICATIONS

Use of the McLean-Ring Enteral Feeding Tube Set is contraindicated in patients with the following conditions:

- Esophageal varices
- Gastric varices
- INR (international normalized ratio) > 1.3 (at time of insertion and/or expected at time of removal)
- Anticoagulated patients (anticoagulated at time of insertion and/or expected to be anticoagulated at time of removal)
- Pathologic coagulopathies
- History of bleeding disorder(s)
- Small or large bowel obstruction(s)
- Ischemic bowel
- Peritonitis
- Esophageal stricture or obstruction
- Complete gastric obstruction
- Recent nasal, oral, esophageal, or gastric surgery or trauma
- Deviated septum
- Inability to pass the feeding tube through the nares
- Uncooperative patient

WARNINGS

None known

PRECAUTIONS

When removing the McLean-Ring Enteral Feeding Tube through the nose or mouth, proceed very slowly.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with placement and use of a McLean-Ring Enteral Feeding Tube include, but are not limited to:

- Bleeding
- Clogged or leaking feeding tube
- Sinusitis
- Premature displacement of the tube
- Aspiration
- Nasal irritation
- Sore throat

ENGLISH

INSTRUCTIONS FOR USE

Adequate Peristalsis and No Gastric Outlet Obstruction

1. If peristalsis is adequate and there is no gastric outlet obstruction, pass the McLean-Ring feeding tube through the external naris in the usual fashion, utilizing a water-soluble surgical lubricant.
2. Advance the feeding tube until there is adequate length within the stomach to allow the feeding tube to pass through the pylorus and duodenum to the level of the ligament of Treitz.

Inadequate Peristalsis or Partial Gastric Outlet Obstruction

1. If gastric outlet obstruction or inadequate peristalsis prevents the feeding tube from advancing on its own, place it over a wire guide.
2. Using either a steerable catheter or a torque control catheter, under fluoroscopy advance the .038 inch (0.97 mm) diameter wire guide through the pylorus and duodenum. The wire guide enclosed with the feeding tube should be used, since this allows adequate length for effective exchange of the feeding tube.
3. Advance the wire guide to well beyond the level of the ligament of Treitz.
4. Leaving wire guide in place, remove catheter.
5. Pre-lubricate the feeding tube prior to introduction.
6. Advance the feeding tube over the wire guide by holding the end of the wire guide steady and sliding the feeding tube along the wire guide as far as it will pass. If resistance is encountered, advance the feeding tube and wire guide together as a unit. Once some progress has been made, withdraw a small amount of wire guide while maintaining slight tension on the end of the feeding tube. This allows the feeding tube to advance slowly through the lumen of the bowel, following the course of the wire guide. With this technique, the ligament of Treitz can be reached in almost all patients.

Pre-existing Gastrostomy or Jejunostomy

1. If the feeding tube is to be placed through a pre-existent gastrostomy or jejunostomy tube, assess the internal diameter of the indwelling tube to assure that the tip of the feeding tube will pass through it.
2. Under fluoroscopy, pass a wire guide through the gastrostomy or jejunostomy tube or stoma to the desired level in the bowel. Advance the feeding tube as described above.

NOTE: It is advisable to vigorously flush the feeding tube with water following the infusion of any of the enteral solutions. Some of these solutions may harden within the feeding tube lumen, causing concretions which may make subsequent infusions difficult.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.