COOK

MEDICAL

EN
3

McLean-Ring Nasojejunal Feeding Tube Set

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Mode d'emploi

T _ M R T _ R E V 6

MCLEAN-RING NASOJEJUNAL FEEDING TUBE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The McLean-Ring Nasojejunal Feeding Tube Set consists of:

- 1 0.038" TFE-coated, stainless steel Coons interventional wire guide
- 1 9.5 French Ultrathane® feeding tube

The feeding tube has a weighted distal end to facilitate passage. The feeding tube endhole allows placement over a wire guide if ineffective peristalsis or obstruction exists.

INTENDED USE

The McLean-Ring Nasojejunal Feeding Tube Set is intended for enteral feeding and medication in adult and pediatric populations. The device is intended to be used in the treatment of the following pediatric subgroups: children and adolescents.

CONTRAINDICATIONS

Use of the McLean-Ring Nasojejunal Feeding Tube Set is contraindicated in patients with the following conditions:

- Esophageal varices
- Gastric varices
- INR (international normalized ratio) > 1.3 (at time of insertion and/or expected at time of removal)
- Anticoagulated patients (anticoagulated at time of insertion and/or expected to be anticoagulated at time of removal)
- Pathologic coagulopathies
- History of bleeding disorder(s)
- Small or large bowel obstruction(s)
- Ischemic bowel
- Peritonitis
- Esophageal stricture or obstruction
- Complete gastric obstruction
- Recent nasal, oral, esophageal, or gastric surgery or trauma
- Deviated septum
- Inability to pass the feeding tube through the nares
- Uncooperative patient
- Not for intravascular use
- Esophageal perforation

WARNINGS

- Improper placement of the feeding tube into the trachea or lungs can lead to serious injury. Always check the final X-ray prior to tube feeding.
- This device was tested to demonstrate a reduced likelihood of a misconnection per Standard AAMI/ ANSI/ISO 80369-1 and AAMI/CN3:2014 Part 3. However, there is a potential to misconnect this device with connectors of other healthcare applications including but not limited to:
 - Intravascular or hypodermic applications (such as female Luer lock fittings)

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in the placement of feeding tubes. Standard techniques for placement of feeding tubes should be employed.
- Do not use a syringe smaller than 30 mL when irrigating tube.
- When removing the McLean-Ring Nasojejunal Feeding Tube through the nose or mouth, proceed very slowly.
- The McLean-Ring Nasojejunal Feeding Tube Set is not intended for use beyond 30 days.
- Use caution with patients who have anomalies or diseases of the nose, throat, or esophagus.
- Any introduction of a wire guide following initial placement, for repositioning or exchange, requires fluoroscopic visualization to prevent perforation.
- Any introduction of a wire guide following initial placement, for repositioning or exchange, requires aseptic technique to prevent contamination.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the McLean-Ring Nasojejunal Feeding Tube Set is MR Unsafe.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with placement and use of a McLean-Ring Nasojejunal Feeding Tube include, but are not limited to:

- Bleeding
- Clogged or leaking feeding tube
- Sinusitis
- Premature displacement of the tube
- Aspiration
- Nasal irritation
- Sore throat
- Manual removal of the feeding tube by the patient
- Bowel erosion and/or intestinal perforation
- Pneumothorax from inadvertent pleural placement
- Pulmonary injury
- Improper placement of the feeding tube into the lung or brain
- Coughing or dyspnea, indicating potential bronchial placement
- Feeding tube coiled in esophagus or posterior pharynx
- Esophageal tear from trauma or tube passing
- Feeding tube dislodgements during therapy, necessitating removal and new feeding tube placement
- Aspiration of stomach contents despite appropriate placement

INSTRUCTIONS FOR USE

Adequate Peristalsis and No Gastric Outlet Obstruction

1. If peristalsis is adequate and there is no gastric outlet obstruction, pass the feeding tube through the external naris in the usual fashion, utilizing a water-soluble surgical lubricant.
2. Advance the feeding tube until there is adequate length within the stomach to allow the feeding tube to pass through the pylorus and duodenum to the level of the ligament of Treitz.

Inadequate Peristalsis or Partial Gastric Outlet Obstruction

1. If gastric outlet obstruction or inadequate peristalsis prevents the feeding tube from advancing on its own, place it over a wire guide.
2. Using either a steerable catheter or a torque control catheter, under fluoroscopy advance the 0.038 inch (0.97 mm) diameter wire guide through the pylorus and duodenum. The wire guide enclosed with the feeding tube should be used, since this allows adequate length for effective exchange of the feeding tube.
3. Advance the wire guide to well beyond the level of the ligament of Treitz.
4. Leaving wire guide in place, remove catheter.
5. Pre-lubricate the feeding tube prior to introduction.
6. Advance the feeding tube over the wire guide by holding the end of the wire guide steady and sliding the feeding tube along the wire guide as far as it will pass. If resistance is encountered, advance the feeding tube and wire guide together as a unit. Once some progress has been made, withdraw a small amount of wire guide while maintaining slight tension on the end of the feeding tube. This allows the feeding tube to advance slowly through the lumen of the bowel, following the course of the wire guide. With this technique, the ligament of Treitz can be reached in almost all patients.

NOTE: The McLean-Ring Nasojejunal Feeding Tube is intended to connect to a catheter tip syringe or oral syringe. If enteral feeding is needed, use a catheter tip syringe. If medications are to be administered, use an oral syringe. Do not use a Luer lock or Luer tip syringe.

Post-Placement Care

1. Irrigate tube after each use. Do not use a syringe smaller than 30 mL when irrigating the tube. In order to prevent feeding tube clogging, always flush feeding lumen with prescribed amount of water after feeding administration, every 4-6 hours during continuous feeding, before administering medication and after medication administration is complete.
NOTE: The amount of water used for flushing will depend on the needs of the patient and clinical condition.
2. Assess the feeding tube daily for damage, clogging, or abnormal discoloration.
3. Device is intended for use for up to 30 days.

Feeding Tube Removal

1. Irrigate the feeding tube prior to removal.
2. If the feeding tube is not occluded, insert a wire guide down the lumen and remove or exchange the feeding tube as needed. If occluded, close the cap of the feeding tube and/or clamp the feeding tube distal to the hub during tube removal to minimize the risk of aspirating fluids into the trachea during withdrawal.
3. Slowly withdraw the feeding tube and discard per standard procedures.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

For Customer support, please contact Cook Medical at:

Cook Medical Inc.
P.O. BOX 4195
Bloomington, IN 47402
U.S.A.
Phone: (812) 339-2235
Alternative Phone: (800) 475-4500 (toll free)
Fax: (800) 554-8335 (toll free)
Email: customersupport@cookmedical.com

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Bloomington, IN 47404 U.S.A.