



COOK**MEDICAL****CE**
0123

EN 4	Multipurpose Drainage Catheter Instructions for Use
BG 7	Универсален дренажен катетър Инструкции за употреба
CS 12	Víceúčelový drenážní katetr Návod k použití
DA 16	Multifunktionsdrænagekateter Brugsanvisning
DE 20	Mehrzweck-Drainagekatheter Gebrauchsanweisung
EL 24	Καθετήρας παροχέτευσης πολλαπλών εφαρμογών Οδηγίες χρήσης
ES 29	Catéter de drenaje multipropósito Instrucciones de uso
ET 33	Mitmeotstarbeline drenaazikateeter Kasutusjuhend
FI 37	Monikäyttöinen dreneerauskatetri Käyttöohjeet
FR 41	Cathéter de drainage multi-usages Mode d'emploi
HR 45	Višenamjenski drenažni kateter Upute za uporabu
HU 49	Multipurpose drenázskatéter Használati utasítás
IT 53	Catetere di drenaggio multiuso Istruzioni per l'uso
LT 58	Daugiafunkcis drenažo kateteris Naudojimo instrukcija
LV 62	Daudzfunkcionāls drenāžas katetrs Lietošanas instrukcija
NL 65	Drainagekatheter voor meerdere doeleinden Gebruiksaanwijzing
NO 70	Generelt drenasjekateter Bruksanvisning
PL 74	Uniwersalny cewnik do drenażu Instrukcja użycia
PT 78	Cateter de drenagem multiusos Instruções de utilização
RO 82	Cateter de drenaj multifuncțional Instrucțiuni de utilizare
SK 86	Viacúčelový drenážny katéter Návod na použitie
SL 90	Večnamenski drenažni kateter Navodila za uporabo
SV 94	Universell dränagekateter Bruksanvisning
TR 98	Çok Amaçlı Drenaj Kateteri Kullanma Talimatı

T _ M U L T I 2 - M _ R E V 1

MULTIPURPOSE DRAINAGE CATHETER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Multipurpose Drainage Catheter Ultrathane®, Multipurpose Drainage Catheter Set Ultrathane® Neff Access Set with Nitinol Wire Guide, Dawson-Mueller Drainage Catheter Ultrathane®, and Multipurpose Small Pigtail Drainage Catheter - Ultrathane® are manufactured in a range of French sizes, lengths and sideport quantities. Locking Loop Catheters are available with one of two types of locking mechanisms:

- Mac-Loc® Locking Loop
- Cook-Cope Type Loop

When ordered with the suffix –HC, the distal part of the catheter, sideports included, has a hydrophilic coating.

Performance Characteristics

These devices are manufactured in 5, 6, 6.3, 7, 8.5, 10.2, 12, 14, or 16 French; 15, 25, 35, 40, 45, 50, or 60 cm length; and 0.018, 0.035, or 0.038 inch end hole diameter, which allows the catheter to accept wire guides 0.038 inches or smaller. Sideports allow fluid drainage facilitated by gravity. The tapered distal tip of the catheter facilitates over-the-wire trackability and is designed with a pigtail shape in order to maintain placement in the body.

Device Compatibility

These devices are designed to be compatible with components such as stiffeners, fixation device, connecting tubes and drainage bags with standard Luer Lock fittings, as well as standard size wire guides (0.038 inches or smaller).

Qualitative and Quantitative Information

Device Material		Weight (grams)
Polyurethane	Shaft Tubing	Up to 7.082
Hydrogel	Hydrophilic Coating (-HC and -RH devices only)	Up to 0.012
Silicone Lubricant	Silicone Coating (not included on -HC and -RH devices)	Up to 0.021
Mac-Loc® Locking Loop Mechanism		
Nylon	Suture	Up to 0.089
Cope-Loop Type Locking Loop Mechanism*		
PTFE-Coated Polyester	Suture	Up to 0.043

*The device includes a sleeve made of natural rubber latex. The sleeve is a skin-contacting component. **NATURAL RUBBER LATEX** may cause allergic reactions.

Expected Lifetime of the Device

The expected lifetime of the device is 3 years shelf-life and 90 days indwell time.

Patient Population

The target population of the device is adult patients who require placement of a catheter to allow for percutaneous drainage, particularly those requiring drainage related to nephrostomy or cholecystostomy and abscess or biliary drainage.

Intended User

These devices are intended for use by physicians or other properly licensed providers trained and experienced in diagnostic and interventional techniques. Drainage procedures that require the use of these devices should only be performed by licensed and trained clinical practitioners, such as physicians, and licensed and trained healthcare providers. Maintenance and care of the device may be performed by laypersons (e.g. patients or lay caregivers), nurses, and technicians.

Contact with Body Tissue

These devices are externally communicating devices contacting tissue. They facilitate percutaneous drainage from tissues, organs, or anatomic spaces to outside the body.

Operating Principle

These devices are placed with imaging guidance to identify and facilitate control of the distal tip of the catheter during placement. Drainage occurs through the sideports and is facilitated by gravity. Deploying the locking mechanism maintains placement of the catheter in the body through retention of the pigtail loop.

INTENDED USE

These devices, which consist of drainage catheters and components for access, insertion, and/or securement, are intended to gain percutaneous access, either by direct stick or Seldinger access technique, and to serve as a conduit for long-term fluid drainage from an organ or cavity to outside the body in a variety of applications (e.g., nephrostomy or cholecystostomy procedures and abscess or biliary drainage).

INDICATIONS FOR USE

These devices are indicated for percutaneous drainage in adult patients for treatment of:

- Biliary tract obstruction, infection, stones, leak, and fistula
- Urinary tract obstruction, infection, stones, leak, and fistula
- Abscess/fluid collections (e.g., pelvic, peritoneal, pleural, subcutaneous)

CLINICAL BENEFITS

- Percutaneous drainage of fluid to outside the body
- Percutaneous drainage in a minimally invasive manner (i.e., avoiding surgical drainage)

CONTRAINDICATIONS

- Pericardial drainage
- Uncorrectable coagulopathy (in the context of percutaneous nephrostomy and percutaneous transhepatic biliary drainage)

WARNINGS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Minimize catheter manipulation during application and removal of the fixation device.
- See product label for latex indication. **NATURAL RUBBER LATEX** may cause allergic reactions.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Do not use fixation device where loss of adherence could occur, such as with a confused patient or non-adherent skin, or if a patient has known tape or adhesive allergies.

PRECAUTIONS

- Standard techniques for placement of percutaneous drainage catheters should be employed. Caution should be used if a safe, appropriate, and acceptable pathway to the target site is not present.
- Manipulation of products requires ultrasound, fluoroscopy, or other imaging guidance.
- When inserting a stiffening cannula into a catheter with retention suture, hold suture during cannula insertion to avoid bunching or tangling of suture.
- A TFE-coated wire guide must be used with Ultrathane® catheters.
- Activate the hydrophilic coating, if present, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.
- Catheters should be irrigated on a routine basis to ensure function.
- Patients with indwelling drainage catheters should be evaluated routinely to ensure continuous function of the catheter and adherence of the fixation device, if applicable.
- Traction on the locking suture, if present, should be sufficient to ensure adequate retention of the tip, but should not be overly tight. Verify catheter tip configuration with imaging.
- It is recommended to use a wire guide when removing a Locking Loop catheter.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- Pregnant women have physiologic and anatomic alterations as well as additional procedural risks (e.g., radiation exposure to the fetus, effects of anesthesia on the fetus) that should be considered when using these devices in pregnant patients.
- Avoid fixation device contact with alcohol or acetone. Both can weaken bonding of components and the fixation device anchor pad adherence to the skin.

MR SAFETY INFORMATION

Nonclinical testing has demonstrated that this device is **MR Conditional**. A patient with this device may be safely scanned after placement under the following conditions:

- Static magnetic field of 3.0 T and 1.5 T only
- Maximum magnetic field spatial gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

For Mac-Loc Locking Loop Mechanism:

Under the scan conditions provided above, this device is expected to produce a maximum temperature rise of less than 2.1°C after 15 minutes of continuous scanning in nonclinical testing. The image artifact extends approximately 40.1 mm from the device as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

For Cook-Cope Type Locking Loop Mechanism:

Under the scan conditions provided above, this device is expected to produce a maximum temperature rise of less than 2.1°C after 15 minutes of continuous scanning in nonclinical testing.

The image artifact extends approximately 40.1 mm from the device as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

POTENTIAL ADVERSE EVENTS

- Air embolism
- Arrhythmia, transient
- Biliary events
- Bleeding and/or hemorrhage
- Catheter-related events (damage, dislodgement, fracture, kinking, leakage, malposition, migration, obstruction, occlusion)
- Death
- Fever
- Fistula formation
- Fluid/electrolyte imbalance
- Hypotension and/or shock
- Infectious and/or inflammatory events
- Injury and/or perforation of organs and/or structures (tissues, viscera, vessels, nerves)
- Pain
- Sepsis and/or septic shock
- Tumor seeding of catheter tract
- User injury
- Vasovagal reaction

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Catheter Placement

1. Under imaging guidance, perform standard techniques for placement of percutaneous drainage catheters, either by Seldinger access or trocar access.
2. Once catheter is in desired location, remove any wire guides, trocars, or stiffeners, allowing the catheter to form its configuration.
3. For Locking Loop catheters, lock the catheter in place using appropriate technique for the locking mechanism type, as described below.

For Mac-Loc® Locking Loop Mechanism:

- a. Stabilize the Mac-Loc catheter hub assembly with one hand and pull back on the drawstring to form the distal catheter loop configuration. **(1)**
- b. While maintaining traction on the drawstring, push the locking cam lever down until a distinct “snap” is felt. The distal loop of the catheter is now locked into position. **(2)**
- c. Trim off the excess drawstring. **(3)**

For Cook-Cope Type Locking Loop Mechanism:

- a. Pull suture end tight to form loop configuration in catheter, and tie securely.
- b. Trim off excess suture, and slide latex sleeve over suture to prevent leakage.

Unlocking Catheter Loop

For Mac-Loc® Locking Loop Mechanism:

- a. While stabilizing the Mac-Loc catheter hub assembly with one hand, position a small, blunt object (approximately the shape and size of a ball point pen or small forceps) into the Mac-Loc release notch.
- b. Pry upward until the locking cam lever is free. (4)
NOTE: For catheter exchange, advance the distal end of a wire guide into the locked loop configuration of the catheter before unlocking the Mac-Loc assembly. Release the Mac-Loc as described above. Advance the wire guide through the catheter end hole. Catheter exchange can now be performed.

For Cook-Cope Type Locking Loop Mechanism:

- a. Advance wire guide into catheter.
- b. Uncover suture by folding back the latex sleeve, and release suture.
- c. Withdraw catheter over wire guide.
NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

Applying the Catheter Fixation Device

- 1. Bring the fixation device anchor pad into position. Stabilize catheter between fingers. Place catheter into fixation device retainer.
- 2. Close lid, peel off paper backing, and place on skin.
- 3. Suture the fixation device to the skin if desired or deemed necessary.

Removing the Catheter Fixation Device

- 1. To open retainer, press tab gently with thumb and lift lid.
- 2. Remove the fixation device anchor pad adhesive with alcohol.
NOTE: The fixation device should be replaced at least every 7 days.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use. Please inform the patient of the **Multipurpose Drainage Catheter - Patient Implant Information** booklet. You can use the basic unique device identifier for this device seen in **Table 1** below to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: ec.europa.eu/tools/eudamed.

Table 1: Basic Unique Device Identification

Part Number Prefix	Part Number Suffix	Basic Unique Device Identifier (BUDI)
ULT...	...MDC	0827002CINC050-321_A3YZ
ULT...	...MCL	0827002CINC050-321_A2YX
ULT...	...CLDM or ...CLDMS	0827002CINC050-321_A-1K3
ULT...	...CLM or ...CLMS	0827002CINC050-321_A-1K3
CLM...	...NPAS-NT	0827002CINC050-321_A-1K3

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2024
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.