COOK®

MEDICAL

CE

0123

EN	Multipurpose Drainage Catheter
5	Instructions for Use
CS	Víceúčelový drenážní katetr
7	Návod k použití
DA	Multifunktionsdrænagekateter
9	Brugsanvisning
DE	Mehrzweck-Drainagekatheter
11	Gebrauchsanweisung
EL	Καθετήρας παροχέτευσης πολλαπλών εφαρμογών
13	Οδηγίες χρήσης
ES	Catéter de drenaje multipropósito
16	Instrucciones de uso
FR	Cathéter de drainage multi-usages
18	Mode d'emploi
HU	Multipurpose drenázskatéter
20	Használati utasítás
IT	Catetere di drenaggio multiuso
23	Istruzioni per l'uso
NL	Drainagekatheter voor meerdere doeleinden
25	Gebruiksaanwijzing
NO	Generelt drenasjekateter
27	Bruksanvisning
PL	Uniwersalny cewnik do drenażu
29	Instrukcja użycia
PT	Cateter de drenagem multiusos
32	Instruções de utilização
SV	Universell dränagekateter
34	Bruksanvisning

T _ M U L T I 2 _ R E V 1

Fig. 3


Fig. 2



Fig. 1


Fig. 4

MULTIPURPOSE DRAINAGE CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Multipurpose Drainage Catheters are constructed from Ultrathane® or polyethylene and come in a range of French sizes, lengths and sideport quantities. Locking Loop Catheters are available with one of two types of locking mechanisms:

- Mac-Loc® Locking Loop
- Cook-Cope Type Loop

When ordered with the suffix -HC, the distal part of the catheter, sideports included, has an AQ® hydrophilic coating.

INTENDED USE

Multipurpose drainage catheters are intended for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary and abscess), either by direct stick or Seldinger access technique.

CONTRAINdications

- Bleeding diathesis and uncontrolled hypertension
- Anticoagulant use

WARNINGS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- See product label for latex indication. **NATURAL RUBBER LATEX** may cause allergic reactions.

PRECAUTIONS

- These products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous drainage catheters should be employed.
- Manipulation of products requires ultrasound, fluoroscopy, or other imaging guidance.
- When inserting a stiffening cannula into a catheter with retention suture, hold suture during cannula insertion to avoid bunching or tangling of suture.
- A TFE-coated wire guide must be used with Ultrathane® catheters.
- Activate the hydrophilic coating, if present, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.
- Catheters should be irrigated on a routine basis to ensure function.
- Patients with indwelling drainage catheters should be evaluated routinely to ensure continuous function of the catheter.
- Traction on the locking suture, if present, should be sufficient to ensure adequate retention of the tip, but should not be overly tight. Verify catheter tip configuration by fluoroscopy.
- It is recommended to use a wire guide when removing a Locking Loop catheter.
- The Peel-Away® Pigtail Straightener, if present, is not to be used as a vascular introducer sheath.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Arrhythmia (transient), including tachycardia
- Bleeding (including chronic hematuria, post-procedural, intraprocedural, puncture site, hemoperitoneum)
- Bile leak/biloma
- Catheter dislocation/dislodgement
- Catheter malposition/transcolonic catheter placement
- Catheter occlusion/obstruction
- Damage to adjacent organs
- Death (due to uncontrolled hypotension)
- Hypoxia
- Inflammation/infection (including cellulitis at puncture site, fever/chills within 24 hours, rigors, sepsis during catheter exchange, septic shock leading to death, urinary tract infection, peritonitis)
- Mechanical catheter damage
- Pain (including intraprocedural pain, post-procedural pain)
- Perihepatic fluid collection
- Small bowel injury
- Urinary leak

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Ultrathane percutaneous drainage catheter is **MR Conditional**. A patient with this device may be safely scanned after placement under the following conditions:

- Static magnetic field of 3.0 T and 1.5 T only
- Maximum magnetic field spatial gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

For Mac-Loc Locking Loop Mechanism:

Under the scan conditions provided above, the Ultrathane percutaneous drainage catheter is expected to produce a maximum temperature rise of less than 2.1°C after 15 minutes of continuous scanning in nonclinical testing.

The image artifact extends approximately 40.1 mm from the Ultrathane percutaneous drainage catheter as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

For Cook-Cope Type Locking Loop Mechanism:

Under the scan conditions provided above, the Multipurpose Drainage Catheter is expected to produce a maximum temperature rise of less than 2.1°C after 15 minutes of continuous scanning in nonclinical testing.

The image artifact extends approximately 40.1 mm from the Multipurpose Drainage Catheter as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)

Fax: 209-669-2450

Web: www.medicalert.org

INSTRUCTIONS FOR USE

Catheter Placement

1. Under fluoroscopic control, perform standard techniques for placement of percutaneous drainage catheters, either by Seldinger access or trocar access.
2. Once catheter is in desired location, remove any wire guides, trocars, or stiffeners, allowing the catheter to form its configuration.

3. For Locking Loop catheters, lock the catheter in place using appropriate technique for the locking mechanism type, as described below.

For Mac-Loc® Locking Loop Mechanism:

- a. Stabilize the Mac-Loc catheter hub assembly with one hand and pull back on the drawstring to form the distal catheter loop configuration. (**Fig. 1**)
- b. While maintaining traction on the drawstring, push the locking cam lever down until a distinct "snap" is felt. The distal loop of the catheter is now locked into position. (**Fig. 2**)
- c. Trim off the excess drawstring. (**Fig. 3**)

For Cook-Cope Type Locking Loop Mechanism:

- a. Pull suture end tight to form loop configuration in catheter, and tie securely.
- b. Trim off excess suture, and slide latex sleeve over suture to prevent leakage.

Unlocking Catheter Loop

For Mac-Loc® Locking Loop Mechanism:

- a. While stabilizing the Mac-Loc catheter hub assembly with one hand, position a small, blunt object (approximately the shape and size of a ball point pen or small forceps) into the Mac-Loc release notch.
- b. Pry upward until the locking cam lever is free. (**Fig. 4**)

NOTE: For catheter exchange, advance the distal end of a wire guide into the locked loop configuration of the catheter *before unlocking the Mac-Loc assembly*. Release the Mac-Loc as described above. Advance the wire guide through the catheter end hole. Catheter exchange can now be performed.

For Cook-Cope Type Locking Loop Mechanism:

- a. Advance wire guide into catheter.
- b. Uncover suture by folding back the latex sleeve, and release suture.
- c. Withdraw catheter over wire guide.

NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd