

COOK**MEDICAL****CE**
0123

EN 10	Retracta® Detachable Embolization Coils Instructions for Use
BG 14	Отделящи се спирали за емболизация Retracta® Инструкции за употреба
CS 19	Odpoutatelné embolizační spirálky Retracta® Návod k použití
DA 23	Retracta® adskillelige emboliseringscoils Brugsanvisning
DE 27	Retracta® ablösbare Embolisationsspiralen Gebrauchsanweisung
EL 32	Αποσπώμενα σπειράματα εμβολισμού Retracta® Οδηγίες χρήσης
ES 37	Espiraless de embolización largables Retracta® Instrucciones de uso
FI 41	Irrotettavat Retracta®-embolisaatiokoilit Käyttöohjeet
FR 45	Coils détachables Retracta® Mode d'emploi
HU 50	Retracta® leválasztható embolizációs spirálok Használati utasítás
IT 54	Spirali per embolizzazione a rilascio controllato Retracta® Istruzioni per l'uso
LT 59	„Retracta®“ atskiriamosios embolizavimo spiralės Naudojimo instrukcija
LV 63	Retracta® atvienojamās embolizācijas spirāles Lietošanas instrukcija
NL 67	Retracta® loskoppelbare embolisatiecoils Gebruiksaanwijzing
NO 72	Retracta® utløsbare emboliseringsspiraler Bruksanvisning
PL 76	Odłączalne spirale embolizacyjne Retracta® Instrukcja użycia
PT 80	Espiraiss de embolização amovíveis Retracta® Instruções de utilização
RO 85	Bobine de embolizare detaşabile Retracta® Instrucţiuni de utilizare
SK 89	Odpojiteľné embolizačné cievky Retracta® Návod na použitie
SL 93	Snemljive embolizacijske tuljave Retracta® Navodila za uporabo
SV 97	Retracta® löstagbara emboliseringsspiraler Bruksanvisning
TR 101	Retracta® Ayrılabilir Embolizasyon Sarmalları Kullanma Talimatı


T _ M W C E R - M _ R E V O

- 1 Straight coil length
- 2 Catheter inner diameter
- 1 Права дължина на спиралата
- 2 Вътрешен диаметър на катетъра
- 1 Délka rovné spirálky
- 2 Vnitřní průměr katetru
- 1 Lige coillængde
- 2 Katetrets indvendige diameter
- 1 Länge der geraden Spirale
- 2 Innendurchmesser des Katheters
- 1 Μήκος σπειράματος σε ευθεία θέση
- 2 Εσωτερική διάμετρος καθετήρα
- 1 Longitud de la espiral recta
- 2 Diámetro interior del catéter
- 1 Koilin pituus suorana
- 2 Katetrin sisäläpimitta
- 1 Longueur de coil droite
- 2 Diamètre interne du cathéter
- 1 Egyenes spirálhossz
- 2 Katéter belső átmérője
- 1 Lunghezza della spirale diritta
- 2 Diametro interno del catetere

- 1 Tiesios spirālės ilgis
- 2 Kateterio vidinis skersmuo
- 1 Taisnas spirālės garums
- 2 Katetra iekšējais diametrs
- 1 Lengte rechte coil
- 2 Binnendiameter katheter
- 1 Rett spirallengde
- 2 Kateterets indre diameter
- 1 Długość prostej spirali
- 2 Wewnętrzna średnica cewnika
- 1 Comprimento da espiral reta
- 2 Diâmetro interno do cateter
- 1 Lungimea bobinei drepte
- 2 Diametrul interior al cateterului
- 1 Dĺžka rovnej cievky
- 2 Vnútorný priemer katétra
- 1 Dolžina poravnane tuljave
- 2 Notranji premer katetra
- 1 Spiralens längd i utsträckt tillstånd
- 2 Kateterns innerdiameter
- 1 Düz sarmal uzunluğu
- 2 Kateter iç çapı

3

- 1 Safety Clip
- 2 Loading Cartridge
- 3 Torque Device
- 4 Luer Lock Adapter
- 1 Предпазен клипс
- 2 Касета за зареждане
- 3 Усукващо устройство
- 4 Луеров заключващ адаптер
- 1 Bezpečnostní svorka
- 2 Nasazovací kazeta
- 3 Otáčecí zařízení
- 4 Adaptér Luer lock
- 1 Sikkerhedsclips
- 2 Ladehylster
- 3 Momenthåndtag
- 4 Luer lock-adapter
- 1 Sicherheitsclip
- 2 Ladekartusche
- 3 Torquer
- 4 Luer-Lock-Adapter
- 1 Κλιπ ασφάλισης
- 2 Φυσίγγιο τοποθέτησης
- 3 Συσκευή ροπής στρέψης
- 4 Προσαρμογέας ασφάλισης Luer
- 1 Clip de seguridad
- 2 Cartucho de carga
- 3 Dispositivo de torque
- 4 Adaptador Luer Lock

- 1 Turvapidike
- 2 Latauskasetti
- 3 Vääntölaite
- 4 Luer-lukkosovitin
- 1 Clip de securitate
- 2 Cartouche de chargement
- 3 Torqueur
- 4 Adaptateur Luer Lock
- 1 Biztosítókapocs
- 2 Betöltőbetét
- 3 Forgatóeszköz
- 4 Luer-záras adapter
- 1 Clip di sicurezza
- 2 Dispositivo di caricamento
- 3 Dispositivo di torsione
- 4 Adattatore Luer Lock
- 1 Apsauginis spaustukas
- 2 Užtaisymo kasetė
- 3 Sukimo įtaisas
- 4 Luerio fiksatoriaus adapteris
- 1 Drošības skava
- 2 Spirāles ielādēšanas kasetne
- 3 Griezes momenta ierīce
- 4 Luera savienojuma adapters
- 1 Veiligheidsclip
- 2 Laadhuls
- 3 Torsie-instrument
- 4 Luerlock-adapter

- 1 Sikkerhetsklemme
- 2 Ladepatron
- 3 Vridningsenhet
- 4 Luer-lock-adapter

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RETRACTA® DETACHABLE EMBOLIZATION COILS

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Retracta® Detachable Embolization Coil is composed of two main components: a fibered platinum embolization coil and a delivery wire. It is designed to be deployed into the target location by the delivery wire which is connected to the embolization coil via a threaded mechanism on the proximal end of the coil. The coil interlocks with a corresponding mechanism on the distal end of the delivery wire allowing the wire to be fully repositioned prior to being permanently deployed by detaching from the coil when the wire is turned 8 to 10 counterclockwise rotations.

Performance Characteristics

- Coil design and soft platinum material allow device to form a tight occluding mass.
- Delivery system allows the coil to be fully advanced and completely repositioned before the coil is finally deployed.
- Delivery system provides safe delivery of embolization coils when correct positioning is especially critical.
- Radiopaque platinum material allows for visibility under fluoroscopy.
- Synthetic fibers spaced along coil increase thrombogenicity.
- Wide range of diameters: 4-20 mm (even numbers); Lengths: 7, 14 cm.
- MR conditional at 3T or less.
- Coil design maintains radial force to prevent migration when oversized by 10-20% for arteries or 30-50% for veins.
- Obstructs or reduces blood flow in vasculature.
- Can be used alone or in combination with other platinum or Inconel alloy pushable and detachable embolization coils.
- Multiple coils may be used to achieve complete occlusion.
- Length is identified on product label by **Fig. 1**.
- The lifetime of the implant is intended to be permanent. Testing has been completed to assure that the device will perform as intended for the lifetime of the patient or duration of device implantation.

Device Compatibility

- 0.035" Retracta Coils are compatible with angiographic catheters with an 0.035" inner diameter.
- Catheter inner diameter is identified on product labels by **Fig. 2**.

Qualitative and Quantitative Information

Device Materials		Weight (grams)
Platinum-tungsten alloy	Wire for coiling	Up to 0.568
Nylon	Fiber	Up to 0.001

Expected Lifetime of the Device

The expected lifetime of Retracta Detachable Embolization Coils is permanent, with an indwell duration spanning the life of the patient.

Patient Population

The target population for Retracta Detachable Embolization Coils is adult patients who require vascular embolization procedures that include cessation or reduction of blood flow in peripheral arterial and venous vessels.

Intended User

This device is intended for use by physicians who are experienced in interventional techniques. Embolization procedures that require the use of the Embolization Devices should only be performed by licensed clinical practitioners.

Contact with Body Tissue

The Retracta Embolization Coils are for single use, categorized as implant devices in permanent contact with circulating blood and peripheral vasculature.

Operating Principle

The Retracta Detachable Embolization Coil is a platinum coil with spaced synthetic fibers, and is attached to a delivery system. The Retracta Detachable Embolization Coil system allows the coil to be advanced and fully repositioned prior to deployment. Once in the vessel, coil design maintains radial force to prevent migration. The coil and its synthetic fibers induce thrombogenicity

to form a tight occluding mass to obstruct or reduce blood flow in the vasculature. Multiple coils may be used to achieve complete occlusion.

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross-sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion. The combination of the coaxial technique and either the anchor or scaffold technique significantly enhances stability of coil deployment.

- **Coaxial technique:** The use of an outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (**Fig. 9**)
- **Anchor technique:** The anchor technique provides safe and distal occlusion when there is a question about instability of coils. At least 2 cm of a coil is advanced into the side branch, which is normally sacrificed. The rest of the coil is then deployed just proximal to that side branch, and additional coils are packed. (**Fig. 10**)
- **Scaffold technique:** The scaffold technique is used for high-flow vessels when there is concern about migration of a softer coil. A high radial force coil is placed initially. Then, several high radial force coils or soft coils may be packed within the scaffold. (**Fig. 11**)

In general, the first coil selected should have a diameter that is 20% larger, or at least 2 mm oversized, than the vessel that is being occluded.

INTENDED USE

This permanent embolic device is a fibered, metallic coil that is intended to cause cessation or reduction of blood flow in peripheral arteries and veins by means of physical occlusion and thrombosis induced by the coil and its synthetic fibers.

INDICATIONS FOR USE

The Retractable Detachable Embolization Coil is indicated for arterial and venous embolization in the peripheral vasculature. They may be used for treatment of vascular abnormalities, treatment of hemorrhage, devascularization of benign and malignant tumors, devascularization of benign tissue, flow redistribution, and endoleak management.

CLINICAL BENEFITS

Successful adjunctive (e.g., pre-surgical), palliative, or curative

- treatment of vascular abnormalities,
- treatment of hemorrhage,
- devascularization of benign or malignant tumors,
- devascularization of benign tissue,
- flow redistribution, or
- endoleak management

by complete or partial vessel occlusion via a minimally invasive procedure.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that single and multiple Retractable Detachable Embolization Coils are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions:

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1700 gauss/cm (17 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Retractable Detachable Embolization Coils are not expected to result in a temperature rise of more than 3.0 °C after 15 minutes of continuous scanning.

The image artifact extends approximately 5.7 mm from the Retractable Detachable Embolization Coil as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

Absolute contraindications to percutaneous transcatheter embolization:

General:

- Active systemic infection
- Clinically manifest hyperthyroidism

In the context of chemoembolization (TACE) or radioembolization (TARE) for hepatic malignancy:

- Active hepatic infection
- Acute variceal bleeding in previous 2 weeks
- Advanced cirrhosis
- Arteriovenous fistulas, untreatable
- Bilirubin >5 mg/dL
- Capecitabine treatment within 2 months before TARE treatment with resin spheres
- Child-Pugh C liver disease or Child-Pugh ≥ 8
- Contraindications for arteriography and/or transarterial interventions
- Extrahepatic metastases
- Inability to catheterize hepatic artery
- Insufficient functional liver reserve (for TARE)
- Large tumor burden (e.g., bilobar tumor affecting a large portion of the liver, >50% liver replaced by tumor)
- Liver failure, fulminant
- Lung shunting, severe (TARE)
- Nontarget embolization of the gastrointestinal tract, anticipated (TARE)
- Portal vein obstruction (TACE)
- Reversed portal flow

WARNINGS

- Neurological deficit, ischemic stroke, spinal infarct, or ischemic infarct can occur from occlusion of normal vessels by emboli.
- The devices are made of a 92% platinum, 8% tungsten alloy; the nitinol delivery wire is coated with Inconel, a nickel-chromium alloy. The devices should not be implanted in patients who have allergies to any of the coil or delivery components.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- The Retractable Detachable Embolization Coil is not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Warning: The user should consider potential risks associated with radiation exposure.

PRECAUTIONS

- This product is intended for use by healthcare providers experienced in arterial and venous vessel embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Access-related and catheter-angiographic side effects (access-site nerve injury, contrast agent-induced nephropathy, contrast reaction, embolism [e.g., air, atheromatous debris, thrombus], hematoma, infection, nausea/vomiting, vessel occlusion, pain, pseudoaneurysm/arteriovenous fistula, subintimal passage, thrombosis, vasospasm, vasovagal syncope, vessel dissection, vessel injury, vessel perforation)
- Allergic reaction to coil components
- Bile leak/bile duct injury/biloma
- Bleeding/hemorrhage
- Coil migration
- Embolization target ischemia/infarction/inflammation (i.e., the clinical effects, whether desired or not, that result from vessel occlusion in the intended target area)
- Fistula
- Infection
- Liver insufficiency/failure
- Nontarget embolization causing ischemia or infarction
- Pleural effusion, ascites

- Postembolization syndrome (systemic signs and symptoms in addition to anatomy-specific signs, symptoms and abnormal laboratory values of postembolization syndrome that may affect particular patient populations)
- Procedure-related mortality
- Renal dysfunction (e.g., increased renal laboratory values, renal insufficiency/failure)
- Thrombosis

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

This device is accompanied by an implant card, that should be given to the patient after it has been completed by the Healthcare Professional.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Perform an angiogram and measure the diameter of the vessel to be occluded.
2. Firmly grasping the loading cartridge between thumb and forefinger, introduce the metal end of the loading cartridge into the base of the catheter hub. Lock the loading cartridge onto the catheter hub by turning the Luer lock adapter clockwise. **(Fig. 5)**
3. Unlock the white safety clip and advance at least 30 cm of the delivery wire into the catheter so that the coil is fully loaded in the catheter. **(Fig. 6)**
4. Unlock the loading cartridge from the catheter hub by turning the Luer lock adapter counter-clockwise.
5. Remove the delivery wire holder and the metal loading cartridge from the catheter hub while holding the delivery wire stationary. Remove the torque device and reserve it for use later in the procedure.
6. Under fluoroscopic visualization, slowly advance the delivery wire until the entire length of the coil exits the distal end of the catheter. Ensure that the junction remains positioned just inside the catheter tip. **(Fig. 7)**
NOTE: Advancing the delivery wire slowly allows the junction to be seen more easily and reduces the risk of damaging it.
NOTE: If significant resistance is encountered during coil advancement, do not continue advancing. Retract the delivery wire slightly, then gently re-advance it. If there is still significant resistance, withdraw the delivery wire from the catheter and try using a new coil with a shorter length.
NOTE: Do not turn the delivery wire counterclockwise during advancement; the coil may be unintentionally detached.
7. Verify correct position of the coil fluoroscopically. If coil position or placement is not satisfactory, the coil may be retracted into the catheter and re-deployed so long as there is no significant resistance.
NOTE: It may be possible to perform a test injection of contrast media using a Tuohy-Borst Sidearm Adapter while the delivery wire is in the catheter.
NOTE: If the size of the coil is not correct, gently remove the entire delivery wire and coil. Do not use the coil again.
8. If the coil position is correct, use the torque device to turn the delivery wire counterclockwise 8-10 times, until coil detachment can be either felt or visualized under fluoroscopy. **(Fig 8)**
NOTE: It is recommended that the junction remain just inside the tip of the catheter.
NOTE: Do not advance the delivery wire after the coil is detached.
9. Gently remove the delivery wire after coil detachment.
10. Additional Retractable or pushable coils may be required to achieve permanent occlusion of the vessel.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC016-031_A2YR) – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: (ec.europa.eu/tools/eudamed).

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2024
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd