COOK

MEDICAL

CE
0088

EN
9

Retracta™ Detachable Embolization Coils

Instructions for Use

DA
11

Retracta™ adskillelige emboliseringscoils

Brugsanvisning

DE
13

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

- 1 Vävning
- 2 Hylsa
- 3 Kateter

9

- 1 Weaving
- 2 Sheath
- 3 Catheter
- 4 High radial force coil
- 5 High radial force coil deployed (scaffold)
- 6 Platinum coil (soft coil)
- 7 Platinum coil deployed, cross-sectional occlusion completed

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RETRACTA™ DETACHABLE EMBOLIZATION COIL

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Retractable Detachable Embolization Coil is a platinum coil with spaced synthetic fibers, and is attached to a delivery wire. (**Figs. 1 & 2**) The Retractable Detachable Embolization Coil is designed to be delivered under fluoroscopy to the target vessel. The Retractable Detachable Embolization Coil system allows the coil to be advanced and fully repositioned before the coil is finally deployed.

INTENDED USE

The Retractable Detachable Embolization Coil is intended for arterial and venous embolization in the peripheral vasculature. The coil delivery system provides safe delivery of embolization coils when correct positioning is especially critical.

CONTRAINDICATIONS

None known

WARNINGS

- Positioning of embolization coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- The Retractable Detachable Embolization Coil is not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.

PRECAUTIONS

- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
- This product is intended for use by physicians trained and experienced in arterial and venous vessel embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.

PRODUCT RECOMMENDATIONS

The following catheters are recommended for use with Retractable Detachable Embolization Coil:

- HNB(R)4.0-35
- HNB(R)5.0-35
- HNB(R)5.0-38
- SCBR-5.0-38
- SCBR-4.0-38

The PTBYC-RA is recommended for use with the Retractable Detachable Embolization Coil to perform a test injection of contrast media.

MRI INFORMATION

Nonclinical testing has demonstrated that the Retractable Detachable Embolization Coil is **MR Conditional** according to ASTM F2503. A patient with this coil can be scanned safely after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient of 1600 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

In nonclinical testing, the Retractable Detachable Embolization Coil produced a maximum temperature rise of 1.1 °C (scaled to an SAR of 2.0 W/kg) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5).

Image Artifact

MR image quality may be compromised if the area of interest is within approximately 9 mm of the position of the Retractable Detachable Embolization Coil, as found during non clinical testing using the sequence: T1-weighted,

spin echo and gradient echo pulse sequence in a MR 3.0 Tesla system (General Electric Healthcare Excite, HDx). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic embolization coil.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

INSTRUCTIONS FOR USE

1. Perform an angiogram and measure the diameter of the vessel to be occluded.
2. Firmly grasping the loading cartridge between thumb and forefinger, introduce the metal end of the loading cartridge into the base of the catheter hub. Lock the loading cartridge onto the catheter hub by turning the Luer lock adapter clockwise. (**Fig. 3**)
3. Unlock the white safety clip and advance at least 30 cm of the delivery wire into the catheter so that the coil is fully loaded in the catheter. (**Fig. 4**)
NOTE: Advancing the delivery wire slowly allows the junction to be seen more easily and reduces the risk of damaging it.
NOTE: If significant resistance is encountered during coil advancement, do not continue advancing. Retract the delivery wire slightly, then gently re-advance it. If there is still significant resistance, withdraw the delivery wire from the catheter and try using a new coil with a shorter length.
NOTE: Do not turn the delivery wire counterclockwise during advancement; the coil may be unintentionally detached.
4. Unlock the loading cartridge from the catheter hub by turning the Luer lock adapter counter-clockwise.
5. Remove the delivery wire holder and the metal loading cartridge from the catheter hub while holding the delivery wire stationary. Remove the torque device and reserve it for use later in the procedure.
6. Under fluoroscopic visualization, slowly advance the delivery wire until the entire length of the coil exits the distal end of the catheter. Ensure that the junction remains positioned just inside the catheter tip. (**Fig. 5**)
NOTE: Advancing the delivery wire slowly allows the junction to be seen more easily and reduces the risk of damaging it.
NOTE: If significant resistance is encountered during coil advancement, do not continue advancing. Retract the delivery wire slightly, then gently re-advance it. If there is still significant resistance, withdraw the delivery wire from the catheter and try using a new coil with a shorter length.
NOTE: Do not turn the delivery wire counterclockwise during advancement; the coil may be unintentionally detached.
7. Verify correct position of the coil fluoroscopically. If coil position or placement is not satisfactory, the coil may be retracted into the catheter and re-deployed so long as there is no significant resistance.
NOTE: It may be possible to perform a test injection of contrast media using a Tuohy-Borst Sidearm Adaptor while the delivery wire is in the catheter.
NOTE: If the size of the coil is not correct, gently remove the entire delivery wire and coil. Do not use the coil again.
8. If the coil position is correct, use the torque device to turn the delivery wire counterclockwise 8-10 times, until coil detachment can be either felt or visualized under fluoroscopy. (**Fig 6**)
NOTE: It is recommended that the junction remain just inside the tip of the catheter.
NOTE: Do not advance the delivery wire after the coil is detached.
9. Gently remove the delivery wire after coil detachment.
10. Additional Retractable or pushable coils may be required to achieve permanent occlusion of the vessel.

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion. The combination of the coaxial technique and either the anchor or scaffold technique significantly enhances stability of coil deployment.

- Coaxial technique: The use of outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (**Fig. 7**)
- Anchor technique: The anchor technique provides safe and distal occlusion when there is a question about instability of coils. At least 2 cm of a coil is advanced into the side branch, which is normally sacrificed. The rest of the coil is then deployed just proximal to that side branch, and additional coils are packed. (**Fig. 8**)
- Scaffold technique: The scaffold technique is used for high-flow vessels when there is concern about migration of a softer coil. A high radial force coil is placed initially. Then, several high radial force coils or soft coils may be packed within the scaffold. (**Fig. 9**)

In general, the first coil selected should have a diameter that is 20% larger, or 2 mm oversized, than the vessel that is being occluded.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2014
Bloomington, IN 47404 U.S.A.