EN 3	Mandril Wire Guides Instructions for Use
CS 3	Mandrénové vodící dráty Návod k použití
DA 4	"Mandril"-kateterledere Brugsanvisning
DE 5	Mandrin-Führungsdrähte Gebrauchsanweisung
EL 6	Συρμάτινοι οδηγοί mandril Οδηγίες χρήσης
ES 7	Guías con mandril Instrucciones de uso
FR 8	Guides mandrins Mode d'emploi
HU 9	Szár-vezetődrót Használati utasítás
IT 10	Guide a mandrino Istruzioni per l'uso
NL 11	Mandrijnvoerdraden Gebruiksaanwijzing
NO 12	Mandrengledevaiere Bruksanvisning
PL 13	Prowadniki typu mandryn Instrukcja użycia
PT 14	Fios guia de mandril Instruções de utilização
SV 15	Mandrinledare Bruksanvisning

MANDRIL WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Mandril Guides are made with a core mandril of stainless steel or nitinol material that runs the entire length of the product. Attached at the distal end is a coil of stainless steel, platinum, or palladium of varying length. Mandril Wire Guides come in a variety of lengths, diameters, and tip configurations. Some devices are coated with TFE coating to ease coaxial use with other devices. Some devices are provided with two usable ends. Refer to the product label for product specifications.

INTENDED USE

Mandril Wire Guides are used to facilitate the placement of devices during diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Use medical imaging when you manipulate the wire guide. Do not advance or torque the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching a syringe with heparinized saline or sterile water to the fitting of the wire guide holder. Inject enough solution to wet the wire surface entirely.

NOTE: If flushing through the wire guide holder is not possible, remove the wire guide from the holder and place it in a bowl of heparinized saline or sterile water, or wet the wire guide surface over the entire length using gauze that has been moistened with heparinized saline solution.

2. Carefully remove the wire guide from the holder.
3. If needed, insert a wire guide insertion tool (provided) through the valve assembly or hub of the guiding sheath or other interventional device. Insert the tip of the wire guide through the insertion tool and advance the wire guide to the desired location.
4. Attach a torque device to the wire guide (if provided).
5. Standard wire guide techniques may now be employed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd