
COOK**MEDICAL****CE**
0123

EN 7	Nester® Embolization Coils and Microcoils Instructions for Use
BG 11	Спирали и микроспирали за емболизация Nester® Инструкции за употреба
CS 15	Embolizační spirálky a mikrospirálky Nester® Návod k použití
DA 19	Nester® emboliseringscoils og microcoils Brugsanvisning
DE 23	Nester® Embolisationsspiralen und Microcoils Gebrauchsanweisung
EL 28	Σπειράματα και Microcoils εμβολισμού Nester® Οδηγίες χρήσης
ES 33	Espirales y microespirales de embolización Nester® Instrucciones de uso
FI 38	Nester®-embolisaatiokoilit ja -mikrokoilit Käyttöohjeet
FR 42	Coils et microcoils d'embolisation Nester® Mode d'emploi
HU 46	Nester® embolizációs spirálok és mikrospirálok Használati utasítás
IT 51	Spirali e microspirali per embolizzazione Nester® Istruzioni per l'uso
LT 55	„Nester®“ embolizavimo spiralės ir mikrospiralės Naudojimo instrukcija
LV 59	Nester® embolizācijas spirāles un mikrospirāles Lietošanas instrukcija
NL 64	Nester® embolisatiecoils en Microcoils Gebruiksaanwijzing
NO 68	Nester® emboliseringsspiraler og mikrospiraler Bruksanvisning
PL 72	Spirale i mikrospirale embolizacyjne Nester® Instrukcja użycia
PT 77	Espirais e microcoils de embolização Nester® Instruções de utilização
RO 81	Bobine de embolizare și microbobine Nester® Instrucțiuni de utilizare
SK 86	Embolizačné cievky a mikrocievky Nester® Návod na použitie
SL 90	Embolizacijske tuljave in mikrotuljave Nester® Navodila za uporabo
SV 94	Nester® emboliseringsspiraler och mikrospiraler Bruksanvisning
TR 98	Nester® Embolizasyon Sarmalları ve Mikro Sarmallar Kullanma Talimatı

T _ N E C 2 - M _ R E V O

1


2

- 1 Straight coil length
2 Catheter inner diameter
- 1 Права дължина на спиралата
2 Вътрешен диаметър на катетъра
- 1 Délka rovné spirálky
2 Vnitřní průměr katetru
- 1 Lige coillængde
2 Katetrets indvendige diameter
- 1 Länge der geraden Spirale
2 Innendurchmesser des Katheters
- 1 Μήκος σπειράματος σε ευθεία θέση
2 Εσωτερική διάμετρος καθετήρα
- 1 Longitud de la espiral recta
2 Diámetro interior del catéter
- 1 Koilin pituus suorana
2 Katetrin sisäläpimitta
- 1 Longueur de coil droite
2 Diamètre interne du cathéter
- 1 Egyenes spirálhossz
2 Katéter belső átmérője
- 1 Lunghezza della spirale diritta
2 Diametro interno del catetere

- 1 Tiesios spirālės ilgis
2 Kateterio vidinis skersmuo
- 1 Taisnas spirāles garums
2 Katetra iekšējais diametrs
- 1 Lengte rechte coil
2 Binnendiameter katheter
- 1 Rett spirallengde
2 Kateterets indre diameter
- 1 Długość prostej spirali
2 Wewnętrzna średnica cewnika
- 1 Comprimento da espiral reta
2 Diâmetro interno do cateter
- 1 Lungimea bobinei drepte
2 Diametrul interior al cateterului
- 1 Dĺžka rovnej cievky
2 Vnútroň priemer katétra
- 1 Dolžina poravnane tuljave
2 Notranji premer katetra
- 1 Spiralens längd i utsträckt tillstånd
2 Kateterens innerdiameter
- 1 Düz sarmal uzunluğu
2 Kateter iç çapı

3

4


5

6

- 0,018 Inch Microcoils
Микроспирали 0,018 inch
Mikrospirálky 0,018 inch
0,018 inch microcoils
0,018 Inch Microcoils
Microcoils 0,018 inch
Microespirales de 0,018 inch
0,018 inch -kokoiset mikrokoilit
Microcoils de 0,018 inch
0,018 inch méretű mikrospirálok
Microspirali da 0,018 inch


- 0,018 inch mikrospirālės
0,018 inch mikrospirāles
Microcoils van 0,018 inch
Mikrospiraler på 0,018 inch
Mikrospirale o średnicy 0,018 inch
Microcoils de 0,018 inch
Microbobine de 0,018 inch
0,018 inch mikrocievky
Mikrotuljave velikosti 0,018 inch
Mikrospiraler, 0,018 inch
0,018 inch Mikro Sarmallar

7


8

0.035 Inch Coils	0,035 inch spirālēs
Спирали 0,035 inch	0,035 inch spirāles
Spirálky 0,035 inch	Coils van 0,035 inch
0,035 inch coils	Spiraler på 0,035 inch
0,035 Inch Spiralen	Spirale 0,035 inch
Σπειράματα 0,035 inch	Espiraies de 0,035 inch
Espiraies de 0,035 inch	Bobine de 0,035 inch
0,035 inch -kokoiset koilit	Cievky veľkosti 0,035 inch
Coils de 0,035 inch	Tuljave velikosti 0,035 inch
0,035 inch méretű spirálok	Spiraler, 0,035 inch
Spirali da 0,035 inch	0,035 inch Sarmallar


9

1 Weaving	1 Audimas
2 Sheath	2 Vamzdēlis
3 Catheter	3 Kateteris
1 Преплитане	1 Aušana
2 Защитна обвивка	2 Ievadapvalks
3 Катетър	3 Katetrs
1 Pohyb tam a zpět	1 Wevende beweging
2 Sheath	2 Sheath
3 Objednací číslo katetru	3 Katheterlengte
1 AVævning	1 Fletting
2 Sheath	2 Hylse
3 Kateterlængde	3 Kateterlengde
1 Hin- und Herbewegung	1 Ruch tam i z powrotem
2 Schleuse	2 Koszulka
3 Katheter	3 Cewnik
1 Κίνηση με ελιγμούς	1 Movimento circular
2 Θηκάρι	2 Bainha
3 Καθετήρας	3 Cateter
1 Zigzagueo	1 Împletire
2 Vaina	2 Teacă
3 Catéter	3 Cateter
1 Mutkittelu ristiin rastiin	1 Tkanie
2 Holkki	2 Puzdro
3 Katetri	3 Katéter
1 Mouvement de va-et-vient circulaire	1 Prepletanje
2 Gaine	2 Tulec
3 Cathéter	3 Kateter
1 Kõrkõrõs ide-oda mozgatás	1 Vævning
2 Hűvely	2 Hylsa
3 Katéter	3 Kateter
1 Movimento alternato	1 Dokuma
2 Introduttore	2 Kilif
3 Catetere	3 Kateter

10

- | | |
|---------------------------------------|-------------------------|
| 1 Weaving | 1 Audimas |
| 2 Sheath | 2 Vamzdelis |
| 3 Catheter | 3 Kateteris |
| 1 Преплитане | 1 Aušana |
| 2 Защитна обвивка | 2 Ievadapvalks |
| 3 Катетър | 3 Katetrs |
| 1 Pohyb tam a zpět | 1 Wevende beweging |
| 2 Sheath | 2 Sheath |
| 3 Objednací číslo katetru | 3 Katheterlengte |
| 1 Vævning | 1 Fletting |
| 2 Sheath | 2 Hylse |
| 3 Kateterlængde | 3 Kateterlengde |
| 1 Hin- und Herbewegung | 1 Ruch tam i z powrotem |
| 2 Schleuse | 2 Koszulka |
| 3 Katheter | 3 Cewnik |
| 1 Κίνηση με ελιγμούς | 1 Movimento circular |
| 2 Θηκάρι | 2 Bainha |
| 3 Καθετήρας | 3 Cateter |
| 1 Zigzaguelo | 1 Împletire |
| 2 Vaina | 2 Teacă |
| 3 Catéter | 3 Cateter |
| 1 Mutkittelu ristiin rastiin | 1 Tkanie |
| 2 Holkki | 2 Puzdro |
| 3 Katetri | 3 Katéter |
| 1 Mouvement de va-et-vient circulaire | 1 Prepletanje |
| 2 Gaine | 2 Tulec |
| 3 Cathéter | 3 Kateter |
| 1 Kõrkõrõs ide-oda mozgatás | 1 Vävning |
| 2 Hűvely | 2 Hylsa |
| 3 Katéter | 3 Kateter |
| 1 Movimento alternato | 1 Dokuma |
| 2 Introduttore | 2 Kilif |
| 3 Catetere | 3 Kateter |

11

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NESTER® EMBOLIZATION COILS AND MICROCOILS

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Nester® 0.035" Embolization Coils and 0.018" Microcoils are made of platinum with spaced synthetic fibers, and are supplied preloaded in a loading cartridge. They are designed to be delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter.

Performance Characteristics

- Coil design and soft platinum material allow device to form a tight occluding mass.
- Radiopaque platinum material allows for visibility under fluoroscopy.
- Synthetic fibers spaced along coil increase thrombogenicity.
- Wide range of diameters (2-20 mm) and lengths (2-20 cm) available.
- MR conditional at 3 T and 1.5 T.
- Device comes pre-loaded in a loading cartridge with an MLLA adapter that can connect to a standard angiographic catheter hub.
- Coils can be deployed using a standard wire guide and/or contrast/saline flush.
- Coil design maintains radial force to prevent migration when oversized by 10-20% for arteries or 30-50% for veins.
- Obstructs or reduces blood flow in vasculature.
- Can be used alone or in combination with other platinum or Inconel alloy pushable and detachable embolization coils.
- Multiple coils may be used to achieve complete occlusion.
- Length is identified on product label by **Fig. 1**.
- The lifetime of the implant is intended to be permanent. Testing has been completed to assure that the device will perform as intended for the lifetime of the patient or duration of device implantation.

Device Compatibility

- 0.035" Nester Coils are compatible with 0.035" wire guides and angiographic catheters with an 0.035" inner diameter.
- 0.018" diameter Nester Microcoils are compatible with 0.018" wire guides and microcatheters with an inner diameter between 0.018" and 0.025".
- Catheter inner diameter is identified on product labels by **Fig. 2**.

Qualitative and Quantitative Information

Device Materials		Weight (grams)
Platinum-tungsten (a metal alloy)	Wire for coiling	Up to 0.819
Nylon	Fiber	Up to 0.008

Expected Lifetime of the Device

The expected lifetime of the Nester Embolization Coil is permanent, with an indwell duration spanning the life of the patient.

Patient Population

The target population for Nester Embolization Coils is adult patients who require vascular embolization procedures that induce cessation or reduction of blood flow in peripheral arterial and venous vessels.

Intended User

This device is intended for use by physicians who are experienced in interventional techniques. Embolization procedures that require the use of the Embolization Devices should only be performed by licensed clinical practitioners.

Contact with Body Tissue

The Nester Embolization Coils are for single use, categorized as implant devices in permanent contact with circulating blood and peripheral vasculature.

Operating Principle

Nester Embolization Coils are made of platinum with spaced synthetic fibers, and are supplied preloaded in a loading cartridge. They are designed to be delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter. Once in the vessel, coil design maintains radial force to prevent migration. The coil and its synthetic fibers induce

thrombogenicity to form a tight occluding mass to obstruct or reduce blood flow in the vasculature. Multiple coils may be used to achieve complete occlusion.

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross-sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion. The combination of the coaxial technique and either the anchor or scaffold technique significantly enhances stability of coil deployment.

- **Coaxial technique:** The use of an outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (**Fig. 9**)
- **Anchor technique:** The anchor technique provides safe and distal occlusion when there is a question about instability of coils. At least 2 cm of a coil is advanced into the side branch, which is normally sacrificed. The rest of the coil is then deployed just proximal to that side branch, and additional coils are packed. (**Fig. 10**) **NOTE:** Use a coil longer than 10 cm for anchor technique.
- **Scaffold technique:** The scaffold technique is used for high-flow vessels when there is concern about migration of a softer coil. A high radial force Inconel coil is placed initially. Then, several Inconel coils or platinum coils (soft coils) may be packed within the scaffold. (**Fig. 11**)

In general, the first coil selected should have a diameter that is 20% larger, or at least 2 mm oversized, than the vessel that is being occluded.

INTENDED USE

This permanent embolic device is a fibered, metallic coil that is intended to cause cessation or reduction of blood flow in peripheral arteries and veins by means of physical occlusion and thrombosis induced by the coil and its synthetic fibers.

INDICATIONS FOR USE

Nester Embolization Coils and Microcoils are indicated for arterial and venous embolization in the peripheral vasculature. They may be used for treatment of vascular abnormalities, treatment of hemorrhage, devascularization of benign and malignant tumors, devascularization of benign tissue, flow redistribution, and endoleak management.

CLINICAL BENEFITS

Successful adjunctive (e.g., pre-surgical), palliative, or curative

- treatment of vascular abnormalities,
- treatment of hemorrhage,
- devascularization of benign or malignant tumors,
- devascularization of benign tissue,
- flow redistribution, or
- endoleak management

by complete or partial vessel occlusion via a minimally invasive procedure.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that single and multiple Nester Embolization Coils and Microcoils are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1700 gauss/cm (17 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Nester Embolization Coils and Microcoils are not expected to result in a temperature rise of more than 3.0°C after 15 minutes of continuous scanning.

The image artifact extends approximately 5.7 mm from the Nester Embolization Coils and Microcoils as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

Absolute contraindications to percutaneous transcatheter embolization:

General:

- Active systemic infection
- Clinically manifest hyperthyroidism

In the context of chemoembolization (TACE) or radioembolization (TARE) for hepatic malignancy:

- Active hepatic infection
- Acute variceal bleeding in previous 2 weeks
- Advanced cirrhosis
- Arteriovenous fistulas, untreatable
- Bilirubin >5 mg/dL
- Capecitabine treatment within 2 months before TARE treatment with resin spheres
- Child-Pugh C liver disease or Child-Pugh ≥ 8
- Contraindications for arteriography and/or transarterial interventions
- Extrahepatic metastases
- Inability to catheterize hepatic artery
- Insufficient functional liver reserve (for TARE)
- Large tumor burden (e.g., bilobar tumor affecting a large portion of the liver, >50% liver replaced by tumor)
- Liver failure, fulminant
- Lung shunting, severe (TARE)
- Nontarget embolization of the gastrointestinal tract, anticipated (TARE)
- Portal vein obstruction (TACE)
- Reversed portal flow

WARNINGS

- Neurological deficit, ischemic stroke, spinal infarct, or ischemic infarct can occur from occlusion of normal vessels by emboli.
- The devices are made of a 92% platinum, 8% tungsten alloy. The devices should not be implanted in patients who have allergies to any of the coil components.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Nester Embolization Coils and Microcoils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Warning: The user should consider potential risks associated with radiation exposure.

PRECAUTIONS

- This product is intended for use by healthcare providers experienced in arterial and venous vessel embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
- If using a 0.018 inch Microcoil, ensure that the delivery catheter has an internal diameter (ID) of 0.018 to 0.025 inch.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Access-related and catheter-angiographic side effects (access-site nerve injury, contrast agent-induced nephropathy, contrast reaction, embolism [e.g., air, atheromatous debris, thrombus], hematoma, infection, nausea/vomiting, vessel occlusion, pain, pseudoaneurysm/arteriovenous fistula, subintimal passage, thrombosis, vasospasm, vasovagal syncope, vessel dissection, vessel injury, vessel perforation)
- Allergic reaction to coil components
- Bile leak/bile duct injury/biloma
- Bleeding/hemorrhage
- Coil migration
- Embolization target ischemia/infarction/inflammation (i.e., the clinical effects, whether desired or not, that result from vessel occlusion in the intended target area)
- Fistula

- Infection
- Liver insufficiency/failure
- Nontarget embolization causing ischemia or infarction
- Pleural effusion, ascites
- Postembolization syndrome (systemic signs and symptoms in addition to anatomy-specific signs, symptoms and abnormal laboratory values of postembolization syndrome that may affect particular patient populations)
- Procedure-related mortality
- Renal dysfunction (e.g., increased renal laboratory values, renal insufficiency/failure)
- Thrombosis

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

This device is accompanied by an implant card, that should be given to the patient after it has been completed by the Healthcare Professional.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

For 0.018 Inch Microcoils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Push the loading cartridge completely into the delivery catheter. **(Fig. 3)**
3. Advance the Luer lock connector fitting toward the catheter hub and lock into place. **(Fig. 4)**
4. Use the pusher stylet to load the Nester Embolization Microcoil into the delivery catheter. **(Fig. 5)** **NOTE:** The stylet must be pushed as far as possible into the loading cartridge to ensure proper loading. Remove the pusher stylet and loading cartridge.
5. To obtain secure placement of the Nester Embolization Microcoil, we recommend the use of the pusher stylet to push the Microcoil further into the delivery catheter. **(Fig. 6)**
6. Once coil is positioned into delivery catheter, verify position of catheter tip prior to deployment. Coil design permits delivery into the target vessel by saline/contrast flush or by push technique using an appropriately sized wire guide or pusher.
7. Perform final angiogram to confirm coil position within target vessel.

For 0.035 Inch Coils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by turning Luer lock adapter clockwise. **(Fig. 7)**
3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 centimeters of the angiographic catheter. **(Fig. 8)** Remove the wire guide and loading cartridge.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.
5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC016-031_A1YP) – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: (ec.europa.eu/tools/eudamed).

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd