- EN**
3 **NGage® Nitinol Stone Extractor**
Instructions for Use
- BG**
5 **Нитинолов екстрактор на камъни NGage®**
Инструкции за употреба
- CS**
8 **Nitinolový extraktor kamenů NGage®**
Návod k použití
- DA**
11 **NGage® stenekstraktor af nitinol**
Brugsanvisning
- DE**
13 **NGage® Nitinol-Steinextraktor**
Gebrauchsanweisung
- EL**
16 **Εξαγωγέας λίθων από νιτινόλη NGage®**
Οδηγίες χρήσης
- ES**
19 **Extractor de cálculos de nitinol NGage®**
Instrucciones de uso
- ET**
22 **Nitinoolist kivieemaldi NGage®**
Kasutusjuhend
- FI**
25 **Nitinolista valmistettu NGage®-kivenpoistolaite**
Käyttöohjeet
- FR**
27 **Extracteur de calculs en nitinol NGage®**
Mode d'emploi
- HR**
30 **Nitinolski ekstraktor kamenaca NGage®**
Upute za uporabu
- HU**
33 **NGage® nitinol kőeltávolító**
Használati utasítás
- IT**
35 **Estrattore di calcoli in nitinol NGage®**
Istruzioni per l'uso

- LT**
38
NGage® nitinolo akmenų šalinimo priemonė
Naudojimo instrukcija
- LV**
41
NGage® nitinola akmeņu izvilcējs
Lietošanas instrukcija
- NL**
43
NGage® steenextractor van nitinol
Gebruiksaanwijzing
- NO**
46
NGage®-stenekestaktor i nitinol
Bruksanvisning
- PL**
49
Nitynolowy ekstraktor złożów **NGage®**
Instrukcja użycia
- PT**
52
Extrator de cálculos em nitinol **NGage®**
Instruções de utilização
- RO**
54
Extractor de calculi din nitinol **NGage®**
Instrucțiuni de utilizare
- SK**
57
Nitinolový extraktor kameňov **NGage®**
Návod na použitie
- SL**
60
Nitinolni ekstraktor kamnov **NGage®**
Navodila za uporabo
- SV**
62
NGage® stenudragare av nitinol
Bruksanvisning
- TR**
65
NGage® Nitinol Taş Ekstraktörü
Kullanma Talimatı

NGAGE® NITINOL STONE EXTRACTOR

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The NGage® Nitinol Stone Extractors are manufactured in 1.7 and 2.2 Fr sizes and in a length of 115 cm. The modified basket is designed to open to a 10 millimeter diameter, and the non-modified baskets open to an 8 mm diameter. The 115 cm shaft houses a 3-wire nitinol basket, and a handle facilitates opening and closing the basket. A lever on the handle slides to control the open or closed position of the basket or extractor. The NGage® Nitinol Stone Extractor includes an adapter and an inserter.

Performance Characteristics

- Nitinol wire basket retains shape after multiple actuations and stone removals.
- Device shaft consisting of inner stainless-steel coil to increase torque response.
- Length of product allows targets in the kidney to be manipulated and removed.
- Stainless steel braid is incorporated into the sheath material for resistance to kinking.
- Ergonomic handle design allows for use with one hand with open and close directions labeled.
- Tuohy-Borst Adapter resists backflow of fluid from the working channel of the endoscope around the extractor.
- Inserter to assist insertion of the extractor into the working channel of the endoscope.

Device Compatibility

- These devices are compatible with endoscope working channels with a minimum inner diameter of 3.0 Fr (1 mm) and a maximum length of 80 cm.
- The Tuohy-Borst Adapter can be connected to any standard female luer lock fitting or bulb-type endoscopic fitting.
- The inserter tip is compatible with the Tuohy-Borst Adapter and allows for insertion of the extractor.
- These devices are compatible with sterile saline and sterile water.

Patient Population

These devices are not intended to be used in any specific target population unless outlined on the labelling, in the instructions, and/or in promotional materials. Patient population could include any person requiring treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues comprising the urinary tract, which may include the urethra, bladder, ureters, and kidneys.

Operating Principle

The NGage Nitinol Stone Extractor is a single-use device. The design enables the physician to engage, reposition, release and extract stones in the urinary tract. An optional Tuohy-Borst adapter is provided that can be connected to the endoscope working channel to resist fluid backflow around the extractor. An optional inserter is also provided that may be used to facilitate insertion of the extractor into the Tuohy-Borst adapter. The device is constructed from nitinol, which is flexible, kink resistant, and durable. The device is manufactured in 1.7 and 2.2 Fr sizes and in a length of 115 cm. The basket is designed to open in either an 8 or 10 mm diameter. A handle facilitates opening and closing of the basket. Sliding the thumb lever in the "OPEN" direction opens the basket, while sliding the thumb lever in the "CLOSE" direction closes the basket.

INTENDED USE

The NGage Nitinol Stone Extractor is intended for stone manipulation and removal in the urinary tract through the working channel of an endoscope during cystoscopic or ureteroscopic procedures by actuating a nitinol wire basket.

DEVICE AND CLINICAL BENEFITS

To manipulate and remove stones from the urinary tract through the working channel of an endoscope.

INDICATIONS FOR USE

These devices are used to assist in extraction of stones (urolithiasis) from the urinary tract.

CONTRAINDICATIONS

No contraindications.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- User should be familiar with and have experience in urological endoscopic surgery.
- This device is conductive. Avoid contact with any electrified instrument.
- Close the basket before removing from the tray.
- Assess stone prior to instrument deployment to ensure that the stone is not too large to be removed through anatomy.
- If resistance is encountered while attempting to remove stone, release the stone. Do not exert excessive force on the device.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding e.g. Hematoma, Hematuria, Hemorrhage
- Death
- Deep venous thrombosis
- Fever
- Fluid overload
- Hypotension
- Hypothermia
- Inflammation e.g. Systemic inflammatory response syndrome
- Infection e.g. Bacteremia, Bacteriuria, Peri-renal Abscess, Pneumonia, Pyelonephritis, Sepsis, Urinary Tract Infection
- Nausea
- Pain
- Pulmonary Embolism
- Tissue, organ or structure injury e.g. Arteriovenous fistula, False tract, Infundibular stenosis, Pleural injury, Pseudoaneurysm, Renal collecting system injury, Selective arterial embolization, Surrounding organs/tissues injury, Urinary fistula
- Ureteral injury e.g. Ureteral avulsion, Ureteral mucosal injury, Ureteral mucosal tear, Ureteral obstruction, Ureteral perforation, Ureteral stricture
- Urinary tract issues e.g. Dysuria, Extrarenal stone migration, Fluid collection, Hydronephrosis, Steinstrasse, Urinary frequency, Urinary leakage, Urinary retention, Urinary urgency

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Move the thumb lever on the handle to the "Close" position to close the basket and remove the extractor from its packaging.
2. The inserter may be placed securely into the Tuohy-Borst Adapter on the scope, as desired.
3. Then carefully insert the retracted extractor into the scope.
4. Under direct visualization, locate the stone and advance the extractor to a location just proximal to it. Slide the thumb lever on the handle to the "Open" position to deploy the basket near the stone and maneuver the device until it is within the basket.
5. Once correctly positioned, draw the basket down around the stone, and slide the thumb lever into the "Close" position.
6. Maintaining position of the stone within the basket, retract the scope and basket.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd