COOK**MEDICAL****CE**
0123

EN 3	Neff Percutaneous Access Sets Instructions for Use
CS 4	Soupravy pro perkutánní přístup Neff Návod k použití
DA 6	Neff-sæt til perkutan adgang Brugsanvisning
DE 8	Neff perkutane Punktions-Sets Gebrauchsanweisung
EL 11	Σετ διαδερμικής πρόσβασης Neff Οδηγίες χρήσης
ES 13	Equipos de acceso percutáneo Neff Instrucciones de uso
ET 16	Perkutaanse juurdepääsu komplektid Neff Kasutusjuhend
FI 18	Perkutaaniseen sisäänvientiin tarkoitettut Neff-setit Käyttöohjeet
FR 20	Sets d'accès percutané Neff Mode d'emploi
HR 22	Kompleti za perkutani pristup Neff Upute za uporabu
HU 24	Neff perkután punkciós szettek Használati utasítás
IT 27	Set per accesso percutaneo Neff Istruzioni per l'uso
LT 29	„Neff“ perkutaninės prieigos rinkiniai Naudojimo instrukcija
LV 31	Neff perkutānās piekļuves komplekti Lietošanas instrukcija
NL 33	Neff perkutane toegangssets Gebruiksaanwijzing
NO 36	Neff perkutane tilgangssett Bruksanvisning
PL 38	Zestawy do dostępu przezskórnego Neff Instrukcja użycia
PT 40	Conjuntos de acesso percutâneo Neff Instruções de utilização
RO 43	Seturi de acces percutanat Neff Instrucțiuni de utilizare
SK 45	Súpravy na perkutánný prístup Neff Návod na použitie
SL 47	Kompleti za perkutani dostop Neff Navodila za uporabo
SV 49	Neff set för perkutan åtkomst Bruksanvisning
TR 51	Neff Perkütan Erişim Setleri Kullanma Talimatı

T _ N P A S - M _ R E V O

NEFF PERCUTANEOUS ACCESS SETS

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Neff Percutaneous Access Sets are comprised of three main components: the introducer set, the wire guide, and the disposable Chiba needle or trocar needle. The wire guide mandril is either nitinol (-NT) or stainless steel.

The introducer sets consist of a radiopaque sheath, an inner dilator, and a blunt stainless steel stiffening cannula. Some outer sheaths (-RH) have a platinum/iridium marker band and are designed with a hydrophilic coating. The stiffening cannula, dilator and sheath all have a transitional fit and are compatible with Luer-locking hubs.

The nitinol or stainless steel mandril wire guide is designed with a flexible spring coil tip manufactured from either platinum or stainless steel. The mandril is 0.018 inch in diameter and has a length of 60 cm.

The 22-gauge disposable Chiba needle or 21-gauge trocar needle is manufactured from stainless steel. Both needles are manufactured as a two part needle with a needle cannula and stylet in lengths of 15 or 20 cm. The Chiba needle has a lancet bevel and the trocar needle has a trocar tip. The Chiba 20-centimeter length needle has an echo-tip characterized by dimpling along the end of the needle.

Performance Characteristics

- Marker Bands
- Radiopaque Sheath
- Tapered Distal Tip of Introducer Set
- Hydrophilic Coating on Sheath
- Flexible Tip Wire Guide
- Sheath accepts a 0.038 inch working wire guide

Device Compatibility

- Sheath compatible with 0.018 and 0.038-inch diameter wire guide
- Compatible with 4 Fr devices

Patient Population

This device is used in adult patients with non-vascular medical conditions requiring percutaneous access for diagnostic or interventional procedures as directed by a physician or licensed healthcare provider.

Intended User

This device is intended for use by physicians and other properly licensed providers trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous catheters should be employed.

Contact with Body Tissue

This device is externally communicating contacting non-vascular tissue.

Operating Principle

Under imaging guidance (e.g., Ultrasound or CT scan) the needle of the device provides percutaneous access to an intended target anatomic site to facilitate a working channel that is used to introduce a 0.018-inch diameter wire guide into the target site. The introducer assembly is then placed over the 0.018-inch diameter wire guide. Once the introducer assembly is in place, the stiffening cannula and dilator are removed. The 0.018-inch diameter wire guide may be removed or left in place as a 0.038-inch diameter wire guide is placed through the introducer sheath.

INTENDED USE

This device is comprised of an introducer set, wire guide, and an access needle and is intended for single-micropuncture percutaneous access to facilitate placement of a working wire guide or introducing catheter for non-vascular interventional radiology procedures.

INDICATIONS FOR USE

This device is indicated for percutaneous access (for diagnosis and/or medical treatment) in non-vascular interventional radiology procedures.

CLINICAL BENEFITS

Gain percutaneous access to a non-vascular lumen, duct, viscus, abscess, fluid collection, or cavity.

CONTRAINDICATIONS

- Uncorrectable coagulopathy
- Bowel ischemia (for enteral access procedures)
- Peritonitis (for enteral access procedures)

- Portal hypertension or gastric varices (for enteral access procedures)

WARNINGS

- The mandril wire guide should not be withdrawn through the needle once placed beyond the needle tip.
- The Neff Percutaneous Access Sets should not be used in patients who are allergic to any of its materials. Nitinol wire guides should not be used in patients with nickel allergy.
- This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- Imaging guidance should be used for the placement of this device.
- This product must be used for non-vascular procedures only.
- Activate the hydrophilic coating, if present, by wetting the sheath with sterile water or saline. For best results, keep sheath surface wet during placement.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a cobalt alloy (or stainless steel alloy containing cobalt), which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to, the following:

- Bleeding, hematoma
- Fever
- Fistula
- Harms associated with device fracture
- Harms associated with wire guide separation, breakage, or unraveling
- Infection
- Inflammation, edema, necrosis, or scarring
- Injury, laceration or perforation of nerve, organ, tissue, or vessel
- Leakage of air, bile, blood, fluid, gastrointestinal contents, lymph, or urine
- Pain
- Sepsis, hypotensive shock, and death
- Tumor seeding of needle tract
- User injury
- Vasovagal reaction

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

NOTE: This procedure should be performed under imaging guidance.

1. Insert needle.
2. Advance the 0.018 inch wire guide through the needle.
3. Withdraw needle, leaving wire guide in place.
4. Advance the assembled introducer (sheath, dilator, and stiffening cannula) over the wire guide while maintaining control of the wire guide.
5. Remove the dilator and stiffening cannula, leaving the sheath in place. The 0.018 wire guide may be removed.
6. Advance a wire guide (up to 0.038 inch diameter) through the introducer sheath.
7. Remove the introducer sheath, leaving the wire guide in place.
8. Proceed with planned intervention.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING information

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC050-181_B12B) – to find the Summary of Safety and Clinical Performance (SSCP) on the Eudamed website (ec.europa.eu/tools/eudamed).

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2024
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd