COOK

MEDICAL

CE
0088

EN
3

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Bruksanvisning

T _ N P A S _ R E V O

APRIMA ACCESS SET WITH TRANSITIONLESS TECHNOLOGY

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Aprima Access Set includes an access needle, mandril wire guide and introducer sheath assembly, which consists of an outer introducer sheath, a sheath dilator and a stiffening cannula.

INTENDED USE

The Aprima Access Set is intended for single-puncture percutaneous access to facilitate placement of an .038 inch (0.97 mm) diameter working wire guide for interventional radiology procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous catheters should be employed.
- The mandril wire guide should not be withdrawn through the needle once placed beyond the needle tip.
- Imaging guidance should be used for the placement of this device.

INSTRUCTIONS FOR USE

1. Place the needle percutaneously from skin puncture site to the target anatomy.
2. Once the needle tip position is confirmed, remove the stylet from the percutaneous puncture needle.
3. Advance the mandril wire guide through the needle cannula into the target anatomy, confirming position with imaging guidance.
4. Maintaining position of the wire guide, remove the needle cannula.
5. Advance the introducer assembly over the wire guide until the distal tip is within target anatomy.
6. Maintaining position of the outer sheath, remove the wire guide, stiffening cannula and dilator. **NOTE:** If desired, the wire guide can remain in position as a safety wire.
7. Insert a .038 inch wire guide of choice through the lumen of the sheath and advance until the distal tip is within the target anatomy.
8. Maintaining position of wire guide, remove the introducer sheath.
9. Utilize wire guide for placement of further diagnostic or interventional devices.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.