EN
3

NSnare™ Stent Retriever

Instructions for Use

BG
5

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NSNARE™ STENT RETRIEVER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The NSnare™ Stent Retrievers are manufactured in 2.4 and 4.5 Fr sizes and in lengths of 115 and 65 cm, respectively. The retrieval loop (snare) is manufactured from flexible nitinol wire which measures 3.3 cm long and 2.5 cm wide in the opened position. A handle facilitates opening and closing with a thumb lever which slides to control the open and closed position of the wire loop. A Tuohy-Borst Adapter and inserter are included with the 2.4 Fr retriever only.

Performance Characteristics

- Nitinol wire loop retains shape after multiple actuations and object retrievals.
- Device shaft consisting of inner stainless-steel coil to increase torque response.
- Stainless steel braid is incorporated into the sheath material for resistance to kinking.
- Length of 2.4 Fr retriever allows targets in the kidney to be manipulated and removed.
- Ergonomic handle design allows for use with one hand with labeled open and close directions.
- Tuohy-Borst Adapter may be included to resist backflow of fluid from the working channel of the endoscope around the retriever.
- Inserter may be included to assist insertion of the retriever into the working channel of the endoscope.

Device Compatibility

- The 2.4 Fr retriever is compatible with endoscope working channels with a minimum inner diameter of 3.0 Fr (1 mm) and a maximum length of 80 cm.
- The 4.5 Fr retriever is compatible with endoscope working channels with a minimum inner diameter of 5.1 Fr (1.7 mm) and a maximum length of 45 cm.
- The Tuohy-Borst Adapter may be included and can be connected to any standard female luer lock fitting or bulb-type endoscopic fitting.
- The inserter tip is compatible with the Tuohy-Borst Adapter and allows for insertion of the closed retriever, if applicable.
- These devices are compatible with sterile saline and sterile water.

Patient Population

These devices are not intended to be used in any specific target population unless outlined on the labelling, in the instructions, and/or in promotional materials. Patient population could include any person requiring treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues comprising the urinary tract, which may include the urethra, bladder, ureters, and kidneys.

Operating Principle

The NSnare Stent Retrievers are introduced into the urinary tract through the working channel of an endoscope. An optional Tuohy-Borst adapter is provided with the 2.4 Fr NSnare Stent Retriever that can be connected to the endoscope working channel to resist fluid backflow around the retriever. An optional inserter is also provided with the 2.4 Fr device that may be used to facilitate insertion of the retriever into the Tuohy-Borst adapter. The retriever features an expandable wire loop

designed to ensnare stents or other foreign objects. The snare is controlled using the thumb lever on the handle. Sliding the thumb lever in the "OPEN" direction opens the snare, while sliding the thumb lever in the "CLOSE" direction closes the snare. Once the retriever is in position near the stent/object, the snare is opened, looped around it and closed. The snare can then be removed from the patient, along with the endoscope, and the process repeated as needed.

INTENDED USE

The NSnare Stent Retrievers are used to ensnare and retrieve stents and foreign objects from the urinary tract through the working channel of an endoscope during cystoscopic or ureteroscopic procedures by actuating a nitinol wire loop.

DEVICE AND CLINICAL BENEFITS

To ensnare and retrieve stents and other objects (foreign bodies) from the urinary tract through the working channel of an endoscope.

INDICATIONS FOR USE

These devices are used to assist in extraction of foreign bodies from the urinary tract. Foreign bodies may include fractured pieces of other medical devices such as a ureteral stent, catheter, or laser fiber.

CONTRAINDICATIONS

No contraindications.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- User should be familiar with and have experience in urological endoscopic surgery.
- Assess foreign bodies prior to instrument deployment to ensure that object is not too large to be removed through anatomy.
- If resistance is encountered while attempting to remove a foreign body, release the object. Do not exert excessive force on the device.
- This device is conductive. Avoid contact with any electrified instrument.
- Enclose the device in the sheath before removing from the tray/holder.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding e.g. Hematoma, Hematuria, Hemorrhage
- Death
- Deep venous thrombosis
- Fever
- Fluid overload
- Hypotension
- Hypothermia
- Inflammation e.g. Systemic inflammatory response syndrome
- Infection e.g. Bacteremia, Bacteriuria, Peri-renal Abscess, Pneumonia, Pyelonephritis, Sepsis, Urinary Tract Infection
- Nausea
- Pain
- Pulmonary Embolism
- Tissue, organ or structure injury e.g. Arteriovenous fistula, False tract, Infundibular stenosis, Pleural injury, Pseudoaneurysm, Renal collecting system injury, Selective arterial embolization, Surrounding organs/tissues injury, Urinary fistula

- Ureteral injury e.g. Ureteral avulsion, Ureteral mucosal injury, Ureteral mucosal tear, Ureteral obstruction, Ureteral perforation, Ureteral stricture
- Urinary tract issues e.g. Dysuria, Extrarenal stone migration, Fluid collection, Hydronephrosis, Steinstrasse, Urinary frequency, Urinary leakage, Urinary retention, Urinary urgency

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Move the thumb lever on the handle to the "Close" position to close the wire loop and remove the retriever from its packaging.
2. If included, the inserter may be placed securely into the Tuohy-Borst Adapter on the scope, as desired.
3. Then carefully insert the retracted retriever into the scope.
4. Under direct visualization, locate the target object and advance the retriever to a location just proximal to it. Slide the thumb lever on the handle to the "Open" position to deploy the wire loop near the target object and maneuver the device until it is within the wire loop.
5. Once correctly positioned, draw the wire loop down around the target object, and slide the thumb lever into the "Close" position.
6. Maintaining position of the target object within the wire loop, retract the scope and wire loop.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd