EN **NTRAP® Stone Entrapment and Extraction Device**

3 Instructions for Use

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NTRAP® STONE ENTRAPMENT AND EXTRACTION DEVICE

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The NTRAP® Stone Entrapment and Extraction Devices have a 2.8 Fr diameter and in a length of 145 cm. The braided nitinol wire extends from the cannula and forms the interwoven nitinol basket with a 7-millimeter diameter and 3.5 mm depth. A Uni-Torquer handle is located on the proximal end of the basket to facilitate opening and closing the basket. An inserter is also included.

Performance Characteristics

- Nitinol wire basket retains shape after multiple actuations and stone removals.
- Decreased material thickness at shaft end to increase flexibility and increase scope deflection and irrigation flow.
- Device shaft consisting of inner stainless-steel coil to increase torque response.
- Length of product allows targets in the kidney to be manipulated and removed.
- Stainless steel braid incorporated into sheath material for resistance to kinking.
- Includes inserter to assist insertion of the extractor into the working channel of the endoscope.

Device Compatibility

- The extractor is compatible with endoscope working channels with a minimum inner diameter of 3 Fr (1 mm) and a maximum length of 80 cm.
- The inserter tip allows for insertion of the extractor.
- These devices are compatible with sterile saline and sterile water.

Patient Population

This device is not intended to be used in any specific target population unless outlined on the labelling, in the instructions, and/or in promotional materials. Patient population could include any person requiring treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended User

This device is intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

This device interacts with the structures and tissues comprising the urinary tract, which may include the urethra, bladder, ureters, and kidneys.

Operating Principle

The NTrap Nitinol Stone Entrapment and Extraction Device has a basket manufactured from nitinol wire with a 2.8 Fr diameter shaft. The device is 145 cm in length, and a unitorquer handle facilitates maneuverability and the opening and closing of the basket. The device also includes an optional inserter that may be used to facilitate insertion of the extractor into the scope.

INTENDED USE

This device is used as an endoscopic entrapment and extraction device for calculi and other foreign bodies in the urinary tract, and to minimize stone migration during laser, ultrasonic, electrohydraulic or pneumatic lithotripsy.

DEVICE AND CLINICAL BENEFITS

To manipulate and remove stones from the urinary tract through the working channel of an endoscope.

INDICATIONS FOR USE

This device is used to assist in extraction of stones (urolithiasis) or other objects (foreign bodies) from the urinary tract. Foreign bodies can include fractured pieces of other medical devices such as a ureteral stent, catheter, or laser fiber.

CONTRAINDICATIONS

No contraindications.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- Users should be familiar with and experienced in urological endoscopic surgery.
- Assess stone or other foreign bodies prior to instrument deployment to ensure that object is not too large to be removed through the anatomy.
- If resistance is encountered while attempting to remove stone or other foreign bodies, release the object. Do not exert excessive force on the device.
- Due to the asymmetric nature of the NTRAP, do not torque or rotate the device in vivo.
- When used in conjunction with holmium: YAG laser, do not fire the laser directly upon the device.
- This device is conductive. Avoid contact with any electrified instrument.
- Enclose device in sheath before removing from tray/holder.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding e.g. Hematoma, Hematuria, Hemorrhage
- Death
- Deep venous thrombosis
- Fever
- Fluid overload
- Hypotension
- Hypothermia
- Inflammation e.g. Systemic inflammatory response syndrome
- Infection e.g. Bacteremia, Bacteriuria, Peri-renal Abscess, Pneumonia, Pyelonephritis, Sepsis, Urinary Tract Infection
- Nausea
- Pain
- Pulmonary Embolism
- Tissue, organ or structure injury e.g. Arteriovenous fistula, False tract, Infundibular stenosis, Pleural injury, Pseudoaneurysm, Renal collecting system injury, Selective arterial embolization, Surrounding organs/tissues injury, Urinary fistula
- Ureteral injury e.g. Ureteral avulsion, Ureteral mucosal injury, Ureteral mucosal tear, Ureteral obstruction, Ureteral perforation, Ureteral stricture
- Urinary tract issues e.g. Dysuria, Extrarenal stone migration, Fluid collection, Hydronephrosis, Steinstrasse, Urinary frequency, Urinary leakage, Urinary retention, Urinary urgency

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Retract the basket into the sheath, remove the stone extractor from its outer package, and place into the sterile field.
2. The inserter may be used to facilitate insertion of the extractor into the scope.
3. Then, carefully insert the retracted extractor into the scope.
4. Using direct vision, advance the device sheath beyond the stone or foreign object prior to deployment.
5. Deploy the device by advancing the handle forward while holding the sheath in position. The device may now be used to minimize stone migration during laser, ultrasonic, electrohydraulic or pneumatic lithotripsy.

NOTE: When using laser or electrohydraulic lithotripsy, a gap of at least 1 cm should be maintained between the stone and the NTRAP to prevent device damage.

NOTE: The NTRAP, in the open position, may be used as an extraction device to sweep the ureter, collecting stone fragments into the bladder.

6. It may be necessary to gently advance the basket the release any remaining stone fragments prior to closing.
7. Upon completion of the procedure, close the basket by pushing the sheath forward while holding the operating handle in position.

NOTE: The operating handle may be detached to facilitate scope removal.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2023
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd