EN
3

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NCIRCLE®, NCOMPASS®, AND NFORCE® STONE EXTRACTORS

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

NCircle® Nitinol Tipless Helical Stone Extractors, NCircle® Nitinol Tipless Stone Extractors, NCompass® Nitinol Stone Extractors, and The NForce® Nitinol Helical Stone Extractors are 3- or 4- wire nitinol baskets available in a range of French sizes, lengths, and basket diameters. An inserter and endoscopic cap or Tuohy-Borst Adapter may be included.

Performance Characteristics

- Nitinol wire basket retains shape after multiple actuations and stone removals.
- Device shaft consisting of inner stainless-steel coil to increase torque response.
- Length of product allows targets in the kidney to be manipulated and removed.
- Stainless steel braid incorporated into sheath material for resistance to kinking.
- Ergonomic handle design allows for use with one hand with labeled open and close directions.
- Includes Tuohy-Borst Adapter or endoscopic cap to resist backflow of fluid from the working channel of the endoscope around the extractor.
- Includes inserter to assist insertion of the extractor into the working channel of the endoscope.

Device Compatibility

- The 1.5, 1.7, 2.2, and 2.4 Fr devices are compatible with endoscope working channels with a minimum inner diameter of 3.0 Fr (1 mm) and a maximum length of 80 cm.
- The 3.0 and 3.2 Fr devices are compatible with endoscope working channels with a minimum inner diameter of 3.6 Fr (1.2 mm) and a maximum length of 80 cm for 115 cm length devices and 45 cm for 65 cm length devices.
- The 4.5 Fr devices are compatible with endoscope working channels with a minimum inner diameter of 5.1 Fr (1.7 mm) and a maximum length of 45 cm.
- The endoscopic cap or Tuohy-Borst Adapter may be included and can be connected to any standard female Luer lock fitting or bulb-type endoscopic fitting.
- The inserter tip is compatible with the endoscopic cap or Tuohy-Borst Adapter and allows for insertion of the extractor, if applicable.
- These devices are compatible with sterile saline and sterile water.

RPN	FR	Length (CM)	Basket Diameter / Configuration
NTSE-015115	1.5	115	1.0 cm
NCT4-017115	1.7	115	9 mm
NTSE-022115-UDH	2.2	115	1.0 cm
NFHSES-022115-3W-UDH	2.2	115	1.0 cm / 3 Wire
NTSE-022115-UDH-MB	2.2	115	1.8 cm
NCT4-024115	2.4	115	9 mm
NTSED-024115-UDH	2.4	115	1.0 cm

RPN	FR	Length (CM)	Basket Diameter / Configuration
NCT4-024115-MB	2.4	115	1.5 cm
NTSED-024115-UDH-MB	2.4	115	1.8 cm
NTHSES-030115-UDH, NTSE-030115-UDH	3	115	1.0 cm
NTSE-030115-UDH-MB	3	115	1.8 cm
NTSE-030065	3	65	1.0 cm
NFHSES-032115-3W-UDH	3.2	115	1.0 cm / 3 Wire
NTHSE-045065-UDH	4.5	65	1.5 cm
NTSE-045065-UDH	4.5	65	1.8 cm

Patient Population

These devices are not intended to be used in any specific target population unless outlined on the labelling, in the instructions, and/or in promotional materials. Patient population could include any person requiring treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues comprising the urinary tract, which may include the urethra, bladder, ureters, and kidneys.

Operating Principle

The NCompass® Nitinol Stone Extractor is a single-use device used for stone manipulation and removal in the urinary tract. At full deployment, the unique multi-wire geometry provides a traditional four-wire basket. Partial closure redesigns the basket to provide the weave of a 16-wire basket. The duality of the NCompass® extractor allows for both routine capture and extraction of intact stones as well as an ability to capture multiple small fragments following intracorporeal lithotripsy. These devices are manufactured in 1.7 and 2.4 Fr diameters and in a length of 115 cm.

The NForce® Nitinol Helical Stone Extractor is a single-use device used for stone manipulation and removal in the urinary tract. The basket is constructed with nitinol Delta Wire. The triangular shape of Delta Wire creates a three-wire basket with the largest possible wire mass. The greater mass of the Delta Wire combined with the inherent memory qualities of nitinol alloy, creates a basket with increased radial force which permits greater radial dilation in the ureter. These devices are manufactured in 2.2 and 3.2 Fr diameters and in a length of 115 cm.

The NCircle® Nitinol Tipless Stone Extractor and NCircle® Delta Wire Nitinol Tipless Stone Extractor is a single-use device used for stone manipulation and removal in the urinary tract. The nitinol memory wire allows the basket to retain its shape following extreme torsion. The unique tipless design aids in calyceal stone capture by allowing the basket to open at the base of the calyx and encircle the stone. The tipless design reduces the likelihood of mucosal damage and perforation. The flexible sheath permits continuous deflection of a flexible ureteroscope during stone manipulation. These devices are manufactured in 1.5, 2.2, 2.4, 3.0, and 4.5 Fr diameters and in lengths of 115 cm and 65 cm, respectively.

The NCircle® Nitinol Tipless Helical Stone Extractor is a single-use device that is used for stone manipulation and removal in the urinary tract. The nitinol memory wire allows the helical basket to retain its shape following extreme torsion. The unique

tipless design aids in calyceal stone capture by allowing the basket to open at the base of the calyx and encircle the stone reducing the likelihood of mucosal damage and perforation. The flexible sheath permits continuous deflection of a flexible ureteroscope during stone manipulation. The helical basket provides torque to aid in stone capture. These devices are manufactured in 3.0 and 4.5 Fr diameters and in lengths of 115 and 65 cm, respectively.

The Uni-Dex Handle of the devices facilitates ergonomic, one-handed use. Sliding the thumb lever in the "OPEN" direction opens the basket, while sliding the thumb lever in the "CLOSE" direction closes the basket. An optional Tuohy-Borst adapter or endoscopic cap may be provided that can be connected to the endoscope working channel to resist fluid backflow around the extractor. An optional inserter may also be provided that may be used to facilitate insertion of the extractor into the adapter or cap.

INTENDED USE

These devices are intended for stone manipulation and removal in the urinary tract through the working channel of an endoscope during cystoscopic or ureteroscopic procedures by actuating a nitinol wire basket.

DEVICE AND CLINICAL BENEFITS

To manipulate and remove stones from the urinary tract through the working channel of an endoscope.

INDICATIONS FOR USE

These devices are used to assist in extraction of stones (urolithiasis).

CONTRAINDICATIONS

No contraindications.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- These devices are designed for the removal of stones no larger than their maximum basket diameter according to the table located in the device compatibility section of this instruction.

PRECAUTIONS

- This device is conductive. Avoid contact with any electrified instrument.
- Enclose the device in the sheath before removing from the tray/holder.
- User should be familiar with and have experience in urological endoscopic surgery.
- Assess stone prior to instrument deployment to ensure that the stone is not too large to be removed through anatomy.
- If resistance is encountered while attempting to remove stone, release the stone. Do not exert excessive force on the device.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding e.g. Hematoma, Hematuria, Hemorrhage
- Death
- Deep venous thrombosis
- Fever
- Fluid overload
- Hypotension
- Hypothermia
- Inflammation e.g. Systemic inflammatory response syndrome
- Infection e.g. Bacteremia, Bacteriuria, Peri-renal Abscess, Pneumonia, Pyelonephritis, Sepsis, Urinary Tract Infection

- Nausea
- Pain
- Pulmonary Embolism
- Tissue, organ or structure injury e.g. Arteriovenous fistula, False tract, Infundibular stenosis, Pleural injury, Pseudoaneurysm, Renal collecting system injury, Selective arterial embolization, Surrounding organs/tissues injury, Urinary fistula
- Ureteral injury e.g. Ureteral avulsion, Ureteral mucosal injury, Ureteral mucosal tear, Ureteral obstruction, Ureteral perforation, Ureteral stricture
- Urinary tract issues e.g. Dysuria, Extrarenal stone migration, Fluid collection, Hydronephrosis, Steinstrasse, Urinary frequency, Urinary leakage, Urinary retention, Urinary urgency

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Move the thumb lever on the handle to the "Close" position to close the basket and remove the extractor from its packaging.
2. If included, the instrument insertion adapter may be placed securely into the endoscopic cap or Tuohy-Borst Adapter on the scope, as desired.
3. Then carefully insert the instrument into the scope.
4. Under direct visualization, locate the stone and advance the extractor to a location just proximal to the stone. Slide the thumb lever on the handle to the "Open" position to deploy the basket near the stone and maneuver the device until the stone is within the basket.
5. Once correctly positioned, draw the basket down around the stone and slide the thumb lever into the "Close" position.
6. Maintaining position of the stone within the basket, retract the scope and basket.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd