Cope Nephroureterostomy Stents

EN
3

Instructions for Use

Cope nefroureterostomi-stents

DA
5

Brugsanvisning

Cope Nephroureterostomie-Stents

DE
7

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ZH
25

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COPE NEPHROURETEROSTOMY STENTS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cope Nephroureterostomy Stent is constructed of Ultrathane® and comes in two variations. Products with an “-A” suffix in their order numbers have sideports in the loops **only**, and are used for bypassing lacerated ureters. Products with a “-B” suffix in their order numbers have sideports in the loops and in the entire distal segment, and are used for all other stenting purposes. Both “A” and “B” stents are available with one of two locking mechanisms: the Mac-Loc® Locking Loop or the Cook-Cope Type Locking Loop.

INTENDED USE

The Cope Nephroureterostomy Stent is intended for internal drainage from the ureteropelvic junction to the bladder while maintaining external access to the stent.

CONTRAINDICATIONS

None known

WARNINGS

Stents with Cook-Cope Type Locking Loop mechanisms contain **NATURAL RUBBER LATEX**. Possible allergic reactions should be considered.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in placement of nephroureterostomy stents. Standard techniques should be employed.
- Manipulation of product requires fluoroscopic control.
- A TFE-coated wire guide must be used with this product.
- Activate hydrophilic coating, if present, by wetting surface of device with sterile water or saline. For best results, maintain wetted condition of device during placement.

INSTRUCTIONS FOR USE

Stent Placement

1. Using standard percutaneous access technique, establish wire guide position well into the bladder.
2. Over the wire guide, introduce the stent/stiffening cannula into the kidney collecting system.
3. After establishing proper proximal and distal position, push the stent off the stiffening cannula over the wire, making sure the distal pigtail forms within the bladder.
4. Verify final position of the stent via fluoroscopy, then carefully withdraw the wire guide.

NOTE: Contrast medium in the collecting system can be diluted by injection of saline solution.

5. Form the proximal pigtail loop in the kidney using appropriate technique for the locking mechanism type as described below, then fix the catheter at the skin surface and apply dressings.

For Mac-Loc® Locking Loop Mechanism:


Fig. 1

- a. Stabilize the Mac-Loc hub assembly with one hand and pull back on the suture to form the distal stent loop configuration. (**Fig. 1**)


Fig. 2

- b. While maintaining traction on the suture, push the locking cam lever down until a distinct “snap” is felt. The distal loop of the stent is now locked into position. (**Fig. 2**)


Fig. 3

- c. Trim off the excess suture. (**Fig. 3**)

For Cook-Cope Type Locking Loop Mechanism:

- Pull suture end tight to form loop configuration in stent, and tie securely.
- Trim off excess suture, and slide latex sleeve over suture to prevent leakage.

Unlocking Catheter Loop

For Mac-Loc® Locking Loop Mechanism:

- While stabilizing the Mac-Loc hub assembly with one hand, position a small, blunt object (approximately the shape and size of a ball point pen or small forceps) into the Mac-Loc release notch.


Fig. 4

- Pry upward until the locking cam lever is free. (**Fig. 4**)
- Withdraw stent.

NOTE: For stent exchange, advance the distal end of a wire guide into the locked loop configuration of the stent before unlocking the Mac-Loc assembly. Release the Mac-Loc as described above. Advance the wire guide through the distal end of the stent until you have an adequate length of wire guide in the bladder. Stent exchange can now be performed.

For Cook-Cope Type Locking Loop Mechanism:

- Advance wire guide through stent, ensuring excess wire guide is positioned within the bladder.
- Uncover suture by folding back the latex sleeve, and release suture.
- Withdraw stent over wire guide.

NOTE: For stent exchange, insert new stent over wire guide. For removal, withdraw wire guide following removal of stent.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.