- EN**
3
Cope Nephroureterostomy Stents
Instructions for Use
- CS**
5
Nefroureterostomické stenty Cope
Návod k použití
- DA**
7
Cope nefroureterostomi-stents
Brugsanvisning
- DE**
10
Cope Nephroureterostomie-Stents
Gebrauchsanweisung
- EL**
12
Ενδοπροσθέσεις νεφροουρητηροστομίας Cope
Οδηγίες χρήσης
- ES**
15
Stents de nefroureterostomía Cope
Instrucciones de uso
- FR**
18
Endoprothèses pour néphro-urétérostomie de Cope
Mode d'emploi
- HU**
20
Cope nephroureterostomiás sztentek
Használati utasítás
- IT**
23
Stent per nefroureterostomia Cope
Istruzioni per l'uso
- NL**
25
Cope nefro-ureterostomiestents
Gebruiksaanwijzing
- NO**
28
Cope nefroureterostomistenter
Bruksanvisning
- PL**
30
Stenty do nefroureterostomii Cope
Instrukcja użycia
- PT**
33
Stents de nefroureterostomia Cope
Instruções de utilização
- SV**
35
Cope stentar för nefrouretärostomi
Bruksanvisning


COPE NEPHROURETEROSTOMY STENTS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions for use before using this device.

DEVICE DESCRIPTION

The Cope Nephroureterostomy Stent is constructed of Ultrathane® and comes in two variations. Products with an "A" suffix in their order numbers have sideports in the loops **only**, and are used for bypassing lacerated ureters. Products with a "B" suffix in their order numbers have sideports in the loops and in the entire distal segment, and are used for all other stenting purposes. Both "A" and "B" stents are available with the Mac-Loc® Locking Loop.

Cope Nephroureterostomy Stent

Flexible Stiffening Cannula

INTENDED USE

The Cope Nephroureterostomy Stent is delivered percutaneously and is intended to establish internal drainage from the ureteropelvic junction to the bladder while maintaining external access to the stent. The Cope Nephroureterostomy Stents are not intended to remain indwelling more than 90 days.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during catheterization.

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in placement of nephroureterostomy stents. Standard techniques should be employed.
- This product is intended for use in adult patients only.
- Manipulation of product requires fluoroscopic control.
- A PTFE-coated wire guide must be used with this product.
- Activate hydrophilic coating, if present, by wetting surface of device with sterile water or saline. For best results, maintain wetted condition of device during placement.
- Where long-term use is indicated, it is recommended that indwelling time not exceed 90 days and that the physician evaluate the catheter before this time has expired.

POTENTIAL ADVERSE EVENTS

The complications that may result from the use of these devices include, but are not limited to:

- Extravasation
- Urinary tract infection
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Abrasions
- Sepsis
- Hemorrhage
- Bowel injury
- Iatrogenic renal artery injury
- Pleural complications (including but not limited to Pneumothorax)
- Hematuria
- Fever
- Encrustation
- Stent dislodgement
- Stent obstruction
- Lower urinary tract symptoms
- Pain

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Cope Nephroureterostomy Stent is **MR Conditional**. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the Cope Nephroureterostomy Stent is expected to produce a maximum temperature rise of less than 2.5° C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the Cope Nephroureterostomy Stent extends approximately 46 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

INSTRUCTIONS FOR USE

Stent Placement

1. Using standard percutaneous access technique, establish wire guide position through the renal pelvis and ureter into the bladder using a support wire guide in a diameter appropriately sized to the endhole of the catheter.
2. Determine the appropriate catheter size using imaging techniques.
3. Flush the catheter with sterile normal saline or water.
4. Insert the stiffening cannula into the catheter and lock it in place.
5. Over the wire guide, introduce the stent/stiffening cannula into the kidney collecting system.
6. After establishing proper proximal and distal position, push the stent off the stiffening cannula over the wire, making sure the distal pigtail forms within the bladder.
7. Remove the stiffening cannula from the stent, leaving the wire guide in place.
8. Verify final position of the stent via fluoroscopy, then carefully withdraw the wire guide.

NOTE: Contrast medium in the collecting system can be diluted by injection of saline solution.

9. Form the proximal pigtail loop in the kidney using the appropriate technique for the locking mechanism type, as described below.
10. Fix the catheter at the skin surface and apply dressings.
11. The catheter may be capped for internal drainage or may be connected to an external drainage bag.
12. Monitor urinary output and perform clinical evaluation for any abnormal decrease in output.

For Mac-Loc® Locking Loop Mechanism:


Fig. 1

- a. Stabilize the Mac-Loc hub assembly with one hand and pull back on the suture to form the distal stent loop configuration. (**Fig. 1**)


Fig. 2

- b. While maintaining traction on the suture, push the locking cam lever down until a distinct "snap" is felt. The distal loop of the stent is now locked into position. (**Fig. 2**)


Fig. 3

- c. Trim off the excess suture. (**Fig. 3**)

Unlocking Catheter Loop

- While stabilizing the Mac-Loc hub assembly with one hand, position a small, blunt object (approximately the shape and size of a ball point pen or small forceps) into the Mac-Loc release notch.


Fig. 4

- Pry upward until the locking cam lever is free. (Fig. 4)
- Withdraw stent.

NOTE: For stent exchange, advance the distal end of a wire guide into the locked loop configuration of the stent before unlocking the Mac-Loc assembly. Release the Mac-Loc as described above. Advance the wire guide through the distal end of the stent until you have an adequate length of wire guide in the bladder. Stent exchange can now be performed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd