## Park Blade Septostomy Catheter

**EN**  
**3**

Instructions for Use

## Park kateter til bladseptostomi

**DA**  
**5**

Brugsanvisning

## Park-Klingenseptostomiekatheter

**DE**  
**7**

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A. Borgschroef  
B. Afdichting

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## PARK BLADE SEPTOSTOMY CATHETER

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner.

### DEVICE DESCRIPTION

The Park Blade Septostomy Catheter is a radiopaque polyethylene catheter with a distal tip section made of stainless steel tubing. The metal tubing contains a blade that is linked to a lever, the distal portion of which pivots at the tip. The proximal portion of the blade is linked to a solid wire guide that passes through the entire catheter and exits the hub. Advancing the wire guide extends the blade and the lever through a slit in the metal tubing to form a triangle.

### INTENDED USE

The Park Blade Septostomy Catheter is intended to enlarge interatrial openings. The device can also be used when balloon atrial septostomy is insufficient or unsuccessful, particularly in older infants and children with a thickened atrial septum. Normally, the existing interatrial opening is used for this procedure. However, if the interatrial septum is intact, the procedure can be performed in conjunction with a transseptal technique.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

### WARNINGS

Under no circumstances should rapid catheter withdrawal through the septum be attempted, as is required for balloon septostomy. A gentle, but firm, force should be maintained while withdrawing the catheter from the left atrium to the right atrium. Discontinue traction when a sudden decrease in resistance is felt.

### PRECAUTIONS

- Do not use this product if the left atrium is smaller than 2 cm in diameter.
- Assess to determine patient does not have bilaterally thrombosed femoral or iliac veins.
- Assess to determine patient does not have an absent hepatic segment of the inferior vena cava.
- Assess to determine patient does not have other conditions that prohibit femoral venous approach.
- Manipulation of product requires biplane fluoroscopic control.
- Do not inject contrast medium through the Park Blade Septostomy Catheter. Contrast medium can cause the blade mechanism to malfunction.
- In the case of an unusually thick interatrial septum where the interatrial opening is very small, or when the transseptal approach is used for an intact atrial septum, the first withdrawal of the Park Blade Septostomy Catheter must be conducted with the blade only partially extended to prevent undue stress to the blade lever component and linkage.
- If resistance is met or the blade cannot be fully extended, check to ensure that the catheter tip is not positioned in the left atrial appendage or in the pulmonary vein. If so, pull back the pin vise to collapse the blade into the catheter, then withdraw the entire catheter system slightly and repeat extension of the blade.

### PRODUCT RECOMMENDATIONS

#### Park Blade Septostomy Sizing Recommendations

| Size of Interatrial Opening | Patient Age    | Patient Weight | Suggested Blade Catheter |
|-----------------------------|----------------|----------------|--------------------------|
| <5 mm                       | Under 5 months | <4 kg          | PBS-100 (small)          |
| 5-8 mm                      | Young children | 4-10 kg        | PBS-200 (medium)         |
| >8 mm                       | Older children | >10 kg         | PBS-300 (large)          |

#### Introducer Sheath Size Recommendations

| Order Number | French Size | Length | Blade Length | Recommended Sheath Size |
|--------------|-------------|--------|--------------|-------------------------|
| PBS-100      | 5.7         | 65 cm  | 9.4 mm       | 7 Fr                    |
| PBS-200      | 5.7         | 65 cm  | 13.4 mm      | 7 Fr                    |
| PBS-300      | 7.3         | 85 cm  | 20 mm        | 9 Fr                    |

### INSTRUCTIONS FOR USE

#### Device Preparation

1. Test the Park Blade Septostomy Catheter before use by advancing and retracting the pin vise to ensure that the blade mechanism works properly. (**Fig. 1**)

2. Pull the pin vise to fold the blade into the slit of the stainless steel tubing.
3. Loosen the pin vise and slide it forward until there is no gap between the pin vise and the Y-adapter attached to the catheter hub.
4. Tighten the pin vise, ensuring no gap remains between the Y-adapter and the pin vise, to prevent inadvertent protrusion of the blade during catheter manipulation.

### **Determining Size of Interatrial Opening**

5. Using standard technique, introduce an angled introducer sheath (e.g., Mullins-type) into the right atrium. Refer to the Product Recommendations Section of these Instructions for Use for introducer sheath size recommendations.
6. Advance a balloon septostomy catheter through the introducer sheath into the left atrium (via the right atrium).
7. Inflate the balloon to a diameter larger than the interatrial opening and gently pull it against the interatrial septum. Gradually deflate the balloon, while maintaining traction, until it passes through the defect.
8. Use the size of the balloon as it passes across the atrial septum, viewed fluoroscopically, to determine the true size of the interatrial opening.
9. Withdraw the balloon catheter.

### **Septostomy Procedure**

10. Position the introducer sheath in the left atrium and perform pressure measurements. **NOTE:** It is recommended that left atrial angiocardiology be performed to evaluate the size of the left atrium prior to blade septostomy. The left atrium must be larger than 2 cm in diameter.
11. Insert the appropriately-sized Park Blade Septostomy Catheter through the introducer sheath into the left atrium. Refer to catheter sizing recommendations in the Product Recommendations Section of these Instructions for Use.
12. Confirm fluoroscopically the position of the Park Blade Septostomy Catheter tip in the left atrium. Maintaining catheter position, withdraw the tip of the introducer sheath to the level of the mid-inferior vena cava.
13. Loosen the pin vise and slide it backward until its distal edge is aligned with the hash mark on the wire guide then tighten the pin vise. **CAUTION:** Once the pin vise is tightened, advancing it will cause the blade to extend. Use caution to ensure that you do not inadvertently advance the pin vise prematurely.
14. Under fluoroscopic control, advance the pin vise gently toward the Y-adapter to extend the blade.  
**WARNING: If resistance is met or if the blade cannot be fully extended, check to ensure the catheter tip is not positioned in the left atrial appendage or in the pulmonary vein. If so, pull back the pin vise to collapse the blade into the catheter, then withdraw the entire catheter system slightly and repeat extension of blade.**
15. Once the blade is fully extended, hold the Y-adapter and pin vise together with the thumb and index finger and slowly rotate the catheter as needed so that the blade faces anteriorly. Confirm anterior orientation of the blade under lateral fluoroscopy.
16. Using both hands to maintain anterior catheter orientation, slowly withdraw the entire catheter system (including introducer sheath) to the right atrium. Once resistance is met at the atrial septum, maintain gentle but firm traction until a sudden decrease in resistance is felt. **CAUTION:** The introducer sheath must be withdrawn simultaneously with the septostomy catheter to prevent possible cutting of the distal tip of the sheath.  
**WARNING: Under no circumstances should rapid catheter withdrawal through the septum be attempted, as is required for balloon septostomy.**
17. Repeat the procedure, with slight rotation of the blade between each pass until resistance is no longer felt upon withdrawal. To establish a satisfactory opening, two or three passages of the blade across the interatrial septum are typically sufficient.
18. A balloon atrial septostomy may now be performed through the introducer sheath to further enlarge the interatrial opening. **CAUTION:** The introducer sheath must be withdrawn simultaneously with the balloon catheter during the procedure to prevent damage to the sheath.

### **HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.