EN 2	Percutaneous Entry Set Instructions for Use
CS 3	Souprava pro perkutánní vstup Návod k použití
DA 4	Perkutant punktursæt Brugsanvisning
DE 5	Perkutaner Zugangs-Set Gebrauchsanweisung
EL 7	Σετ διαδερμικής εισόδου Οδηγίες χρήσης
ES 8	Equipo de acceso percutáneo Instrucciones de uso
FR 9	Set de ponction percutanée Mode d'emploi
HU 11	Perkután bemeneti készlet Használati utasítás
IT 12	Set per accesso percutaneo Istruzioni per l'uso
NL 13	Percutane-toegangsset Gebruiksaanwijzing
NO 14	Perkutant inngangssett Bruksanvisning
PL 16	Zestaw dojścia przezskórnego Instrukcja użycia
PT 17	Conjunto de acesso percutâneo Instruções de utilização
SV 18	Perkutant införingsset Bruksanvisning

PERCUTANEOUS ENTRY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using device.

DEVICE DESCRIPTION

Set contains:

- 22 gage Skinny Needle® with Chiba Tip, 20 cm long
- 18 gage Disposable Two-Part Trocar needle, 20 cm long
- .038 (0.97mm) PTFE-coated, stainless steel Safe-T-J® wire guide, 80 cm long
- Radiopaque dilators 6 through 16 French, 20 cm long

INTENDED USE

The Percutaneous Entry Set is used to establish percutaneous tract into the renal pelvis for catheter placement or stone manipulation.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

WARNINGS

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for percutaneous access should be employed.
- If resistance is encountered during advancement or withdrawal, stop. Do not continue without first determining cause of resistance and taking remedial action.
- If urine cannot be confirmed, stop, determine fluoroscopic position and reintroduce the access needle if needed.

POTENTIAL ADVERSE EVENTS

- Renal perforation
- Tissue Injury
- Hemorrhage
- Edema

INSTRUCTIONS FOR USE

1. Perform fluoroscopy or ultrasound to provide visualization of the target.
2. After prepping the skin in the usual manner, using fluoroscopy and/or ultrasound, make a needle puncture with the 22 gage needle to provide direct access to the intended area of interest.
3. Perform an antegrade pyelogram to provide visualization of the target collecting system, to gauge depth, and to distend the system for subsequent puncture.

NOTE: A urine specimen may be taken directly from the renal pelvis before contrast injection if desired.

4. If the needle position is suitable for nephrostomy, move on to step 5. The needle is withdrawn by grasping the needle at the skin so the depth of the collecting system can be gauged.

NOTE: In the case of small, partially obstructed systems, the needle can be left in place and contrast infused during the next step.

5. Using fluoroscopy and/or ultrasound, make a needle puncture with the 18 gage needle at the depth previously gauged, providing direct access to the intended area of interest. Free flow of urine indicates proper positioning. **If desired, an antegrade pyelogram can be performed through the 18 gage needle to provide visualization of the target collecting system.**

6. The wire guide is placed through the cannula of the 18 gage needle to achieve access into the collecting system and allow placement of subsequent devices.
7. While maintaining wire guide position, withdraw the needle.
8. The tract can now be sequentially dilated by passing the dilators over the wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd