- EN**
3 **Percutaneous Neonatal Pigtail Nephrostomy Set**
Instructions For Use
- CS**
5 **Perkutánní neonatální nefrostomická pigtailová souprava**
Návod k použití
- DA**
7 **Perkutant neonatalt nefrostomisæt med grisehale**
Brugsanvisning
- DE**
9 **Perkutanes neonatales Pigtail-Nephrostomie-Set**
Gebrauchsanweisung
- EL**
12 **ΣΕΤ διαδερμικής νεφροστομίας νεογνών με σπειροειδές άκρο**
Οδηγίες χρήσης
- ES**
14 **Equipo pigtail de nefrostomía percutánea neonatal**
Instrucciones de uso
- FR**
16 **Set pigtail pour néphrostomie percutanée chez les nouveau-nés**
Mode d'emploi
- HU**
19 **Perkután újszülöttgyógyászati pigtail nephrostomiás készlet**
Használati utasítás
- IT**
21 **Set con catetere pigtail neonatale per nefrostomia percutanea**
Istruzioni per l'uso
- NL**
23 **Percutane neonatale pigtail-nefrostomieset**
Gebruiksaanwijzing
- NO**
26 **Perkutant neonatalt pigtail-nefrostomiset**
Bruksanvisning
- PL**
28 **Zestaw noworodkowego przezskórnego cewnika nefrostomijnego typu „pigtail”**
Instrukcja użycia
- PT**
30 **Conjunto de cateter pigtail para nefrostomia percutânea neonatal**
Instruções de utilização
- SV**
32 **Set för perkutan neonatal pigtailnefrostomi**
Bruksanvisning

T _ P N N S _ R E V O

PERCUTANEOUS NEONATAL PIGTAIL NEPHROSTOMY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions prior to using this device.

DEVICE DESCRIPTION

- 22 gauge Chiba needle, 5 cm long
- 18 gauge trocar needle, 5 cm long
- .038 inch (0.97 mm) diameter PTFE-coated stainless steel wire guide, 30 cm long
- Radiopaque dilators, (5, 6, and 7 French), 10 cm long
- 6 French, 7 cm long radiopaque pigtail catheter
- Retention disc with pull tie
- 10 French clear polyvinylchloride connecting tube, 10 cm long with one-way stopcock

INTENDED USE

This device is intended for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage in neonatal patients.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk for nephrostomy drainage.

WARNINGS

Sterilized by ethylene oxide gas. Do not use if sterile barrier is damaged. This device is intended for one-time use.

PRECAUTIONS

- The maximum indwelling time for this device is four weeks. Periodic evaluation of the catheter is advised.
- This device is intended for use by physicians trained and experienced in percutaneous access and placement of drainage catheters. Standard techniques should be employed.
- Do not forcefully advance any component during placement or removal of the catheter. Carefully remove the components if any resistance is encountered.
- Decompression of the renal pelvis is recommended to alleviate overdistension of collecting system.
- Percutaneous access should be achieved from below the 12th rib to reduce risk of pneumothorax.

POTENTIAL ADVERSE EVENTS:

- Urosepsis or urinary tract infection
- Bleeding, hematuria, or vascular injury
- Pelvis/renal injury
- Catheter occlusion
- Catheter displacement
- Urinoma
- Urinary leakage or site oozing
- Loss of initial access (due to mobility of the anatomy in young children and infants)
- Injury to spleen, liver, or bowel
- Perforation of the collecting system
- Abdominal or flank pain

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated the Percutaneous Neonatal Pigtail Nephrostomy Catheter is **MR Conditional**. A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Percutaneous Neonatal Pigtail Nephrostomy Catheter is not expected to induce RF heating greater than that of biological tissue. The image artifact extends approximately 2.1 mm from the Percutaneous Neonatal Pigtail Nephrostomy Catheter as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 tesla MRI system.

INSTRUCTIONS FOR USE

1. By preliminary plain film, I.V.P., retrograde pyelogram, ultrasound or CAT scan, localize the kidney to be drained. A prone or slightly oblique position is preferred. Prepare and drape the flank in the usual fashion. Under fluoroscopic guidance and/or ultrasound, identify and anesthetize the skin site overlying the collecting system.
2. Introduce the 22 gauge (black hub) needle into the collecting system at a right angle to the skin. Obtain a urine specimen, if desired. Visualization of the pyelocalyceal system may be further enhanced by the injection of contrast medium.

NOTE: If desired, maintain the 22 gauge needle at this location and position throughout the remainder of the procedure. This fixes the kidney in position and allows further injection of the pyelocalyceal system.

3. Rotate the patient into a more lateral position, with the posterior axillary line of the ipsilateral flank in the superior position. Pass the 18 gauge (pink hub) needle vertically through a small incision in the skin (made with a scalpel blade); insert the needle into the appropriate part of the pyelocalyceal system. Free flow of fluid from the needle after removing the obturator confirms a satisfactory location for the tip of the needle cannula.
4. After confirming access with the 18 gauge needle (pink hub) into the collecting system, the 22 gauge needle (black hub) can be removed.
5. Pass the flexible end of the wire guide through the hub of the 18 gauge needle cannula and advance it into the collecting system. Confirm the position fluoroscopically. Maintain the position of the wire guide as the needle cannula is withdrawn over it. **NOTE:** Care must be exercised to avoid withdrawal of the wire guide itself.
6. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided, advancing the dilator over the wire guide.
7. Pass the pigtail end of the catheter over the external end of the wire guide; gradually advance the pigtail end well into the collecting system. Confirm the position fluoroscopically.
8. Holding the shaft of the pigtail catheter securely in position, withdraw the wire guide, to allow the retentive coil to form in the renal pelvis. Confirm the appropriate position of the catheter via fluoroscopic visualization. Then, fully remove the wire guide.

NOTE: To reposition the catheter, advance the wire guide back through the catheter and into the renal pelvis to straighten the catheter pigtail. The catheter can then be repositioned to the appropriate location. Confirm the appropriate position of the catheter via fluoroscopic visualization.

9. Tape or suture the retention disc to the skin; secure the catheter to the disc with the plastic tie. A dressing may be applied.
10. Use the connecting tube to connect the pigtail catheter to a drainage bag. The one-way stopcock is useful for testing the drainage capability of the collecting system, if desired.

Catheter Removal

1. Pass the flexible end of the wire guide through the pigtail catheter and into the renal pelvis.
2. While maintaining wire access in the renal pelvis, remove the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd