- EN**
3 **Pediatric Nephrostomy Stent Set**
Instructions For Use
- CS**
4 **Sada pediatrického nefrostomického stentu**
Návod k použití
- DA**
6 **Pædiatrisk nefrostomistsæt**
Brugsanvisning
- DE**
8 **Pädiatrisches Nephrostomie-Stent-Set**
Gebrauchsanweisung
- EL**
9 **Παιδιατρικό σετ ενδοπρόσθεσης νεφροστομίας**
Οδηγίες χρήσης
- ES**
11 **Equipo de stent de nefrostomía pediátrico**
Instrucciones de uso
- FR**
13 **Set d'endoprothèse pour néphrostomie pédiatrique**
Mode d'emploi
- HU**
15 **Gyermekgyógyászati nefrostomiás sztentkészlet**
Használati utasítás
- IT**
17 **Set con stent per nefrostomia pediatrico**
Istruzioni per l'uso
- NL**
18 **Pediatische nefrostomiestentset**
Gebruiksaanwijzing
- NO**
20 **Pediatrisk nefrostomistentsett**
Bruksanvisning
- PL**
22 **Zestaw pediatrycznego stentu nefrostomijnego**
Instrukcja użycia
- PT**
24 **Conjunto de stent para nefrostomia pediátrica**
Instruções de utilização
- SV**
26 **Stentset för pediatrisk nefrostomi**
Bruksanvisning

T _ P N S S _ R E V O

PEDIATRIC NEPHROSTOMY STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions prior to using this device.

DEVICE DESCRIPTION

- Radiopaque pigtail catheter and stent
- Catheter support with pull tie
- 10 French polyvinylchloride connecting tube, 30 cm long with one-way stopcock

*Set components may vary

NOTE: Catheter and stent accept up to an 0.038 inch diameter wire guide.

INTENDED USE

This device is used as a nephrostomy drainage catheter and ureteral stent in pediatric patients 2 years and older.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk for nephrostomy drainage and/or use of a ureteral stent.

WARNINGS

Sterilized by ethylene oxide gas. Do not use if sterile barrier is damaged. This device is intended for one-time use.

PRECAUTIONS

- The maximum indwelling time for this device is four weeks. Periodic evaluation of the catheter is advised.
- This device is intended for use by physicians trained and experienced in percutaneous access and placement of drainage catheters. Standard techniques should be employed.
- Do not forcefully advance any component during placement or removal of the catheter. Carefully remove the components if any resistance is encountered.
- Decompression of the renal pelvis is recommended to alleviate overdistension of collecting system.

POTENTIAL ADVERSE EVENTS:

- Urosepsis or urinary tract infection
- Bleeding, hematuria, or vascular injury
- Pelvis/renal injury
- Catheter occlusion
- Catheter displacement
- Urinoma
- Urinary leakage or site oozing
- Loss of initial access (due to mobility of the anatomy in young children and infants)
- Abdominal or flank pain

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated the Pediatric Nephrostomy Stent is **MR Conditional**.

A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only

- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Pediatric Nephrostomy Stent is not expected to induce RF heating greater than that of biological tissue. The image artifact extends approximately 2.2 mm from the Pediatric Nephrostomy Stent as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 tesla MRI system.

INSTRUCTIONS FOR USE

1. Perform standard techniques for establishing a nephrostomy tract under fluoroscopic guidance and/or ultrasound, leaving an 0.038 wire guide in the ureter.
2. Pass the stent end of the catheter over the external end of the wire guide. Gradually advance the catheter end well into the collecting system and place the stent segment into the ureter at the desired length.
3. While holding the shaft of the catheter securely in position, slightly withdraw the wire guide enough to allow the retentive coil to form in the renal pelvis. When you have confirmed appropriate position of the catheter via fluoroscopic visualization, fully remove the wire guide.
NOTE: To reposition the catheter, advance the wire guide back through the catheter and into the renal pelvis to straighten the catheter tip. The catheter can then be repositioned to the appropriate location. Confirm the appropriate position of the catheter via fluoroscopic visualization.
4. Tape or suture the retention disc to the skin; secure the catheter to the disc with the plastic tie. A dressing may be applied.
5. Use the connecting tube to connect the pigtail catheter to a drainage bag. The one-way stopcock is useful for testing the drainage capability of the collecting system, if desired.

Catheter Removal

1. Pass a wire guide through the pigtail catheter and stent portion, and down into the bladder.
2. While maintaining wire access, remove the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd.