EN
3

Advance® CS Coronary Sinus Infusion Catheter

Instructions for Use

DA
6

Advance® CS infusionskateter til sinus coronarii

Brugsanvisning

DE
9

Advance® CS Koronarsinus-Infusionskatheter

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Bruksanvisning

ADVANCE® CS CORONARY SINUS INFUSION CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Advance CS Coronary Sinus Infusion Catheter is a double lumen catheter with a balloon near its distal tip. The catheter consists of two independent lumens, which are labeled "DISTAL" and "BALLOON." The distal lumen extends the length of the catheter and is used for the placement of wire guides and/or compatible devices, and infusion. The distal lumen diameter is 0.038 inch (0.97 mm) and is compatible with 0.035 inch (0.89 mm) wire guides. The balloon lumen is used to expand the balloon. Inscribed on the tip of the balloon manifold are the balloon diameter (mm) and the balloon length (cm). Platinum/iridium, radiopaque markers are positioned on the shaft within the balloon to enable visualization of the catheter/balloon under fluoroscopy. Refer to the label for further information.

INTENDED USE

The Advance CS Coronary Sinus Infusion Catheter is intended for temporary occlusion of the coronary sinus for infusion of contrast media, drugs, or therapeutic agents; or possible introduction of devices into the coronary venous system.

PRODUCT RECOMMENDATIONS

Figure 1: Maximum Pressures

	Maximum
Dynamic Pressure	300 psi (20.7 bar)
Static Pressure	100 psi (6.9 bar)

CONTRAINDICATIONS

None known

WARNINGS

- Ensure balloon diameter is appropriate for placement in the intended segment of the coronary sinus based on vessel diameter and anatomy. Do not oversize the balloon. Oversizing the balloon may cause vessel damage or rupture.
- Use the minimum level of inflation required to occlude the vessel. The maximum balloon inflation pressure recommended in the coronary sinus is 2 atm. Use of the Advance CS Catheter in the coronary sinus at balloon pressures above 2 atm has not been evaluated, and may cause vessel damage or rupture.
- Do not exceed rated burst pressure. Rupture of the balloon may occur. Adhere to the balloon inflation parameters in **Figure 2**. Over-inflation may cause rupture of the balloon or damage to the vessel wall. Use of a pressure gauge is recommended to monitor inflation pressures.
- Do not use a power injector for balloon inflation through the balloon lumen. Rupture may occur.
- Do not exceed the distal lumen maximum static or dynamic pressures. Exceeding indicated pressures may result in catheter damage or rupture. Do not use a power injector for injection through the distal lumen.
- If distal lumen occlusion is suspected, do not flush the lumen. Over-pressurization may occur, resulting in catheter damage or rupture.
- To avoid vessel damage, ensure that the catheter tip is centered within the lumen of the coronary sinus and is not in contact with the vessel wall.
- Advancement of the Advance CS Coronary Sinus Infusion Catheter into the distal portion of the coronary sinus may increase the likelihood of vessel damage or rupture due to vessel anatomy.
- Infuse solutions slowly through the lumen into the coronary sinus to minimize the likelihood of vessel injury.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Prior to use, read all package inserts, warnings, precautions, and instructions.
- The Advance CS Coronary Sinus Infusion Catheter has been designed for introduction into the vascular system using percutaneous (Seldinger) technique. It is recommended that an introducer sheath or guiding sheath be used to facilitate placement. Refer to label information for the appropriate introducer sheath size.

- The balloon is constructed of heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- Adhere to the infusion pressure parameters indicated in **Figure 1**.
- The balloon is manufactured from an extra-thinwall, high-strength, minimally compliant material. Particular care should be taken in handling the balloon to prevent damage.
- In heavily scarred access sites, use of an introducer sheath is recommended.
- Manipulate the catheter using fluoroscopic control.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the expiration date specified on the package.

POTENTIAL ADVERSE EVENTS

- Air embolism
- Aneurysm
- Arrhythmias
- Arteriovenous fistula
- Cardiac vessel related complications (pericardial effusion, cardiac tamponade)
- Death
- Drug reactions, allergic reaction to contrast medium
- Hemorrhage or hematoma
- Hypo/hypertension
- Infection and pain at insertion site
- Systemic embolization
- Vascular thrombosis
- Vessel dissection, perforation, rupture, or injury
- Vessel spasm

INSTRUCTIONS FOR USE

Balloon Preparation

1. Choose an appropriate balloon size based on vessel diameter and anatomy.
2. Upon removal from package, inspect the catheter to ensure no damage has occurred during shipping.
3. Remove the protective sleeve from the balloon. Do not discard the protective sleeve after initial introduction of the catheter into the vascular system; it may be helpful in wrapping the balloon prior to subsequent insertions.
4. Prepare the balloon lumen with standard contrast-saline mixture as follows:
 - a) Fill a syringe of appropriate size with 1:1 mixture of contrast medium and normal saline.
 - b) Attach the syringe to the hub of the catheter marked "BALLOON" and pull back to apply negative pressure.
 - c) Release pressure, allowing negative pressure to draw mixture into balloon lumen.
 - d) Detach the syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
 - e) Prepare the inflation device in standard manner and purge to remove all air from the syringe and tubing.
 - f) Attach the inflation device to the balloon lumen, ensuring no air bubbles remain at connection.
 - g) Pull negative pressure on the inflation device.

Catheter Introduction, Balloon Inflation, and Infusion

5. Flush the catheter lumen labeled "DISTAL" using heparinized saline solution.
6. Apply negative pressure to the lumen labeled "BALLOON" prior to introduction. Advance the catheter while rotating counter-clockwise over a pre-positioned 0.035 inch (0.89 mm) wire guide.
7. Under fluoroscopy, advance the balloon into the coronary sinus. Carefully position the balloon using both the distal and proximal radiopaque balloon markers. **NOTE:** If resistance is met while advancing the catheter, determine the cause and proceed with caution.
8. Inflate the balloon to the desired pressure. Adhere to the recommended balloon inflation pressures (See **Figure 2**). **NOTE:** To avoid vessel damage, ensure that the catheter tip is centered within the lumen of the coronary sinus and is not in contact with the vessel wall.

Figure 2: Balloon Inflation Pressures Chart

Pressure (ATM) (BAR)	Balloon Diameter (mm)				
	5	6	7	8	9
1	5.22	5.94	6.79	7.57	8.54
2	5.30	6.08	6.97	7.77	8.78
3	5.44	6.32	7.25	8.09	9.17
4	X	X	X	X	X
5	X	X	X	X	X
6	X	X	X	X	X
7	X	X	X	X	X
8	X	X	X	X	X
9	X	X	X	X	X
10	X	X	X	X	X
11	X	X	X	X	X
12	X	X			
13	X				
14	X				
15	X				

 = Nominal and Maximum Recommended Balloon Pressure

 = Rated Burst Pressure

- Carefully withdraw the wire guide.
- Perform angiography by infusing contrast through the distal lumen of the infusion catheter to confirm coronary sinus occlusion. Do not exceed the maximum dynamic or static distal lumen pressures (**Figure 1**). **NOTE:** Repositioning of the balloon and/or modifying the balloon inflation pressure may be used to achieve satisfactory occlusion. Prior to repositioning the balloon or modifying the balloon pressure, the wire guide must be reinserted.
- Attach the syringe containing the infusion solution to the distal lumen port and infuse the solution into the coronary venous circulation. Do not exceed the maximum dynamic or static distal lumen pressures (**Figure 1**). To avoid vessel injury, infuse the solution slowly through the lumen into the coronary sinus. **NOTE:** If desired, 2 mL of saline may be used to completely flush the distal lumen before or after infusion.
- If balloon pressure is lost and/or balloon rupture occurs, deflate the balloon and remove the balloon and sheath as a unit.

Balloon Deflation and Withdrawal

- After infusion is complete, remove the syringe and advance the wire guide back through the distal lumen.
- Completely deflate the balloon by pulling vacuum on the inflation device or syringe. Allow adequate time for the balloon to deflate. **NOTE:** Balloons with large diameters and/or longer lengths may require longer deflation times.
- Maintain vacuum on the balloon and withdraw the catheter. Upon catheter withdrawal, a gentle counter-clockwise rotation of the catheter will assist in balloon rewrap, minimizing trauma to the percutaneous entry site.
- If resistance is met during withdrawal, apply negative pressure with a larger syringe before proceeding. If resistance continues, remove the balloon and sheath as unit.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages and intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, and cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.