EN 3	PTFE Wire Guides Instructions for Use
BG 5	PTFE телени водачи Инструкции за употреба
CS 8	Vodicí dráty PTFE Návod k použití
DA 10	PTFE-kateterledere Brugsanvisning
DE 12	PTFE-Führungsdrähte Gebrauchsanweisung
ES 15	Guías de PTFE Instrucciones de uso
ET 18	PTFE-juhtetraadid Kasutusjuhend
FI 20	PTFE-johtimet Käyttöohjeet
FR 22	Guides en PTFE Mode d'emploi
HR 25	Žice vodilice od PTFE-a Upute za uporabu
HU 27	PTFE vezetődrótok Használati utasítás
IT 30	Guide in PTFE Istruzioni per l'uso
LV 33	PTFE vaditājstīgas Lietošanas instrukcija

NL
35

PTFE-voerdraden

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PL
37

Prowadniki powlekane PTFE

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SV
40

PTFE-ledare

Bruksanvisning

PTFE WIRE GUIDES

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Bentson, Fixed Core, and Fixed Core Curved PTFE Wire Guides are manufactured as fixed core wire guides with polytetrafluoroethylene (PTFE) coating on the outer coils, and may be manufactured with a flexible straight or curved distal tip. These devices have a tapered mandril wire, coil, and a safety wire. These devices are manufactured in a length of 145 or 80 cm and diameters of 0.035 or 0.038 inches.

Performance Characteristics

- Stainless steel coiling material for resistance to kinking.
- These devices identified by part numbers "638854" have a curved proximal end.

Device Compatibility

These devices are compatible with appropriately sized needles, sheaths, dilators, endoscopes, stents and catheters.

Patient Population

These devices are used for adult patients with a urological condition who require a procedure using a wire guide.

Intended Users

These devices are intended for use by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues of the upper and lower urinary tract.

Operating Principle

These devices are inserted into the body through the lumen of a percutaneous entry needle, guiding sheath, or other interventional device. Only the flexible tip is intended as the initially inserted end. A wire guide inserter may be employed to aid in wire guide insertion (e.g. by straightening out a curved tip). After insertion, medical imaging is used to visualize the movement of the wire guide to the desired location within the body. Catheters or other devices can then be inserted over the wire guide and directed to the intended location.

INTENDED USE

These devices are stainless steel fixed core wire guides that assist the insertion, removal and positioning of other devices within the urinary tract. These devices act as a guide to access a portion of the urinary tract during a urological procedure to allow placement, replacement, or exchange of medical devices.

INDICATIONS FOR USE

Multiple medical conditions may result in the need for a urological procedure. Urology wire guides are used as a guide to access a portion of the urinary tract to allow placement, replacement, or exchange of medical devices during urological procedures such as percutaneous nephrostomy (PCN), percutaneous nephrolithotomy (PNL), ureteroscopy (URS), and stent placement or exchange.

DEVICE AND CLINICAL BENEFITS

To assist with urinary tract access for minimally invasive placement, replacement, or exchange of medical devices (without the wire guide itself having a direct diagnostic or therapeutic function).

CONTRAINDICATIONS

None

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Do not force the device during the procedure; proceed with caution if excessive resistance is encountered. Excessive force can lead to the wire guide fracturing, bending, kinking, unraveling or uncoiling. Additionally, excessive force can cause the wire guide to become entrapped.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1%w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement or migration.
- Manipulation of wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access. Over manipulation can cause the wire guide to become entrapped.
- When using wire guide through a metal cannula/needle, use caution as damage may occur to outer coating, such as the coating shearing off the wire guide.

POTENTIAL ADVERSE EVENTS

- Acute Bleeding / hematoma
- Acute pyelonephritis
- Injury to adjacent organs
- Hematuria
- Infection / UTI / urosepsis
- Mucosal tear / laceration / abrasion
- Pain
- Ureteral laceration / avulsion / tear / lesion
- Ureteral perforation / false passage / urinoma
- Ureteral stenosis
- Ureteral stricture
- Urinary tract injury
- Vascular injury

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. Remove the wire guide from the holder.

3. Advance the flexible end of the wire guide into the device and position to the desired anatomical location.

NOTE: Endhole size and length of the device must be taken into consideration to ensure proper fit between wire guide and device.

4. After wire guide use is complete, remove from patient and discard.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2024
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd