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ZH
80-83

T - P V A - R E V . 3

Polyvinyl Alcohol Foam Embolization Particles

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Polyvinyl Alcohol Foam Embolization particles are packaged in individual bottles containing 1 cc of particles.

Order Number	Color Code	Particle Size
PVA-100	black	90-180 microns
PVA-200	green	180-300 microns
PVA-300	purple	300-500 microns
PVA-500	red	500-710 microns
PVA-700	blue	710-1000 microns
PVA-1000	orange	1000-1400 microns
PVA-1500	yellow	1400-2000 microns
PVA-2000	gray	2000-2800 microns

INTENDED USE

Polyvinyl Alcohol Foam Embolization Particles are intended for embolization of the blood supply to hypervascular tumors and arteriovenous malformations, including use in intracranial embolization. The product is intended for use by physicians trained and experienced in embolization procedures in the targeted area. Standard techniques for embolization procedures should be employed.

CONTRAINDICATIONS

- Presence or suspected presence of severe atheromatous disease
- Presence or suspected presence of unfavorable patient anatomy, such as vascular configurations that do not allow superselective catheter placement
- Inadequate vessel diameter to accept emboli
- Vascular resistance peripheral to the feeding vessels that will not allow emboli to be carried into the lesion
- Inappropriate vascular anatomy such as feeding vessels smaller than the distal branches from which they emerge
- Inappropriate vascular anatomy such as extra-to-intracranial anastomoses or shunts
- Inappropriate vascular anatomy such as the presence of collateral vessel pathways that, if embolized, could endanger normal tissues
- Presence or likely onset of hemorrhage
- Presence or likely onset of vasospasm
- Patient intolerance to temporary occlusion of targeted vessels
- Presence of target vessels leading directly to cranial nerves

WARNINGS

- Neurological deficit, ischemic stroke or ischemic infarct can occur from occlusion of normal vessels by emboli.
- Artificial embolization may not occlude all arteries feeding a large arteriovenous malformation. If treatment is incomplete, the possibility of subsequent hemorrhage and/or development of alternative feeding routes may persist.
- The packaging of this product contains natural rubber latex, which may cause allergic reactions.

PRECAUTIONS

- Small contaminating particles found on the angiography table may cause foreign body reactions or actual infection. The physician should use the utmost caution to avoid contaminants during preparation of the device for use.
- Use of angiography for preoperative evaluation, operative control and postoperative follow-up is recommended.
- Use of artificial embolization devices requires careful evaluation of the vascular network associated with the lesion.
- Smaller particles are reported to be more likely to cause cranial nerve palsies and ischemic infarction because of their ability to block vessels at the precapillary level. The physician's experience must be the final judge as to the amount of the device to use, the size of the particles to use, and even whether a treatment should be undertaken.
- The highest degree of caution is recommended in the presence of visible extra-to-intracranial shunts, and in the area of the cranial nerves.
- Performing therapeutic embolizations deliberately to occlude blood vessels is a high-risk procedure. Appropriate facilities should be available for coping with the potential complications of the procedure.
- As with any surgical procedure, strict attention to sterile technique is required.
- Emboli of appropriate size must be selected by the physician, based upon the lesion to be treated and the measurements from angiography.
- As treatment progresses, the vessel will typically accept less emboli. Slowing or termination of emboli acceptance may occur when the vessel or lesion is occluded by prior emboli, or in the presence of vasospasm or severe atheromatous disease. Care must be taken to prevent continued infusion, as continued infusion may result in emboli reflux into normal vasculature, creating the potential of normal tissue ischemic infarction.
- Termination of emboli infusion is recommended prior to complete vessel occlusion.
- The use of tapered tip delivery catheters is not recommended, as emboli may obstruct the catheter tip.

ENGLISH

- Catheter obstruction by emboli may result if particles are inadequately suspended in non-ionic contrast. Care should be taken to ensure that emboli are free-floating in a sufficient volume of contrast.
- In the event of catheter obstruction, remove the catheter from the patient. Do not attempt to clear the catheter by forceful injection, passage of wire guides, or passage of other instruments. Discard the obstructed catheter, as damage to the device may have occurred.
- Only non-ionic contrast agents indicated for use in the anatomy to be embolized are recommended.

POTENTIAL ADVERSE EVENTS

- Catheter tip thrombosis and subsequent dislodgment
- Spasm of the artery adjacent to the catheter tip
- Rupture of a nearby saccular aneurysm
- Passage of emboli into normal vessels adjacent to the treatment/delivery site
- Passage of emboli through the lesion and into normal vasculature, resulting in normal tissue damage
- Reflux of emboli into normal vasculature, including arterial beds such as the internal carotid artery, pulmonary or coronary circulations, resulting in normal tissue damage
- Subarachnoid hemorrhage from recurrent bleeding of a vascular malformation or from rupture of an associated aneurysm. The relationship between embolization and the appearance of aneurysms or other vascular lesions is unclear.
- Foreign body reactions necessitating medical intervention
- Infection necessitating medical intervention
- Capillary bed saturation and tissue damage
- Ischemic stroke or ischemic infarction
- Vessel or lesion rupture resulting in hemorrhage
- Recurrent hemorrhage or vasospasm
- Vessel recanalization requiring re-treatment
- As with any embolization device, patient injury including permanent disability or death may occur as a result of its use.

PRODUCT RECOMMENDATIONS

Suggested Minimum Catheter Internal Diameter Requirements

Order Number	Particle Size	Minimum Catheter I.D.
PVA-100	90-180 microns	0.018 inches (0.46 mm)
PVA-200	180-300 microns	0.018 inches (0.46 mm)
PVA-300	300-500 microns	0.018 inches (0.46 mm)
PVA-500	500-710 microns	0.018 inches (0.46 mm)
PVA-700	710-1000 microns	0.044 inches (1.22 mm)
PVA-1000	1000-1400 microns	0.044 inches (1.22 mm)
PVA-1500	1400-2000 microns	0.044 inches (1.22 mm)
PVA-2000	2000-2800 microns	0.044 inches (1.22 mm)

Sizing parameters listed above may vary, as hydrated PVA is deformable and may pass through catheter lumen of smaller size. Use of PVA outside of the parameters above requires pre-testing and careful monitoring to prevent catheter obstruction. Testing was done based on non-hydrated particle sizes.

INSTRUCTIONS FOR USE

1. Prior to use, quality baseline angiography is essential to determine the vascular supply to the lesion in question. The angiogram should demonstrate the route of catheter entry and should be extensive enough to allow the physician to plan the entire interventional procedure. Extreme care should be taken to search for potentially dangerous collateral pathways.
2. Following percutaneous entry of the vascular system (generally by Seldinger technique), the carrier catheter is introduced as close as possible to the hypervascular lesion. Good runoff can be assured with real-time fluoroscopy.
3. After properly positioning the catheter, the physician should inspect the embolization particles. The appropriate size particles must be chosen based upon the lesion to be treated and the measurements taken from the baseline angiogram. There should not be a break in the package integrity; otherwise, sterility may be compromised. The particles should be visually inspected for uniformity of size.

4. Approximately 10 cc of non-ionic contrast agent should be introduced into the bottle containing the embolization particles. Approximately 10 minutes should be allowed, with intermittent agitation, for adequate mixing of contrast and embolization particles.
5. The physician should then inspect the particles to ensure that no part of the rubber stopper has been inadvertently introduced during the needle puncture of the vial.
6. The physician should also inspect the catheter through which the embolization particles are to be delivered. There should be no tip taper on the delivery catheter.
7. When the physician is ready to introduce the particles, a small amount of contrast agent should be first injected through the delivery catheter to ensure that the tip has not moved and that adequate runoff exists.
8. The particles are then withdrawn from the bottle into a small syringe (3 cc or larger to prevent excess pressurization of the delivery catheter). Intermittent agitation of the contrast/emboli mixture is recommended to ensure adequate suspension. **CAUTION: Catheter obstruction by emboli may result if particles are inadequately suspended in non-ionic contrast. Care should be taken to ensure that emboli are free-floating in a sufficient volume of contrast.**
9. The physician may then infuse the particles suspended in the contrast agent while directly watching the artery's acceptance of the contrast agent. **CAUTION: As the treatment progresses, the artery will typically accept less and less of the contrast agent/embolization particle mixture. The physician must watch for this change and decrease the infusion accordingly. If this change is not appreciated, inadvertent reflux into critical arteries may occur, creating the potential for ischemic infarction.**
10. The physician must use his or her judgement in determining the endpoint of the infusion. Most physicians stop infusing particles before complete occlusion of the vessel, when the acceptance rate or forward flow is approximately 1 cm/sec.
11. Upon completion of the treatment, particles will likely remain present within the catheter. The catheter should be removed from the patient with a gentle suction upon it so that no particles are inadvertently lost during the withdrawal of the catheter.

12. Should the catheter become plugged during the procedure, no attempt should be made to clean the particles from it; rather, the catheter should be removed from the patient and assessed for damage.
13. Patient follow-up after the conclusion of the procedure is recommended to assess the continued level of vascular occlusion. Angiography may be indicated.

HOW SUPPLIED

Supplied sterilized by gamma irradiation in peel-open pouches. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. Discard any unused materials after use. Do not resterilize.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47402-0489 U.S.A.