EN 18	<b>Zenith Renu® AAA Ancillary Graft with the Z-Trak™ Introduction System</b> Instructions for Use
CS 31	<b>Pomocný graft Zenith Renu® AAA se zaváděcím systémem Z-Trak™</b> Návod k použití
DA 44	<b>Zenith Renu® AAA hjælpeprotese med Z-Trak™ indføringssystem</b> Brugsanvisning
DE 57	<b>Zenith Renu® AAA-Hilfsprothese mit Z-Trak™ Einführsystem</b> Gebrauchsanweisung
EL 71	<b>Βοηθητικό μόσχευμα AAA Zenith Renu® με το σύστημα εισαγωγής Z-Trak™</b> Οδηγίες χρήσης
ES 85	<b>Endoprótesis auxiliar para AAA Zenith Renu® con el sistema de introducción Z-Trak™</b> Instrucciones de uso
FR 99	<b>Endoprothèse auxiliaire Zenith Renu® AAA avec système d'introduction Z-Trak™</b> Mode d'emploi
HU 113	<b>Z-Trak™ felvezetőrendszerrel ellátott Zenith Renu® AAA kiegészítő graft</b> Használati utasítás
IT 127	<b>Endoprotesi addominale ausiliaria Zenith Renu® con sistema di introduzione Z-Trak™</b> Istruzioni per l'uso
NL 141	<b>Zenith Renu® AAA hulprothese met het Z-Trak™ introductiesysteem</b> Gebruiksaanwijzing
NO 155	<b>Zenith Renu® AAA hjelpeimplantat med Z-Trak™ innføringssystem</b> Bruksanvisning
PL 168	<b>Stent-graft pomocniczy Zenith Renu® AAA z systemem wprowadzania Z-Trak™</b> Instrukcja użycia
PT 182	<b>Prótese auxiliar AAA Zenith Renu® com o sistema de introdução Z-Trak™</b> Instruções de utilização
SV 196	<b>Zenith Renu® AAA kompletterande graft med Z-Trak™ införingssystem</b> Bruksanvisning

## TABLE OF CONTENTS

<b>1 DEVICE DESCRIPTION</b> .....	<b>18</b>
1.1 Zenith Renu® AAA Ancillary Graft Components.....	18
1.2 Zenith Renu AAA Ancillary Graft Introduction System.....	18
1.3 Zenith AAA Endovascular Graft Ancillary Components.....	18
<b>2 INDICATIONS FOR USE</b> .....	<b>18</b>
<b>3 CONTRAINDICATIONS</b> .....	<b>18</b>
<b>4 WARNINGS AND PRECAUTIONS</b> .....	<b>18</b>
4.1 General.....	18
4.2 Patient Selection, Treatment and Follow-up.....	18
4.3 Pre-Procedure Measurement Techniques and Imaging.....	19
4.4 Device Selection.....	19
4.5 Implant Procedure.....	19
4.6 Molding Balloon Use.....	19
4.7 MRI Information.....	19
<b>5 POTENTIAL ADVERSE EVENTS</b> .....	<b>20</b>
Device-Related Adverse Event Reporting.....	20
<b>6 SUMMARY OF CLINICAL STUDIES</b> .....	<b>20</b>
<b>7 PATIENT SELECTION AND TREATMENT</b> .....	<b>20</b>
7.1 Individualization of Treatment.....	20
<b>8 PATIENT COUNSELING INFORMATION</b> .....	<b>20</b>
<b>9 HOW SUPPLIED</b> .....	<b>20</b>
<b>10 CLINICAL USE INFORMATION</b> .....	<b>21</b>
10.1 Physician Training.....	21
10.2 Inspection Prior to Use.....	21
10.3 Materials Required.....	21
10.4 Materials Recommended.....	21
10.5 Device Sizing Guidelines.....	21
Table 10.5.1 Zenith Renu AAA Main Body Extension Sizing Guide.....	21
Table 10.5.2 Zenith Renu AAA Converter Sizing Guide.....	22
<b>11 DIRECTIONS FOR USE</b> .....	<b>22</b>
General Use Information.....	22
Pre-Implant Determinants.....	22
Patient Preparation.....	22
11.1 Zenith Renu AAA Main Body Extension.....	22
11.1.1 Zenith Renu AAA Main Body Extension Preparation/Flush.....	22
11.1.2 Vascular Access and Angiography.....	22
11.1.3 Zenith Renu AAA Main Body Extension Placement.....	22
11.1.4 Zenith Renu AAA Main Body Extension Proximal (Top Stent) Deployment.....	22
11.1.5 Zenith Renu AAA Main Body Extension Distal (Bottom) Deployment.....	23
11.1.6 Docking of Top Cap.....	23
11.1.7 Molding Balloon Insertion.....	23
11.1.8 Final Angiogram.....	23
11.2 Zenith Renu AAA Converter.....	23
11.2.1 Zenith Renu AAA Converter Preparation/Flush.....	23
11.2.2 Vascular Access and Angiography.....	23
11.2.3 Zenith Renu AAA Converter Placement.....	23
11.2.4 Zenith Renu AAA Converter Proximal (Top Stent) Deployment.....	23
11.2.5 Zenith Renu AAA Converter Distal (Bottom) Deployment.....	24
11.2.6 Docking of Top Cap.....	24
11.2.7 Iliac Leg Placement and Deployment (optional).....	24
11.2.7.1 Iliac Leg Preparation/Flush.....	24
11.2.7.2 Iliac Leg Placement and Deployment.....	24
11.2.8 Contralateral Iliac Occlusion.....	24
11.2.9 Molding Balloon Insertion.....	24
11.2.10 Final Angiogram.....	24
11.2.11 Femoral-to-Femoral Crossover.....	24
<b>12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP</b> .....	<b>24</b>
12.1 General.....	24
Table 12.1 Recommended Imaging Schedule for Patients with a Zenith Renu AAA Ancillary Graft.....	25
12.2 Contrast and Non-Contrast CT Recommendations.....	25
Table 12.2 Acceptable Imaging Protocols.....	25
12.3 Abdominal Radiographs.....	25
12.4 Ultrasound.....	25
12.5 MRI Information.....	25
12.6 Additional Surveillance and Treatment.....	26
<b>13 POST-APPROVAL REGISTRY</b> .....	<b>26</b>
Summary of Post-Approval Study Methods.....	26
Table 13.1 Pre-existing Grafts Requiring Treatment with the Zenith Renu AAA Ancillary Graft.....	26
Table 13.2 Failure Modes of Pre-existing Grafts Requiring Treatment with the Zenith Renu AAA Ancillary Graft.....	27
Table 13.3 Clinical and imaging follow-up data as of February 8, 2011.....	27
Summary of Post-Approval Study Results.....	27
Table 13.4 Incidence of Adverse Events within 30 days, by Category.....	27
Table 13.5 Secondary Interventions.....	28
Table 13.6 Conversions.....	28
Table 13.7 Deaths.....	29
<b>14 PATIENT TRACKING INFORMATION</b> .....	<b>30</b>

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
Fig. 31


Fig. 32

## ZENITH RENU® AAA ANCILLARY GRAFT WITH THE Z-TRAK™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**CAUTION: All contents of the outer pouch (including the introduction system and the endovascular grafts) are supplied sterile, for single use only.**

For the Zenith product line there are four applicable Suggested Instructions For Use (IFU). This IFU describes the suggested Instructions for Use for the Zenith Renu AAA Ancillary Graft (main body extension and converter configurations). For information on other Zenith components, please refer to the following Suggested Instructions for Use:

- Zenith AAA Endovascular Graft (Zenith AAA Endovascular Graft main body, iliac legs, and occluder);
- Zenith Flex® AAA Endovascular Graft (Zenith Flex AAA Endovascular Graft main body and iliac legs);
- Zenith AAA Endovascular Graft Ancillary Components (main body extension, iliac leg extension, converter, and iliac plug); and
- Coda® Balloon Catheter.

### 1 DEVICE DESCRIPTION

#### 1.1 Zenith Renu AAA Ancillary Graft Components

The Zenith Renu AAA Ancillary Graft is available as either a straight component (Zenith Renu AAA main body extension) system or a longer, tapered component (Zenith Renu AAA converter) system for secondary endovascular intervention in patients having received prior endovascular repair of infrarenal abdominal aortic or aortic iliac aneurysms in which there is inadequate proximal fixation or seal. (Figs. 1 and 2) The Zenith Renu AAA converter module may be used independently or in combination with a Zenith AAA Endovascular Graft iliac leg.

These graft modules are an extension of the Zenith product line. They are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z™ stents with braided polyester and monofilament polypropylene suture. The modules are fully stented to provide stability and the expansive force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the necessary attachment and seal of the graft to the landing zones.

The bare suprarenal stent at the proximal end of the graft contains bars that are placed at 3 mm increments for additional fixation of the device. To facilitate fluoroscopic visualization of the stent graft, four gold radiopaque markers are positioned in a circumferential orientation within 2 mm of the most proximal aspect of the graft material.

#### 1.2 Zenith Renu AAA Ancillary Graft Introduction System

Both the Zenith Renu AAA main body extension and the Zenith Renu AAA converter are shipped preloaded onto the Z-Trak Introduction System. (Figs. 3 and 4) The system has a sequential deployment method with built-in features to provide continuous control of the endovascular graft throughout the deployment procedure. The Z-Trak Introduction System is designed for precise positioning and allows readjustment of the final graft position before deployment of the barbed suprarenal stent. Delivery systems are available in 18, 20, and 22 French sizes, depending upon component diameter. Dual trigger-wire release mechanisms lock the endovascular graft onto the introduction system until released by the physician. All systems are compatible with a .035 inch wire guide. For added hemostasis, the Capto® Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The introduction system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

#### 1.3 Zenith AAA Endovascular Graft Ancillary Components

Additional ancillary endovascular components (main body extensions, iliac leg extensions and iliac plugs) are available. Refer to the Zenith AAA Endovascular Graft Ancillary Components Instructions for Use for more information.

### 2 INDICATIONS FOR USE

The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is indicated for secondary endovascular intervention in patients having received prior endovascular repair of infrarenal abdominal aortic or aortic iliac aneurysms in which there is inadequate proximal fixation or seal with:

- Adequate iliac/femoral access compatible with the required introduction systems,
- Adequate proximal fixation site:
  - with a length from the lowest renal artery to the bifurcation of the previously placed endovascular graft of >43 mm for the main body extension and >37 mm for the converter,
  - with a diameter measured outer wall to outer wall of ≥18 mm and ≤32 mm,
  - with an angle <60 degrees relative to the long axis of the aneurysm, and
  - with an angle <45 degrees relative to the axis of the suprarenal aorta.
- Adequate distal fixation site:
  - for the main body extension, distal fixation site within a graft segment of ≤34 mm in diameter and one Cook-Z stent in length, with more overlap length being preferred,
  - for the converter used without an iliac leg, distal fixation site within a graft segment of ≤12 mm in diameter and ≥17 mm (one Cook-Z stent) in length, with more overlap length being preferred,
  - for the converter used in combination with the iliac leg, distal fixation site 7.5 to 20 mm in diameter (measured outer wall to outer wall) and >10 mm in length with 20 to 30 mm being preferred.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Standard techniques for placement of arterial access sheaths, angiographic catheters and wire guides should be employed.

### 3 CONTRAINDICATIONS

The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), polypropylene, or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

### 4 WARNINGS AND PRECAUTIONS

#### 4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

- The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 10.1, Physician Training.**
- The Zenith Renu AAA Ancillary Graft is not intended for primary endovascular treatment of patients with abdominal aortic or aortic iliac aneurysms. It is intended for use in patients in whom an endovascular graft has already been placed.
- The Zenith Renu AAA Ancillary Graft is designed to provide positive proximal fixation but may not address deficiencies in previously implanted endovascular grafts or correct the clinical problem caused by the pre-existing graft.
- The effect of the Zenith Renu AAA Ancillary Graft on the pre-existing graft or the effect of the pre-existing graft on the Zenith Renu AAA Ancillary Graft is not known. Mechanical interactions between the Zenith Renu AAA Ancillary Graft and the pre-existing graft have not been fully evaluated or characterized.
- Monitoring of the Zenith Renu AAA Ancillary Graft and pre-existing graft combination through clinical and imaging follow-up is required after implantation of the Zenith Renu AAA Ancillary Graft for the following reasons:
  - all endovascular grafts require follow-up to ensure continued success of the treatment over time,
  - the Zenith Renu AAA Ancillary Graft is intended to be used with a pre-existing endovascular graft having inadequate proximal fixation or seal. It may not prevent migration of the pre-existing graft,
  - the multiplicity of graft manufacturers, variety of endovascular graft failure modes, and considerable variability in anatomy and health status between patients precluded evaluation of the Zenith Renu AAA Ancillary Graft in every possible clinical situation to determine patient outcome.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at least yearly intervals thereafter. More frequent follow-up may be appropriate for all patients receiving the Zenith Renu AAA Ancillary Graft until stability of the Zenith Renu AAA Ancillary Graft, pre-existing graft, and aneurysm has been established.
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing an enlarging aneurysm, unacceptable fixation length (both the vessel and component overlaps) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft and/or leaks may be required to undergo secondary interventions or surgical procedures.

#### 4.2 Patient Selection, Treatment and Follow-Up

- Use of the Zenith Renu AAA converter (and not the Zenith Renu AAA main body extension) is recommended in the following:
  - length from the lowest renal artery to the bifurcation of the pre-existing graft is ≤43 mm (Refer to section 10.5 for specific device sizing guidelines),
  - length of pre-existing graft extending from the top of the graft to the graft bifurcation is less than one Cook-Z stent in length (Refer to section 10.5 for specific device sizing guidelines),
  - proximal pre-existing graft diameter is ≥34 mm,
  - adequate oversizing for both the pre-existing graft and the aortic neck is not possible due to significant diameter differences between the pre-existing graft and the aortic neck, or
  - pre-existing graft is unstable (e.g., proximal aspect of the pre-existing graft is situated within the aneurysm sac).
- Use of the Zenith Renu AAA converter in combination with an iliac leg graft (and not without an iliac leg graft) is recommended in the following:
  - pre-existing graft is composed of PTFE graft material, or
  - the distal fixation site within the pre-existing graft leg is >12 mm in diameter.
- The Zenith Renu AAA main body extension may be used when all of the following apply:
  - pre-existing graft is not composed of PTFE graft material,
  - length from the lowest renal artery to the bifurcation of the pre-existing graft is >43 mm (Refer to section 10.5 for specific device sizing guidelines), and
  - distal fixation site (within the pre-existing graft body) is ≤34 mm in diameter and at least one Cook-Z stent in length, with more overlap length being preferred.
- Use of the Zenith Renu AAA converter without an iliac leg graft may be considered if all of the following apply:
  - pre-existing graft is not composed of PTFE graft material, and
  - distal fixation site (within the pre-existing graft leg) is ≤12 mm in diameter and at least 17 mm (one Cook-Z stent) in length, with more overlap length preferred.
- Use of the Zenith Renu AAA converter with an iliac leg graft may be considered when the distal fixation site (arterial) is 7.5 to 20 mm in diameter (measured outer wall to outer wall) and >10 mm in length, with 20 to 30 mm being preferred.
- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees for infrarenal neck to axis of AAA or >45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic sealing length (<15 mm); an inverted funnel shape (>10% increase in diameter over 15 mm of proximal aortic sealing length); and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic sealing zone and distal iliac artery interface (if using Zenith Renu AAA converter device). In the presence of anatomical limitations, a longer neck may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the attachment and sealing at the fixation sites. Proximal aortic sealing zones exhibiting these key anatomical elements may be more conducive to graft migration or endoleak.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall), morphology (minimal tortuosity, occlusive disease and/or calcification), and pre-existing graft diameter should be compatible with vascular access techniques and delivery systems of a 16 to 22 French vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization. A vascular conduit technique may be necessary to achieve success in some patients.
- The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks.

Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.

- The safety and effectiveness of the Zenith Renu AAA Ancillary Components with the Z-Trak Introduction System have not been evaluated in the following patient populations:
  - traumatic aortic injury
  - leaking, pending rupture or ruptured aneurysms
  - mycotic aneurysms
  - pseudoaneurysms resulting from previous graft placement
  - uncorrectable coagulopathy
  - indispensable mesenteric artery
  - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
  - concomitant thoracic aortic or thoracoabdominal aneurysms
  - active systemic infections
  - pregnant or nursing females
  - morbidly obese patients
  - less than 18 years of age
  - less than 15 mm in length or greater than 60 degrees angulation of the proximal aortic neck relative to the long axis of the aneurysm
  - a neck diameter measured outer wall to outer wall greater than 32 mm and less than 18 mm
  - iliac distal fixation site less than 10 mm in length and less than 7.5 mm and greater than 20 mm in diameter (measured outer wall to outer wall)
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

#### 4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Renu AAA Ancillary Graft. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during procedural angiography such that the origins of the renal arteries, and particularly the lowest patent renal artery, are well demonstrated prior to deployment of the proximal edge of the graft material (sealing stent) of the main body. Additionally, angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg component.

#### Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

#### Lengths

- Utilizing CT, length measurements should be determined to accurately assess infrarenal proximal neck length as well as planning main body extension sizes, converter configurations and components for the Zenith Renu AAA Ancillary Graft. These reconstructions should be performed in sagittal, coronal, and 3-D.
- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
  - The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
  - After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

#### 4.4 Device Selection

Strict adherence to the Zenith Renu AAA Ancillary Graft IFU sizing guide is strongly recommended when selecting the appropriate device size (Tables 10.5.1 through 10.5.2). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

#### 4.5 Implant Procedure

(Refer to **Section 11, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Renu AAA Ancillary Graft and assure accurate apposition to the aortic wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Renu AAA Ancillary Graft.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Renu AAA Ancillary Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis. Vessel occlusion may occur.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the

covered proximal stent may result in damage to the stent graft and/or vessel injury.

- Inaccurate placement and/or incomplete sealing of the Zenith Renu AAA Ancillary Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith Renu AAA Ancillary Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of a Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g., adequate hydration).
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- The Zenith Renu AAA Ancillary Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional and angiographic devices in the region of the suprarenal stent.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization or rupture of the aneurysm.
- Avoid damaging the graft or disturbing graft positioning after placement in the event of reinstrumentation (secondary intervention) of the graft is necessary.
- Before deployment of the suprarenal stent, verify that the position of the access wire guide extends just distal to the aortic arch.
- Care should be taken not to displace the pre-existing graft during the placement and deployment of the Zenith Renu AAA Ancillary Graft.

#### 4.6 Molding Balloon Use

- **CAUTION: When using molding balloon in the Renu converter, do not inflate molding balloon in vessel outside of graft.** For the converter with a single internal sealing stent, molding balloon inflation within the proximal aspect of the graft component may result in part of the balloon protruding out into the vessel if the molding balloon is not properly selected. Molding balloon position should be carefully monitored during inflation to assure complete positioning within the graft.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

#### 4.7 MRI Information

Nonclinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

##### Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Highest spatial magnetic gradient field of 720 Gauss/cm

Nonclinical evaluation was conducted in a 3.0 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

##### MRI-Related Heating

###### 1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg, for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

###### 3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

##### Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available for seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 100,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners.

Mail: MedAlert Foundation International  
2323 Colorado Avenue  
Turlock, CA 95382

Phone: 888-633-4298 (toll free)  
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

## 5 POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Damage or dislodgement of pre-existing graft
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprostheses: improper placement; incomplete deployment; migration; component separation; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic (local or systemic) complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

## Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Renu AAA Ancillary Graft should be reported to Cook immediately. For customers in the United States, to report an incident, call the Customer Relations Department at 1-800-457-4500 (24-hour) or 1-812-339-2235. For customers outside the United States, please call your distributor.

## 6 SUMMARY OF CLINICAL STUDIES

The safety and effectiveness of the Zenith product line was evaluated in a U.S. clinical study, which did not include the Zenith Renu AAA Ancillary Graft components. For clinical study information regarding patients receiving Zenith products, please refer to the Zenith (Flex) AAA Endovascular Graft Instructions for Use. A copy is available online at [www.cookmedical.com](http://www.cookmedical.com). A post-market registry designed to collect information specific to physician experience with the Zenith Renu AAA Ancillary Graft has been completed. For more information about the Zenith Renu AAA Ancillary Graft registry, refer to **Section 13, POST-APPROVAL REGISTRY**

## 7 PATIENT SELECTION AND TREATMENT

(See **Section 4, WARNINGS AND PRECAUTIONS**)

### 7.1 Individualization of Treatment

Cook recommends that the Zenith Renu AAA Ancillary Graft component diameters be selected as described in **Tables 10.5.1** and **10.5.2**. The length of the Zenith Renu AAA Ancillary Graft should extend from the lowest renal artery and have sufficient overlap with the pre-existing endovascular graft. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits previously described in **Section 6, SUMMARY OF CLINICAL STUDIES**, should be carefully considered for each patient before use of the Zenith Renu AAA Ancillary Graft. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy.
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity).
- Patient's suitability for open surgical repair.
- Patient's anatomical suitability for endovascular repair.
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Renu AAA Ancillary Graft.
- Ability to tolerate general, regional or local anesthesia.
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 16 to 22 French vascular introducer sheath.
- Use of the Zenith Renu AAA converter (and not the Zenith Renu AAA main body extension) is recommended in the following:
  - length from the lowest renal artery to the bifurcation of the pre-existing graft is  $\leq 43$  mm (Refer to section 10.5 for specific device sizing guidelines),
  - length of pre-existing graft extending from the top of the graft to the graft bifurcation is less than one Cook-Z stent (Refer to section 10.5 for specific device sizing guidelines),
  - proximal pre-existing graft diameter is  $> 34$  mm,

- adequate oversizing for both the pre-existing graft and the aortic neck is not possible due to significant diameter differences between the pre-existing graft and the aortic neck, or
- pre-existing graft is unstable (e.g., proximal aspect of the pre-existing graft is situated within the aneurysm sac).
- Use of the Zenith Renu AAA converter in combination with an iliac leg graft (and not without an iliac leg graft) is recommended in the following:
  - pre-existing graft is composed of PTFE graft material, or
  - the distal fixation site within the pre-existing graft leg is  $> 12$  mm in diameter.
- The Zenith Renu AAA main body extension may be used when all of the following apply:
  - pre-existing graft is not composed of PTFE graft material,
  - length from the lowest renal artery to the bifurcation of the pre-existing graft is  $> 43$  mm (Refer to section 10.5 for specific device sizing guidelines), and
  - distal fixation site (within the pre-existing graft body) is  $\leq 34$  mm in diameter and at least one Cook-Z stent in length, with more overlap length preferred.
- Use of the Zenith Renu AAA converter without an iliac leg graft may be considered if all of the following apply:
  - pre-existing graft is not composed of PTFE graft material, and
  - distal fixation site (within the pre-existing graft leg) is  $\leq 12$  mm in diameter and at least 17 mm (one Cook-Z stent) in length, with more overlap length preferred.
- Use of the Zenith Renu AAA converter with an iliac leg graft may be considered when the distal fixation site (arterial) is 7.5 to 20 mm in diameter (measured outer wall to outer wall) and  $> 10$  mm in length, with 20 to 30 mm being preferred.
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.
- Use of the Zenith Renu AAA converter typically requires occlusion of the contralateral iliac artery in conjunction with a standard femoral-to-femoral bypass procedure. For instructions for preparation and deployment of the Zenith AAA iliac plug, refer to Zenith AAA Endovascular Graft Ancillary Components Instructions for Use.

The final treatment decision is at the discretion of the physician and patient.

## 8 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair.
- Potential advantages of open surgical repair.
- Potential advantages of endovascular repair.
- The possibility that subsequent endovascular or open surgical repair of the aneurysm may be required.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of endovascular grafts with secondary endovascular intervention using additional components has not yet been established.** All patients receiving the Zenith Renu AAA Ancillary Graft should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- **All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft component) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise each patient that it is important to seek prompt medical attention if he/she experience signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.
- Procedure-related risks include cardiac, pulmonary, neurologic, bowel and bleeding complications. Device-related risks include occlusion, endoleak, migration, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5, POTENTIAL ADVERSE EVENTS**). The physician should complete the Patient I.D. Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

## 9 HOW SUPPLIED

- The Zenith Renu AAA Ancillary Graft is supplied sterile and pre-loaded in peel-open packages.
- The device is intended for single use only. Do not re-sterilize the device.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook.
- Prior to use, verify correct device(s) (quantity and size) have been supplied for the patient by matching the device(s) to the order prescribed by the physician for that particular patient.
- The Zenith Renu AAA Ancillary Graft is loaded into a 18, 20 or 22 French Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

## 10 CLINICAL USE INFORMATION

### 10.1 Physician Training

**CAUTION: Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.**

**CAUTION: The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques. Additionally, the Zenith Renu AAA Ancillary Graft should only be used by those who have received both general Zenith training and training in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System are outlined below:**

#### Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair
- Knowledge of radiographic image interpretation, device selection, planning and sizing

#### A multidisciplinary team that has combined procedural experience with:

- Femoral cutdown, femoral-to-femoral crossover, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities
- Reinstrumentation of pre-existing endovascular grafts.

#### 10.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook. Prior to use, verify correct device(s) (quantity and size) have been supplied for the patient by matching the device(s) to the order prescribed by the physician for that particular patient.

### 10.3 Materials Required

- A selection of ancillary components including main body extensions, iliac extensions and iliac plugs.
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

### 10.4 Materials Recommended

The following products are recommended for implantation of any component in the Zenith product line. For information on the use of these products, refer to the individual product's suggested instructions for use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
  - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
  - Cook .035 inch wire guides
  - Cook Nimble™ Wire Guides
- Molding Balloons; for example:
  - Cook Coda Balloon Catheter
- Introducer sets; for example:
  - Cook Check-Flo® Introducer Sets
  - Cook Extra Large Check-Flo® Introducer Sets
  - Cook Flexor Balkin Up & Over® Contralateral Introducers
- Sizing catheter; for example:
  - Cook Auros® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
  - Cook Beacon® Tip Angiographic Catheters
  - Cook Beacon Tip Royal Flush Catheters
- Entry needles; for example:
  - Cook single wall entry needles
- Endovascular Dilators; for example:
  - Cook Endovascular Dilator Sets

### 10.5 Device Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. The choice of length should be determined using the length from the lowest renal artery to the bifurcation of the pre-existing graft. Undersizing or oversizing may result in incomplete sealing or compromised flow.

**Table 10.5.1 Zenith Renu AAA Main Body Extension Sizing Guide\***

Intended Aortic Vessel Diameter <sup>1,2,3</sup> (mm)	Graft Diameter (mm)	Intended Proximal Working Length <sup>4</sup> (mm)	Overall Graft Length (mm)	Sealing Stent Length <sup>5</sup> (mm)	Introducer Sheath (Fr)
18-19	22	43-61	43	17	18
		62-80	62	17	18
		≥81	81	17	18
20-21	24	43-61	43	17	18
		62-80	62	17	18
		≥81	81	17	18
22	26	43-61	43	17	18
		62-80	62	17	18
		≥81	81	17	18
23-24	28	43-61	43	17	20
		62-80	62	17	20
		≥81	81	17	20
25-26	30	43-61	43	17	20
		62-80	62	17	20
		≥81	81	17	20
27-28	32	43-61	43	17	20
		62-80	62	17	20
		≥81	81	17	20
29-32	36	54-76	54	22	22
		77-99	77	22	22
		≥100	100	22	22

<sup>1</sup>Maximum diameter along the proximal fixation site.

<sup>2</sup>Round measured aortic diameter to nearest mm.

<sup>3</sup>Additional considerations may affect choice of diameter.

<sup>4</sup>Intended Proximal Working Length is the length from the lowest renal artery to the bifurcation of the pre-existing graft.

<sup>5</sup>Sealing Stent Length is the length of the most distal Cook-Z stent, which is the minimum length of the distal fixation site.

\*All dimensions are nominal.

**NOTE:** There should be a minimum of 2 mm oversizing with respect to the diameter of the pre-existing graft.

**Table 10.5.2 Zenith Renu AAA Converter Sizing Guide\***

Intended Aortic Vessel Diameter <sup>1,2,3</sup> (mm)	Graft Diameter (mm)	Intended Proximal Working Length <sup>4</sup> (mm)	Proximal Graft Length (mm)	Overall Graft Length (mm)	Sealing Stent Length <sup>5</sup> (mm)	Introducer Sheath (Fr)
18-19	22	37-58	37	91, 108, 125	17	18
		≥59	59	113, 130, 147	17	18
20-21	24	37-58	37	91, 108, 125	17	18
		≥59	59	113, 130, 147	17	18
22	26	37-58	37	91, 108, 125	17	18
		≥59	59	113, 130, 147	17	18
23-24	28	37-58	37	91, 108, 125	17	20
		≥59	59	113, 130, 147	17	20
25-26	30	37-58	37	91, 108, 125	17	20
		≥59	59	113, 130, 147	17	20
27-28	32	37-58	37	91, 108, 125	17	20
		≥59	59	113, 130, 147	17	20
29-32	36	52-79	52	99, 116, 133	22	22
		≥80	80	127, 144, 161	22	22

<sup>1</sup>Maximum diameter along the proximal fixation site.  
<sup>2</sup>Round measured aortic diameter to nearest mm.  
<sup>3</sup>Additional considerations may affect choice of diameter.  
<sup>4</sup>Intended Proximal Working Length is the length from the lowest renal artery to the bifurcation of the pre-existing graft.  
<sup>5</sup>Sealing Stent Length is the length of the most distal Cook-Z stent, which is the minimum length of the distal fixation site.  
<sup>6</sup>All dimensions are nominal.

**NOTE:** If the converter is used with an iliac leg, the distal fixation site (arterial) should be 7.5 to 20 mm in diameter (measured outer wall to outer wall) and >10 mm in length, with 20 to 30 mm being preferred. If the converter is used without an iliac leg, the distal fixation site (within the pre-existing graft leg) should be ≤12 mm in diameter and at least 17 mm (one Cook-Z stent) in length, with more overlap length being preferred.

**11 DIRECTIONS FOR USE**

Prior to use of a Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System, review this *Suggested Instructions for Use* booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

**General Use Information**

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of a Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System. The Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is compatible with .035 inch diameter wire guides.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

**Pre-Implant Determinants**

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the delivery system.
2. Angulation of aortic neck, aneurysm and iliac arteries.
3. Quality of the aortic sealing zone.
4. Diameters of infrarenal aortic neck and distal iliac arteries.
5. Distance from renal arteries to the bifurcation of the pre-existing graft.
6. Length from the bifurcation of the pre-existing graft to the internal iliac arteries/attachment site(s), when using a converter with an iliac leg.
7. Consider the degree of vascular calcification.

**Patient Preparation**

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Expose selected common femoral artery using standard surgical technique. When a femoral-to-femoral crossover will be required as part of the procedure, both common femoral arteries should be exposed.
4. Establish adequate proximal and distal vascular control of selected femoral vessel.

**NOTE:** The Zenith Renu AAA Ancillary Graft components are designed to be introduced through an exposed common femoral artery on the selected side of introduction. Angiography at the implant site may be achieved by the use of a straight angiographic catheter inserted from the contralateral side, either by surgical exposure or percutaneous approach.

**NOTE: Section 11.1** describes Zenith Renu AAA main body extension deployment and **Section 11.2** describes Zenith Renu AAA converter deployment.

**11.1 Zenith Renu AAA Main Body Extension (Fig. 3)**

**NOTE:** The deployment sequence for the Zenith Renu AAA main body extension is comparable to the deployment sequence for the Zenith Flex AAA main body component.

**11.1.1 Zenith Renu AAA Main Body Extension Preparation/Flush**

1. Remove gray-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away<sup>®</sup> sheath from the back of the Captor Hemostatic Valve. (Fig. 5) Elevate distal tip of system and flush through the stopcock on the Captor Hemostatic Valve until fluid emerges from the sideport near the tip of the introducer sheath. (Fig. 6) Continue to inject a full 20 ml of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

**NOTE:** Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

**NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pad in saline solution and use to wipe Flexor introducer sheath to activate hydrophilic coating. Hydrate both sheath and dilator liberally.

**11.1.2 Vascular Access and Angiography**

**NOTE:** The Zenith Renu AAA main body extension is designed to be introduced through an exposed common femoral artery on the selected side of introduction. Angiography at the implant site may be achieved by the use of a straight angiographic catheter inserted from the contralateral side, either by surgical exposure, or percutaneous approach.

1. Puncture the selected common femoral artery(ies) using standard technique with an 18 or 19 UT gauge arterial needle. Upon vessel entry, insert:
  - Wire guides – standard .035 inch diameter, 145 cm long, J tip or Bentson Wire Guide
  - Appropriate size sheaths (e.g., 6 or 8 French)
  - Flush catheter (often radiopaque sizing catheters – e.g., centimeter sizing catheter or straight flush catheter)
2. Perform angiography to identify level(s) of renals, bifurcation of pre-existing graft and iliac bifurcations.

**NOTE:** If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

**11.1.3 Zenith Renu AAA Main Body Extension Placement**

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

**NOTE:** Monitor the patient’s coagulation status throughout the procedure.

3. Replace J wire with stiff wire guide (LES) .035 inch, 260 cm long and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
4. Insert the Zenith Renu AAA main body extension delivery system over the wire, into the femoral artery.

**CAUTION: Maintain wire guide position during delivery system insertion.**

**CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).**

5. Advance delivery system until the four gold radiopaque markers (Fig. 8) are just inferior to the most inferior renal orifice.

**NOTE:** To ensure patency of renal arteries, recognize that the proximal graft markers are 2 mm below the proximal edge of the graft material.

6. Repeat the angiogram to verify the four gold radiopaque markers are just inferior to the most inferior renal orifice. Confirm that the position of distal end of Zenith Renu AAA main body extension will provide at least one Cook-Z stent (17 or 22 mm) overlap in the aortic portion of the pre-existing stent-graft.

**NOTE:** To prevent occlusion of one or both iliac arteries, ensure the device is deployed so that the distal portion is above the bifurcation of the pre-existing graft.

7. Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Fig. 9)
8. Use the gripper to stabilize the gray positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first covered stent by withdrawing the sheath while monitoring device location.
9. Without moving the table, decrease magnification to check position of the distal end of the device and location of renal arteries. Proceed with deployment until the distal segment of the Zenith Renu AAA main body extension is fully deployed. (Fig. 10) Stop withdrawing sheath.

**11.1.4 Zenith Renu AAA Main Body Extension Proximal (Top Stent) Deployment**

**NOTE:** The proximal stent is secured within the top cap by a trigger-wire.

1. Perform angiography through an angiographic catheter to verify position of the endovascular graft does not restrict flow to the renal arteries. If necessary, carefully reposition the covered portion of the endovascular graft with respect to the renal arteries. (Repositioning can only take place over a small range of distance at this stage.)

**NOTE:** To ensure patency of renal arteries, recognize that the proximal graft markers are 2 mm below the proximal edge of the graft material.

**CAUTION:** During proximal trigger-wire removal, top cap advancement, and subsequent suprarenal stent deployment, verify that the position of the main body extension wire guide extends just distal to the aortic arch and that support of the system is maximized.

2. Remove the safety lock from the top stent trigger-wire release mechanism. Under fluoroscopy, withdraw and remove the trigger-wire by sliding the top stent trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Fig. 11)

If resistance is felt or system bowing is noticed, the trigger-wire is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation. If unable to remove the top stent trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy:

- a. Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith main body.
  - b. Retighten the pin vise.
  - c. Remove the top stent trigger-wire release mechanism.
3. Loosen the pin vise. (Fig. 12) Control the position of the graft by stabilizing the gray positioner of the introducer.
  4. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the device until the top stent is fully deployed. (Figs. 13 and 14) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

**NOTE:** Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

**NOTE:** Once the barbed suprarenal stent has been deployed, further attempts to reposition graft are not recommended.

**WARNING:** The Zenith Renu AAA Ancillary Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

### 11.1.5 Zenith Renu AAA Main Body Extension Distal (Bottom) Deployment

**NOTE:** The distal stent is still secured by the trigger-wire.

1. Remove the safety lock from the trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Fig. 15)

### 11.1.6 Docking of Top Cap

1. Loosen the pin vise. (Fig. 16)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the gray positioner over the inner cannula until it docks with the top cap. (Figs. 17, 18 and 19)

**NOTE:** If resistance occurs, slightly rotate gray positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and gray positioner through the graft and through the sheath by pulling on the inner cannula. Leave the sheath and wire guide in place.

**NOTE:** Maintain position of sheath and wire guide.

5. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops. (Fig. 20)

### 11.1.7 Molding Balloon Insertion

For warnings, precautions and instructions, please refer to the Cook Coda Balloon Catheter Instructions for Use.

**NOTE:** Prior to molding balloon use, carefully withdraw any accessory devices (e.g., angiographic catheter) from the contralateral limb until the tip of the catheter is below the distal edge of the graft extension.

1. Prepare molding balloon as follows:
  - Flush wire lumen with heparinized saline.
  - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve of the Zenith Renu AAA main body extension introduction system to the level of the renal arteries. Maintain proper sheath position.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

**CAUTION:** Do not inflate balloon in vessel outside of the graft.

5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (Fig. 21)

**CAUTION:** Confirm complete deflation of balloon prior to repositioning.

**CAUTION:** Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

6. Withdraw the molding balloon to the distal segment of the Zenith Renu AAA main body extension and expand.

**CAUTION:** Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

7. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
8. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

### 11.1.8 Final Angiogram

1. Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks.
2. Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.

**NOTE:** If endoleaks or other problems are observed and require intervention, consider use of additional Zenith AAA Endovascular Graft Ancillary Components.

3. Repair vessels and close in standard surgical fashion.

### 11.2 Zenith Renu AAA Converter (Fig.4)

**NOTE:** The deployment sequence for the Zenith Renu AAA converter is comparable to the deployment sequence for the Zenith Flex AAA main body component.

#### 11.2.1 Zenith Renu AAA Converter Preparation/Flush

1. Remove gray-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from the back of the hemostatic valve. (Fig. 5) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introducer sheath. (Fig. 6) Continue to inject a full 20 ml of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

**NOTE:** Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

**NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pad in saline solution and use to wipe Flexor introducer sheath to activate hydrophilic coating. Hydrate both sheath and dilator liberally.

#### 11.2.2 Vascular Access and Angiography

**NOTE:** The Zenith Renu AAA converter is designed to be introduced through an exposed common femoral artery on the selected side of introduction. Angiography at the implant site may be achieved by the use of a straight angiographic catheter inserted from the contralateral side, either by surgical exposure or percutaneous approach.

1. Puncture the selected common femoral artery(ies) using standard technique with an 18 or 19 UT gage arterial needle. Upon vessel entry, insert:
  - Wire guides – standard .035 inch diameter, 145 cm long, J tip or Bentson Wire Guide
  - Appropriate size sheaths (e.g., 6 or 8 French)
  - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
2. Perform angiography to identify level(s) of renals, bifurcation of pre-existing graft and iliac bifurcations.

**NOTE:** If fluoroscope angulation is used with an angulated neck, it may be necessary to perform angiograms using various projections.

#### 11.2.3 Zenith Renu AAA Converter Placement

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

**NOTE:** Monitor the patient's coagulation status throughout the procedure.

3. Replace J wire with stiff wire guide (LES) .035 inch, 260 cm long and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
4. Insert delivery system over the wire, into the femoral artery.

**CAUTION:** Maintain wire guide position during delivery system insertion.

**CAUTION:** To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

5. Advance delivery system until the four gold radiopaque markers (Fig. 22) are just inferior to the most inferior renal orifice.
6. Repeat the angiogram to verify the four gold radiopaque markers are just inferior to the most inferior renal orifice. If the Zenith Renu AAA converter is used independently (i.e., without an iliac leg), confirm that the position of the distal end will provide at least one Cook-Z stent overlap within the iliac portion of the pre-existing stent-graft.

**NOTE:** To ensure patency of renal arteries, recognize that the proximal graft markers are 2 mm below the proximal edge of the graft material.

7. Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Fig. 9)
8. Use the gripper to stabilize the gray positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first covered stent by withdrawing the sheath while monitoring device location.
9. Without moving the table, decrease magnification to check position of the distal end of the device and location of renal arteries. Proceed with deployment until the distal segment of the Zenith Renu AAA converter is fully deployed. (Fig. 23) Stop withdrawing sheath.

#### 11.2.4 Zenith Renu AAA Converter Proximal (Top Stent) Deployment

**NOTE:** The proximal stent is secured within the top cap by a trigger-wire.

1. Perform angiography through an angiographic catheter to verify that the position of the endovascular graft does not restrict flow to the renal arteries. If necessary, carefully reposition the covered portion of the endovascular graft with respect to the renal arteries. (Repositioning can only take place over a small range of distance at this stage.)

**NOTE:** To ensure patency of renal arteries, recognize that the proximal graft markers are 2 mm below the proximal edge of the graft material.

2. Remove the safety lock from the top stent trigger-wire release mechanism. Under fluoroscopy, withdraw and remove the trigger-wire by sliding the top stent trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Fig. 11)

If resistance is felt or system bowing is noticed, the trigger-wire is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation. If unable to remove the top stent trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy:

- a. Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith Renu AAA converter.
  - b. Retighten the pin vise.
  - c. Remove the top stent trigger-wire release mechanism.
3. Loosen the pin vise. (Fig. 12) Control the position of the graft by stabilizing the gray positioner of the introducer.

**CAUTION:** Before deployment of suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch.

4. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the device until the top stent is fully deployed. (Figs. 13 and 24) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

**WARNING:** The Zenith Renu AAA Ancillary Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

### 11.2.5 Zenith Renu AAA Converter Distal (Bottom) Deployment

**NOTE:** The distal stent is still secured by the trigger-wire.

1. Remove the safety lock. Withdraw and remove the trigger-wire by sliding the trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Fig. 15)

### 11.2.6 Docking of Top Cap

1. Loosen the pin vise. (Fig. 16)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the gray positioner over the inner cannula until it docks with the top cap. (Figs. 17, 18 and 25)

**NOTE:** If resistance occurs, slightly rotate gray positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and gray positioner through the graft and through the sheath by pulling on the inner cannula. Leave the sheath and wire guide in place.

**NOTE:** Maintain position of sheath and wire guide.

5. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops. (Fig. 20)

### 11.2.7 Iliac Leg Placement and Deployment (optional)

To aid in graft sealing, an iliac leg may be used to extend the Zenith Renu AAA converter. (Fig. 28)

#### 11.2.7.1 Iliac Leg Preparation/Flush

1. Remove gray-hubbed inner stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from the back of the Captor Hemostatic Valve. (Fig. 29) Elevate distal tip of system and flush through the stopcock on the Captor Hemostatic Valve until fluid emerges from the sideport near the tip of the introducer sheath. (Fig. 30) Continue to inject a full 20 ml of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

**NOTE:** Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

**NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

#### 11.2.7.2 Iliac Leg Placement and Deployment

**NOTE:** Ensure the Captor Hemostatic Valve on the introducer sheath is turned to the open position. (Fig. 9)

1. Utilize the converter graft wire and sheath assembly to introduce the ipsilateral iliac leg graft. Advance dilator and sheath assembly into the converter sheath.

**NOTE:** In tortuous vessels, the position of the internal iliac arteries may alter significantly with the introduction of the rigid wires and sheath systems.

2. Advance slowly until the ipsilateral iliac leg graft overlaps a minimum of one full iliac leg stent (i.e., proximal stent of iliac leg graft) inside the Renu converter. (Fig. 31)
3. Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure internal iliac patency.
4. To deploy, use the gripper to stabilize the iliac leg graft with the gray positioner while withdrawing the iliac leg sheath. (Fig. 32) If necessary, withdraw the converter sheath.
5. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to gray positioner. Tighten pin vise. Maintain sheath position while withdrawing gray positioner with secured inner cannula.
6. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops.
7. Re-check the position of the wire guides. Leave sheath and wire guides in place.
8. Refer to Section 11.2.9 Molding Balloon Insertion for molding balloon use.

#### 11.2.8 Contralateral Iliac Occlusion

Use of the Zenith Renu AAA converter system typically requires occlusion of the contralateral iliac artery. (Fig. 26) For instructions for preparation and deployment of the Zenith AAA iliac plug, refer to the Zenith AAA Endovascular Graft Ancillary Components Instructions for Use.

#### 11.2.9 Molding Balloon Insertion

**NOTE:** To avoid inflation of the molding balloon outside of the proximal aspect of the converter with a single internal sealing stent, use of an appropriately inflated 40 mm Coda balloon is recommended. Distal graft segments may be molded with a 32 mm Coda balloon. Recommended balloon inflation volumes should be strictly observed.

**NOTE:** Prior to molding balloon use, carefully withdraw any accessory devices (e.g., angiographic catheter) from the contralateral limb until the tip of the catheter is below the distal edge of the graft extension.

1. Prepare molding balloon as follows:
  - Flush wire lumen with heparinized saline.
  - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve of the Zenith Renu AAA converter introduction system to the level of the renal arteries. Maintain proper sheath position.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.
5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (Fig. 27)

**CAUTION: Do not inflate molding balloon in vessel outside of graft. For the converter with a single internal sealing stent, molding balloon inflation within the proximal aspect of the graft component may result in part of the balloon protruding out into the vessel if the molding balloon is not properly selected. Molding balloon position should be carefully monitored during inflation to assure complete positioning within the graft.**

6. Withdraw the molding balloon to the distal section of the device and expand.

**CAUTION: Confirm complete deflation of balloon prior to repositioning.**

**CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.**

7. Deflate and remove molding balloon.

**CAUTION: Do not inflate balloon in vessel outside of the iliac plug.**

**CAUTION: Do not position the iliac plug during insertion and expansion of molding balloon.**

8. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
9. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

#### 11.2.10 Final Angiogram

1. Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac artery.
2. Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.

**NOTE:** If endoleaks or other problems are observed and require intervention, consider use of additional Zenith AAA Endovascular Graft Ancillary Components.

#### 11.2.11 Femoral-to-Femoral Crossover

Perform femoral-to-femoral crossover in standard surgical fashion to revascularize the contralateral limb. Repair vessels and close in standard fashion.

## 12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

### 12.1 General

- **The long-term performance of endovascular grafts with secondary endovascular intervention using additional components has not yet been established.** All patients receiving the Zenith Renu AAA Ancillary Graft should receive enhanced follow-up.
- **All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft component) should receive enhanced follow-up.
- In the Renu Registry, proximal Type I endoleaks were identified as one of the primary failure modes for pre-existing grafts, with 96 being identified preoperatively or during the procedure. Of the proximal Type I endoleaks reported, 99% (95/96) resolved without further intervention following Renu implantation; one persisted through one-month follow-up and was converted to open surgical repair. In addition, several new proximal Type I or Type III endoleaks were identified throughout the course of the study.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at least yearly intervals thereafter. More frequent follow-up may be appropriate for all patients receiving the Zenith Renu AAA Ancillary Graft until stability of the Zenith Renu AAA Ancillary Graft, pre-existing graft, and aneurysm has been established. Patients exhibiting endoleak after treatment with the Zenith Renu AAA Ancillary Graft should be monitored closely, particularly those with either Type I or Type III endoleak, which can result in aneurysm growth/rupture.
- The Zenith Renu AAA Ancillary Graft is not intended for primary endovascular treatment of patients with abdominal aortic or aortoiliac aneurysms. It is intended for use in patients in whom an endovascular graft has already been placed.
- Due to the multiplicity of graft manufacturers, variety of endovascular graft failure modes, and considerable variability in anatomy and health status between patients, evaluation of the Zenith Renu AAA Ancillary Graft in every situation was not possible and the clinical outcome is unknown. The Zenith Renu AAA Ancillary Graft is not intended to alter the course of the pre-existing graft. Furthermore, the effect of the Zenith Renu AAA Ancillary Graft on the pre-existing graft or the effect of the pre-existing graft on the Zenith Renu AAA Ancillary Graft is not known. For these reasons, continued long-term monitoring of the Zenith Renu AAA Ancillary Graft and pre-existing graft combination through clinical and imaging follow-up is required after implantation.
- Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in Table 12.1. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. After secondary endovascular intervention, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, imaging should include: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and disease progression. Only patients with renal complications or other factors precluding the use of image contrast media should be monitored using non-contrast CT and duplex ultrasound.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.

**Table 12.1** lists the minimum imaging follow-up for patients with a Zenith Renu AAA Ancillary Graft. Patients requiring enhanced follow-up should have interim evaluations.

**Table 12.1 Recommended Imaging Schedule for Patients With a Zenith Renu AAA Ancillary Graft**

	Angiogram	CT (Contrast and non-contrast)	Abdominal Radiographs
Pre-procedure	X <sup>1</sup>	X <sup>1</sup>	
Procedural	X		
Pre-discharge (within 7 days)		X <sup>2,3,4</sup>	X
1 month		X <sup>2,3,4</sup>	X
3 month		X <sup>2,4,5</sup>	
6 month		X <sup>2,4</sup>	X
12 month (annually thereafter)		X <sup>2,4</sup>	X

<sup>1</sup>Imaging should be performed within 6 months before the procedure.

<sup>2</sup>Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast-enhanced CT scan. With ultrasound, non-contrast CT is still recommended.

<sup>3</sup>Either pre-discharge or 1 month CT recommended.

<sup>4</sup>If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended, see **Section 12.6, Additional Surveillance and Treatment**.

<sup>5</sup>Recommended if endoleak reported at pre-discharge or 1 month or if the pre-existing graft is unstable.

**12.2 Contrast and Non-Contrast CT Recommendations**

- Film sets should include all sequential images at lowest possible slice thickness (≤3 mm). DO NOT perform large slice thickness (>3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch sheets if film is used.
- Both non-contrast and contrast runs are required, with matching or corresponding table positions.
- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast-enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 12.2** lists examples of acceptable imaging protocols.

**Table 12.2 Acceptable Imaging Protocols**

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high-performance MDCT capable of >40 seconds	Spiral CT or high-performance MDCT capable of >40 seconds
Injection volume	n/a	Per institutional protocol
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E., or equivalent
Coverage – start	Diaphragm	1 cm superior to celiac axis
Coverage – finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

**12.3 Abdominal Radiographs**

The following views are required:

- Four films: supine-frontal (AP), cross-table lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise. If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views.

The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

**12.4 Ultrasound**

Ultrasound imaging may be performed in place of contrast CT when patient factors preclude the use of image contrast media. Ultrasound may be paired with non-contrast CT. A complete aortic duplex is to be videotaped for maximum aneurysm diameter, endoleaks, stent patency and stenosis. Included on the videotape should be the following information as outlined below:

- Transverse and longitudinal imaging should be obtained from the level of the proximal aorta demonstrating mesenteric and renal arteries to the iliac bifurcations to determine if endoleaks are present, utilizing color flow and color power angiography (if accessible).
- Spectral analysis confirmation should be performed for any suspected endoleaks.
- Transverse and longitudinal imaging of the maximum aneurysm should be obtained.

**12.5 MRI Information**

Nonclinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

**Static Magnetic Field**

- Static magnetic field of 3.0 Tesla or less.
- Highest spatial magnetic gradient field of 720 Gauss/cm

Nonclinical evaluation was conducted in a 3.0 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

**MRI-Related Heating**

**1.5 Tesla Systems:**

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg, for

15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

**3.0 Tesla Systems:**

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

**Image Artifact**

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available for seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 100,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners.

Mail: MedicAlert Foundation International  
2323 Colorado Avenue  
Turlock, CA 95382  
Phone: 888-633-4298 (toll free)  
209-668-3333 from outside the US  
Fax: 209-669-2450  
Web: www.medicalert.org

## 12.6 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement,  $\geq 5$  mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length
- Separation from pre-existing graft
- Pre-existing graft is unstable

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

### 13 POST-APPROVAL REGISTRY

The safety and effectiveness of the Zenith product line was evaluated in a multi-center U.S. clinical study, which did not include the Zenith Renu AAA Ancillary Graft components. For clinical study information regarding patients receiving Zenith products (other than Zenith Renu AAA Ancillary Graft), please refer to the Zenith (Flex) AAA Endovascular Graft Instructions for Use. A copy is available online at [www.cookmedical.com](http://www.cookmedical.com)

A post-approval study designed to collect information specific to physician experience with the Zenith Renu AAA Ancillary Graft has been completed. A summary of the study follows.

### Summary of Post-Approval Study Methods

#### Study Objective

The purpose of this post-market surveillance registry was to accumulate cases of physician experience with the Zenith Renu AAA Ancillary Graft to confirm device performance. Specifically, the aims were to:

1. Obtain case reports of physician experience with the Zenith Renu AAA Ancillary Graft.
2. Confirm device performance as measured by the incidence of subsequent device-related problems.

#### Study Design

The Renu registry was a non-randomized, single-arm, prospective, post-market registry of cases involving the Zenith Renu AAA Ancillary Graft. There was no limit to the number of institutions allowed to participate, providing at least one physician at the institution had completed the Zenith Training Program. The registry was intended to provide information on a minimum of 100 Renu devices. Physicians were requested to provide registration information as well as pertinent clinical and imaging results obtained during the procedure, post-procedure within 30 days, at 12 months, and annually up to 5 years. The registry utilized a centralized core lab for analysis of available imaging and an independent Clinical Events Committee (CEC) for adverse event adjudication.

## Study Population

Only patients having received prior repair of infrarenal abdominal aortic or aortoiliac aneurysms in which there was inadequate proximal fixation or seal were entered in the registry. Additionally, patient anatomy had to meet the approved indications, i.e.:

- Adequate iliac/femoral access compatible with the required introduction systems,
- Adequate proximal fixation site:
  - with a length from the lowest renal artery to the bifurcation of the previously placed endovascular graft of  $>43$  mm for the main body extension and  $>37$  mm for the converter,
  - with a diameter measured outer wall to outer wall of  $\geq 18$  mm and  $\leq 28$  mm,
  - with an angle  $<60$  degrees relative to the long axis of the aneurysm, and
  - with an angle  $<45$  degrees relative to the axis of the suprarenal aorta.
- Adequate distal fixation site:
  - for the main body extension, distal fixation site within a graft segment of  $\leq 30$  mm in diameter and  $\geq 17$  mm (one Cook-Z<sup>®</sup> stent) in length, with more overlap length being preferred,
  - for the converter used without an iliac leg, distal fixation site within a graft segment of  $\leq 12$  mm in diameter and  $\geq 17$  mm (one Cook-Z<sup>®</sup> stent) in length, with more overlap length being preferred,
  - for the converter used in combination with the iliac leg, distal fixation site 7.5 to 20 mm in diameter (measured outer wall to outer wall) and  $>10$  mm in length with 20 to 30 mm being preferred.

Recipients of the Renu device were primarily male (87%, n=131) with an average age of 77 years (n=147).

### Pre-existing Device Information

**Table 13.1** presents the types of pre-existing grafts that required treatment. The types of pre-existing grafts included 126 AneuRx<sup>®</sup>, 9 Ancure<sup>®</sup>, 6 Excluder<sup>®</sup>, 3 Talent<sup>™</sup>, 2 Vanguard<sup>™</sup>, 2 hand-made grafts, 1 Fortron<sup>™</sup>, 1 Lifepath<sup>™</sup>, and 1 Zenith. The average overall implantation time of pre-existing grafts treated with the Zenith Renu AAA Ancillary Graft was 43.4 months (n=139). The implantation times of 10 AneuRx<sup>®</sup> grafts, 1 Ancure<sup>®</sup> graft, and 1 Talent<sup>™</sup> graft were not provided.

**Table 13.1 Pre-existing Grafts Requiring Treatment with the Zenith Renu AAA Ancillary Graft**

Pre-existing Graft Type	Number	(%)	Implantation Time Prior to Treatment (months) <sup>2</sup>
AneuRx <sup>®</sup>	126	83.4%	42.1 ± 17.1 (n=116)
Ancure <sup>®</sup>	9	6.0%	55.1 ± 8.6 (n=8)
Excluder <sup>®</sup>	6	4.0%	31.3 ± 26.2 (n=6)
Talent <sup>™</sup>	3	2.0%	38.5 ± 0.7 (n=2)
Vanguard <sup>™</sup>	2	1.3%	96.0 ± 0.0 (n=2)
Other <sup>1</sup>	2	1.3%	75.0 ± 14.1 (n=2)
Fortron <sup>™</sup>	1	0.7%	44.0 (n=1)
Lifepath <sup>™</sup>	1	0.7%	38.0 (n=1)
Zenith	1	0.7%	15.0 (n=1)

<sup>1</sup> Hand-made grafts (1 aortouni-iliac and 1 bifurcated).

<sup>2</sup> The implantation times of 10 AneuRx<sup>®</sup> grafts, 1 Ancure<sup>®</sup> graft, and 1 Talent<sup>™</sup> graft were not provided.

The failure modes of the pre-existing grafts as reported by each site through the online registry are provided in **Table 13.2**.

**Table 13.2 Failure Modes of Pre-existing Grafts Requiring Treatment with the Zenith Renu AAA Ancillary Graft<sup>1</sup>**

	Pre-existing Graft									
	All	AneuRx®	Ancure®	Excluder®	Fortron™	Lifepath™	Talent™	Vanguard™	Zenith	Other <sup>2</sup>
<b>Devices Treated</b>	151	126	9	6	1	1	3	2	1	2
Endoleak	108	89	6	4	1	1	2	2	1	2
Proximal Type I	86	74	4	2	1	1	1	1	1	1
Migration	136	120	6	2	1	1	3	2	0	1
Stent Fracture/ Breakage	3	3	0	0	0	0	0	0	0	0
Graft Tear	3	2	1	0	0	0	0	0	0	0
Component Separation	2	1	0	0	0	0	0	1	0	0
Occlusion	2	1	1	0	0	0	0	0	0	0
Kink	7	5	1	0	0	0	0	0	0	1
Other <sup>4</sup>	2	1	0	1	0	0	0	0	0	0

Reported Failure Modes<sup>3</sup>

<sup>1</sup> Failure modes of pre-existing grafts are based on site-reported data as previously reported to the FDA.

<sup>2</sup> Hand-made grafts (1 aortouni-iliac and 1 bifurcated).

<sup>3</sup> Ninety-nine pre-existing grafts were reported as having multiple failure modes. Failure mode(s) of one AneuRx® graft was not provided by one implanting institution. Per Cook representative present at the procedure, the pre-existing graft had both migrated and had a proximal Type I endoleak. These failure modes were confirmed by evaluation of preoperative imaging and have been included in the registry.

<sup>4</sup> One AneuRx® was noted as having a loss of graft integrity discovered during the Renu procedure. One Excluder® was noted as having aneurysm sac growth due to the material.

The most commonly reported failure modes of pre-existing grafts were proximal Type I endoleak (86 cases) and migration (136 cases). Other failure modes included additional endoleak, stent fracture, graft tear or leakage, component separation, occlusion, and kink. More than one failure mode was reported in 99 cases.

#### Data Source

The Zenith Renu AAA Ancillary Graft post-market surveillance registry represents the only data source for this report.

#### Key Variables Assessed for Safety and Effectiveness

The key variables assessed for safety and effectiveness of the Zenith Renu AAA Ancillary Graft registry were as follows:

- Safety: death, rupture, conversion, and other reported adverse events
- Effectiveness: proximal attachment status (i.e., barb separation, stent fracture/breakage, position of the Zenith Renu AAA Ancillary Graft with respect to the fixation site and integrity of the seal with the aorta by the absence of a Type I endoleak as well as the pre-existing endovascular graft by the absence of a Type III endoleak) and need for secondary intervention.

#### Total Number of Enrolled Study Sites and Subjects, Follow-up Rate

The registry included case data for 151 patients from 95 US institutions, with very few institutions performing more than one or two cases. The follow-up rate for patients treated with the Zenith Renu is presented in **Table 13.3**.

**Table 13.3 Clinical and imaging follow-up data as of February 8, 2011**

Visit	Eligible for Follow-Up <sup>1</sup>	Percent of Data Available Total <sup>2</sup>
30-day	148	100% (148/148)
12-month	132	98.5% (130/132)
24-month	105	90.5% (95/105)
36-month	83	80.7% (67/83)
48-month	67	70.1% (47/67)
60-month	23	56.5% (13/23)

<sup>1</sup> Eligible for follow-up = previous eligible for follow-up – (previous death + conversion + lost to follow-up + not due for next visit).

<sup>2</sup> Includes cases with at least one of any of the following submitted on or after the anniversary of the registration date: clinical form, death form, conversion form, imaging, or considered lost to follow-up.

In many cases where follow-up has not been provided, Cook has been notified by the site of potential loss of follow-up due to reasons such as the patient's unwillingness to return for follow-up, the site's inability to contact the patient (e.g., patient moved without providing forwarding address), or the site's inability or unwillingness to provide additional data (e.g., the site unable to identify the patient).

#### Summary of Post-Approval Study Results

##### Final Safety Findings

1. Post-Procedure (within 30 days of the procedure) Rate of Rupture and Subsequent Death:
  - No cases of post-procedure aneurysm rupture with subsequent death were reported within the registry.
2. Post-Procedure (within 30 days of the procedure) Proximal Attachment Status:
  - Two patients were reported to have a failed post-procedure proximal attachment of the Renu (i.e., barb separation, stent fracture/breakage, proximal Type I endoleak, Type III endoleak). Specifically, these two patients were identified with proximal Type I endoleak by the core laboratory and/or site within 30 days of the procedure.

##### 3. Post-Procedure (within 30 days of the procedure) Adverse Events:

- In addition to two intraoperative conversions to open surgical repair, 12 adverse events in 11 cases (9 events in 8 cases were considered device or procedure-related) have been reported within 30 days of the procedure. One of those adverse events (limb occlusion on side opposite of Renu implantation) led to a secondary intervention (additional stent implantation). No other secondary interventions were reported within 30 days of the procedure.
- Adverse event information was requested at each follow-up. **Table 13.4** lists the incidences of adverse events reported intraoperatively and within one month, categorized by organ system. All cases with any follow-up information were considered available for adverse event evaluation.

**Table 13.4 Incidence of Adverse Events within 30 days, by Category**

Category	Intraoperative	Adverse events Reported within 30 Days
Cases available for adverse event evaluation	151	143
Cardiovascular	0	0
Pulmonary	0	1
Renal	0	2
Bowel	0	1
Neurologic	0	0
Vascular	0	3
Wound	0	2
Other	0	3 <sup>1</sup>
<b>Total</b>	<b>0 events (0 cases)</b>	<b>12 events (11 cases)<sup>2</sup></b>

<sup>1</sup> Spontaneous retroperitoneal hematoma (1), allergic reaction (1), and fall with subsequent pubic ramus fracture (1).

<sup>2</sup> One case was reported with both a renal event and an "other" event (spontaneous retroperitoneal hematoma).

No adverse events were considered related to the Zenith Renu AAA Ancillary Graft, although 9 reported events for 8 cases were considered related to the procedure in which the Renu was placed.

##### Final Effectiveness Findings

Proximal Type I endoleaks were identified as one of the primary failure modes for pre-existing grafts, with 96 being identified preoperatively or during the procedure. Of the proximal Type I endoleaks reported, 99% (95/96) resolved without further intervention following Renu implantation; 1 persisted through 1-month follow-up and was converted to open surgical repair.

The majority of new endoleaks at follow-up have been identified as Type II. It is likely many had been previously overlooked, misidentified, or masked by other existing endoleaks during review at earlier time periods. Both proximal Type I endoleaks at 48 months were identified despite previous secondary interventions (graft placement and/or angioplasty) to treat the endoleak. The distal Type I endoleak present at 48 months was not present at 60 months according to the site. Both Type III endoleaks present at 48 months were successfully treated with additional graft implantation. These late endoleaks demonstrate the need for lifelong, regular follow-up to assess the patient's health and the performance of their endovascular graft. Cook also recommends that patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft component) receive enhanced follow-up. If Type I or III endoleak is identified, prompt intervention and additional follow-up post-intervention is recommended.

## Secondary Intervention

Secondary intervention information was requested at each follow-up. Each secondary intervention is listed in **Table 13.5**.

**Table 13.5 Secondary Interventions**

Time from Procedure (months)	Secondary Intervention Type	Reason for Secondary Intervention
1	Placement of additional stent	Limb occlusion on opposite side from Renu delivery
10	Placement of additional iliac leg graft component	Persistent distal Type I endoleak
12 <sup>1</sup>	Coil embolization	Persistent Type II endoleak
19	Axillary-femoral, femoral-femoral bypass	Limb occlusion
22	Placement of Zenith iliac leg graft to cover the endoleak	Persistent Type III endoleak in the limb of the original endograft
24 <sup>2</sup>	Angioplasty with implantation of a Palmaz stent	Proximal Type I endoleak
25 <sup>3</sup>	Coil embolization	Persistent proximal Type I endoleak
29 <sup>3</sup>	Placement of additional stent at neck	Persistent proximal Type I endoleak
30	Placement of Renu converter, TFLE leg extension, and occluder plug with a fem-fem bypass	Device migration <sup>4</sup> with persistent Type III endoleak
34 <sup>2</sup>	Percutaneous angioplasty	Proximal Type I endoleak
38	Zenith stent-graft placed within original graft	Persistent Type III endoleak
38	Proximal Renu cuff implanted	Persistent Type I endoleak
44	Implantation of branched iliac endograft	Iliac artery dilatation at attachment site
45	Two iliac limbs attempted (AneuRx® and Talent®) in difficult anatomy, which caused iliac artery to rupture. Patient underwent conversion to open surgical repair (see Conversion). Two thoracic endografts (both Talent) were implanted.	Distal Type I endoleak
48	Transfemoral catheter-based sac/lumbar embolization with Onyx bioglue	Persistent Type II endoleak
50	Gore bifurcated graft implanted to traverse across the dissociation between the Renu main body extension and the AneuRx® graft	Pre-existing graft separation from Renu and persistent Type III endoleak
52	Treatment with Renu converter	Graft separation and endoleak
60 <sup>1</sup>	Implantation of Renu converter with fem-fem bypass	Persistent proximal Type I and Type II endoleak

<sup>1</sup> Coil embolization at 12 months and implantation of a Renu converter with fem-fem bypass at 60 months were performed on the same patient.

<sup>2</sup> Angioplasty with implantation at 24 months and additional angioplasty at 34 months was performed on the same patient to treat proximal Type I endoleak.

<sup>3</sup> Coil embolization at 25 months and placement of an additional stent at the neck were performed on the same patient to treat proximal Type I endoleak.

When these interventions were not successful, conversion to open surgical repair was performed successfully at 30 months.

<sup>4</sup> Reported by site as migration of the pre-existing graft

All secondary interventions except one have been reported by the sites as successful. The secondary interventions for one persistent proximal Type I endoleak were not successful; therefore, the patient successfully underwent conversion to open surgical repair. Of note, these endoleaks were identified during routine follow-up, highlighting the need for continued follow-up to mitigate the risk of aneurysm expansion and rupture due to the presence of endoleak.

## Conversion

Nine conversions to open repair were reported. Each conversion is listed in **Table 13.6**

**Table 13.6 Conversions**

Months after Procedure	Reason for Conversion
0	Rupture of aortic wall proximal to aneurysm and Renu device
0	Leakage due to incomplete sealing/ persistent blood flow into the aneurysm from patent vessels (proximal Type I endoleak)
3	Leakage due to incomplete sealing of the aneurysm (proximal Type I endoleak) <sup>1</sup>
12	Leakage due to incomplete sealing of the aneurysm, inadequate sealing between the Renu main body extension and the AneuRx®, and aneurysm rupture <sup>2</sup>
12	Inadequate sealing between the Renu main body extension and the AneuRx® graft, and aneurysm rupture <sup>2</sup>
16	Leakage due to migration of pre-existing graft (AneuRx®) and aneurysm rupture
19	Infection of the pre-existing graft (AneuRx®) <sup>3</sup>
30	Leakage due to incomplete sealing of the aneurysm (proximal Type I endoleak)
45	Rupture of external iliac artery during secondary intervention to treat distal Type I endoleak (distal to Renu main body extension)

<sup>1</sup> Proximal Type I endoleak identified intraoperatively, but not treated.

<sup>2</sup> Failure to properly follow the instructions, warnings, and precautions in this Instructions for Use may lead to serious consequences or injury to the patient.

<sup>3</sup> Core lab analysis of pre-Renu™ imaging noted findings that were potentially indicative of infection. An independent CEC adjudicated this case to be unrelated to the Renu endovascular repair.

Two intraoperative conversions and seven late (> 30 day) conversions were reported. The intraoperative conversions to open repair were due to rupture of the aortic wall proximal to the Renu device (1) and proximal Type I endoleak (1). The patient with aortic wall rupture did not survive the conversion. The late conversions were related to a suspected graft infection (1), persistent proximal Type I endoleak (2), inadequate sealing between the Renu main body extension and the AneuRx® graft leading to component separation and aneurysm rupture (3), and rupture of the external iliac artery during secondary intervention to treat distal Type I endoleak (1).

**Mortality**

**Table 13.7** summarizes the 44 deaths that were reported, along with the results from CEC adjudication of each death.

**Table 13.7 Deaths**

Months after Procedure	Age at Registration	Cause of death	CEC Adjudication
0	82	Intraoperative rupture of aorta proximal to aneurysm and Renu device with subsequent conversion <sup>1</sup>	Procedure-related and technique-related
1	90 <sup>2</sup>	Low platelet count, hematological complications <sup>3</sup>	Procedure-related
1	78	Congestive heart failure 35 days post-procedure	Not related <sup>4</sup>
2	79	Wegener's granulomatosis <sup>5</sup>	Not related
3	80	Cardiopulmonary failure <sup>6</sup>	Cause unable to be determined
4	83	Direct cause of death not available to reporting institution <sup>7</sup>	Not related
7	90 <sup>2</sup>	Failure to thrive/old age	Not related
10	73	Cardiorespiratory arrest secondary to hypotension and sepsis	Procedure-related and technique-related
10	79	Myocardial infarction <sup>8</sup>	Cause unable to be determined
11	69	Recurrent cholangiocarcinoma	Not related
12	77	Unknown, information unable to be obtained by reporting institution <sup>9</sup>	Not related
12	80	Metastatic lung cancer	Not related
12	81	Multi-system organ failure following aortic aneurysm rupture and subsequent emergent conversion <sup>10</sup>	Renu-related
13	82	Cancer	Not related
14	69	Pneumonia with fever and septic shock	Not related
16	76	Cardiac arrest following aneurysm rupture and emergent conversion to open repair <sup>11</sup>	Renu-related
17	77	Pulmonary emboli secondary to malignancy	Not related
17	65	Cancer	Not related
18	80	Paralysis secondary to cervical fracture from fall. Renal failure.	Not related
19	75	Ventricular fibrillation, ischemic cardiomyopathy, and GI bleed	Not related
20	73	Congestive heart failure and respiratory failure secondary to congestive heart failure	Not related
20	79	Pulmonary	Not related
21	81	Complications from pneumonia and organ failure	Not related
21	76	Heart attack	Not related
23	85	Unknown, but believed by site to be unrelated to the Renu	Cause unable to be determined
24	Not provided	Unknown <sup>12</sup>	Cause unable to be determined
27	75	Unrelated to aneurysm, patient died at home <sup>13</sup>	Cause unable to be determined
28	67	Cardiac issues related to CHF	Not related
29	82	Cancer	Not related
30	69	Lung cancer	Not related
30	78	Ruptured AAA after patient refused to return for repair of Type II endoleak with aneurysm expansion	Pre-existing graft-related <sup>14</sup>
33	84	Cardiopulmonary arrest	Not related
36	67	Cancer	Not related
36	61	Unknown <sup>15</sup>	Cause unable to be determined
36	71	CVA, respiratory failure, cancer	Not related
36	74	Cardiac causes	Not related
36	77	Respiratory/cardiac arrest	Not related
36	90 <sup>2</sup>	Unknown <sup>16</sup>	Cause unable to be determined
37	82	Unknown – notification came from primary care physician	Cause unable to be determined
40	79	Cardiogenic shock following hip repair of femoral neck fracture	Not related
40	90 <sup>2</sup>	Respiratory failure	Not related
45	78	Cardiopulmonary arrest, respiratory failure, renal failure 2 days after secondary intervention <sup>17</sup>	Pre-existing graft-related <sup>18</sup>
48	72	Urosepsis, prostate cancer, MRSA	Not related
48	87	Cancer	Not related
Unknown	76	Unknown <sup>19</sup>	Cause unable to be determined

<sup>1</sup> The patient failed to recover from conversion to open surgical repair, which was performed to treat a rupture of the aorta proximal to the Renu device. The aorta was ruptured by a spicule of calcium after Renu deployment, either during deployment of a Palmaz stent or ballooning of the Renu device.

<sup>2</sup> To comply with HIPAA regulations, the age of any patient ≥ 90 years old was recorded and reported as 90 years.

<sup>3</sup> The patient was admitted with a low platelet count and an AneuRx® with a proximal Type I endoleak. The patient's aneurysm ruptured prior to the scheduled Renu implantation date, but was able to be treated emergently with the Renu converter.

<sup>4</sup> Previously reported as procedure-related; however, upon additional review of this event, the CEC adjudicated it to be unrelated to endovascular repair.

<sup>5</sup> Wegener's granulomatosis was reported to be a pre-existing comorbidity identified after the procedure in which the Zenith Renu AAA Ancillary Graft was implanted.

<sup>6</sup> Exact cause of death was unknown. The family described the death as related to cardiopulmonary failure; the patient had a documented 10-year history of severe cardiopulmonary disease.

<sup>7</sup> Patient was undergoing evaluation for neuromuscular degeneration. Per institution, there was no indication that death was related to the aneurysm or the endograft.

<sup>8</sup> The patient died after an MI; however, the CEC was unable to adjudicate the death because the patient passed away at home and an autopsy was not performed.

<sup>9</sup> Cause of death was unknown. Per the reporting institution, an autopsy was performed and the death was related to a pre-existing comorbidity. Based on the available information, the CEC determined the death to be unrelated to the endovascular repair.

<sup>10</sup> Rupture with emergent conversion was secondary to separation of the Renu main body extension from the pre-existing AneuRx® graft. A Renu converter had been recommended prior to the procedure.

<sup>11</sup> Rupture with emergent conversion was secondary to migration of the AneuRx® graft with subsequent Type III endoleak. A Renu converter had been recommended prior to the procedure.

<sup>12</sup> Cause of death was unknown. Per the reporting institution, the death was reported to the physician by the family. The CEC was unable to adjudicate the death based on the information provided.

<sup>13</sup> The site noted this death to be unrelated to the aneurysm; however, the CEC was unable to adjudicate the death without confirmation that the site obtained the information from a death certificate or an autopsy was performed.

<sup>14</sup> CEC adjudicated death as related to endovascular intervention, but related to the pre-existing graft, not the Renu.

<sup>15</sup> Cause of death was unknown. Per the reporting institution, the death was related to a pre-existing comorbidity and the graft was fine two months prior to death.

<sup>16</sup> Cause of death was unknown. Per the reporting institution, a family member had called the office with the date of death, but the cause was unknown. The site noted that the patient was being treated for cancer.

<sup>17</sup> Secondary intervention was to treat distal Type I endoleak. Iliac artery rupture during advancement of additional iliac legs led to conversion to open surgical repair. Patient died 2 days following intervention/conversion.

<sup>18</sup> CEC adjudicated death as related to endovascular intervention, but related to the salvage procedure, which was for a distal Type I endoleak of the pre-existing graft (distal to the Renu device).

<sup>19</sup> Cause of death was unknown. Per the reporting institution, information was not available about cause of death or time of death with respect to Renu implantation. The last contact the site had with the patient was at 1-month follow-up.

One intraoperative, one early ( $\leq$  30-day), and 43 late ( $>$  30-day) deaths were reported. One intraoperative, endovascular intervention-related death occurred after the patient failed to recover from conversion to open repair following rupture of the aorta proximal to the Renu device. One early death occurred approximately 2 weeks after the initial procedure; the patient died from hematological complications (including a low platelet count). Of note, the patient was treated emergently with the Renu after being admitted with a low platelet count and after aneurysm rupture. This event was determined to be procedure-related.

Death beyond 30 days of the initial procedure occurred in 43 cases. Twenty-nine cases were determined to be unrelated to endovascular repair. Of the 14 remaining cases, 5 were adjudicated as related to endovascular repair (i.e., procedure, technique, pre-existing graft and/or Renu-related) and 9 were unable to be adjudicated due to insufficient information from the site (no additional information was able to be collected).

Nine adverse events, 1 intraoperative and 2 late conversions with subsequent death, 4 additional conversions (1 intraoperative, 3 late), and 3 additional deaths related to the procedure, technique, and/or Zenith Renu AAA Ancillary Graft were reported. The 6 Renu-related conversions and deaths

were related to either proximal Type I (3 cases) or Type III (3 cases) endoleaks, 4 of which were first identified after 30 days. Of note, no unanticipated adverse events occurred (all procedure-related or device-related adverse events, deaths, and conversions are noted as possibilities in this Instructions for Use).

Of the 151 registered cases of Renu implantation, 89.4% (135/151) had no procedure or Renu-related adverse events, conversions, or deaths. These mid-term post-market registry data confirm that the Renu device may be used during secondary intervention to successfully treat proximal fixation failures.

#### **Study Strength and Weaknesses**

This registry was limited by the presence of multiple co-morbidities found in these cases, which confounded data analysis. The wide variety of failed device types, sizes, and specific circumstances that were encountered also confounded data analysis. Additional limitations of this registry include the lower incidence of AAA in women, with fewer women than men being enrolled in the registry. This registry was not designed to provide information beyond five years post-procedure. The registry had no control group and therefore comparison with alternative treatment was not an aim of this registry.

#### **14 PATIENT TRACKING INFORMATION**

In addition to these *Instructions for Use*, the Zenith Renu AAA Ancillary Graft with the Z-Trak Introduction System is packaged with a *Device Tracking Form* which the hospital staff is required to complete and forward to Cook for the purposes of tracking all patients who receive the Zenith Renu AAA Ancillary Graft (as required by U.S. Federal Regulation).

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© COOK 2022
Bloomington, IN 47404 U.S.A.