
- EN**  
6 **Renal Access and Biopsy Needle Set**  
Instructions for Use
- CS**  
7 **Souprava jehly pro ledvinový přístup a biopsii**  
Návod k použití
- DA**  
8 **Nålesæt til nyreadgang og -biopsi**  
Brugsanvisning
- DE**  
9 **Kanülenset für Nierenzugang und -biopsie**  
Gebrauchsanweisung
- EL**  
11 **Σετ νεφρικής πρόσβασης και βελόνας βιοψίας**  
Οδηγίες χρήσης
- ES**  
12 **Equipo de acceso renal y aguja de biopsia**  
Instrucciones de uso
- FR**  
14 **Le set d'accès et d'aiguille de biopsie rénal**  
Mode d'emploi
- IT**  
15 **Set con ago per accesso e biopsia del rene**  
Istruzioni per l'uso
- NL**  
17 **Niertoegangs- en biopsienaaldset**  
Gebruiksaanwijzing
- PT**  
18 **Conjunto de agulha para acesso e biopsia renal**  
Instruções de utilização
- SV**  
19 **Set för njuråtkomst med biopsinål**  
Bruksanvisning



**Fig. 2**


- 1 Etch Mark
- 1 Leptaná značka
- 1 Skraveringsmærke
- 1 Geätzte Markierung
- 1 Σήμανση χάραξης
- 1 Marca grabada
- 1 Repère gravé
- 1 Incisione di riferimento
- 1 Geëtste markering
- 1 Marca gravada
- 1 Markering

**Fig. 3**

**Fig. 4**

**Fig. 5**

## RENAL ACCESS AND BIOPSY NEEDLE SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Renal Access and Biopsy Needle Set includes a Quick-Core® Biopsy Needle (**Fig. 1**) and a preassembled renal access set consisting of a radiopaque sheath, a Check-Flo® valve adapter and a stiffening cannula. The set also contains one straight catheter and one curved catheter.

### INTENDED USE

The Renal Access and Biopsy Needle Set is intended for transjugular renal biopsy procedures.

The Quick-Core Biopsy Needle is intended for soft tissue biopsy.

### CONTRAINDICATIONS

None known

### WARNINGS

Extreme care must be exercised during manipulation and withdrawal of catheter to prevent pulling catheter apart. Insertion through a synthetic vascular graft should be avoided whenever possible.

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in transjugular renal biopsy techniques. Standard techniques for jugular venous access and catheter positioning should be employed.
- Fluoroscopic guidance is required to ensure accurate needle placement.
- When test-firing the Quick-Core Biopsy Needle, place distal tip of device against sterile solid surface before activating cutting cannula. This will prevent damage to the distal tip.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### POTENTIAL ADVERSE EVENTS

- Hemorrhage
- Puncture site hematoma
- Cardiac arrhythmias
- Arteriovenous fistulas
- Capsule perforation

### INSTRUCTIONS FOR USE

1. Utilizing standard access techniques, introduce an appropriate .035 inch wire guide into the inferior vena cava via jugular vein access. **CAUTION: Continuous cardiac monitoring is recommended while transiting through the right atrium.**
2. Using a selective catheter and wire guide of choice, catheterize the appropriate renal vein. Leave the wire in a safe, distal position, and remove the catheter.
3. Introduce the transjugular introducer sheath and stiffening cannula over the wire guide and advance into the selected renal vein. When introducing these components as a preassembled set, the included straight catheter may be used to facilitate introduction. Once the set is positioned within the renal vein, the catheter and wire guide should be removed. **NOTE:** Care should be taken to prevent damage to the straight catheter during removal through the metal stiffening cannula. Leaving a wire guide through the straight catheter during removal may aid in preventing catheter damage.
4. Prepare the Quick-Core Biopsy Needle by pulling back on the plunger until a firm click is felt, indicating the needle spring is locked into the ready position. (**Fig. 2**)
5. With the stylet fully retracted advance the Quick-Core Biopsy Needle through the renal access set until the etch mark on the proximal cannula is at the hub of the Check-Flo valve adapter. (**Fig. 3**) This will position the needle tip at the distal end of the access set. **CAUTION: For accurate needle placement under fluoroscopy, the X-ray tube must be obliquely positioned so that the X-ray beam is perpendicular to the needle and tangential to the capsule in the area of the renal biopsy. Use the directional arrow on the hub of the stiffening cannula to ensure that the X-ray beam is perpendicular to the needle.**
6. Using the directional arrow on the hub of the stiffening cannula, direct the renal access set in the desired direction for biopsy and apply gentle forward pressure to "tent" the wall of the renal vein adjacent to the biopsy site. (**Fig. 4**)
7. With the stylet fully retracted, so that the specimen notch is completely covered by the cutting cannula, advance the Quick-Core Biopsy Needle out the tip of the renal access set and into the kidney tissue. **Do not advance the stylet until the needle is in position.**

8. While maintaining needle position, advance the stylet to expose the specimen notch within the area to be biopsied. Fire the cutting cannula by fully depressing the plunger to capture tissue within the specimen notch. (**Fig. 5**)
9. Withdraw the Quick-Core Biopsy Needle from the biopsy area and completely out of the renal access set.
10. To remove the tissue specimen, pull back on the plunger until a firm click is felt. This indicates that the cutting cannula is locked into position. Push the stylet forward to expose the tissue specimen and remove the specimen from the notch.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd