**COOK**

MEDICAL

**CE**  
0123

- EN**  
2 **Renal Access Cobra Catheter**  
Instructions for Use
- CS**  
3 **Katetr kobra pro renální přístup**  
Návod k použití
- DA**  
3 **Cobra-kateter til nyreadgang**  
Brugsanvisning
- DE**  
4 **Cobra Katheter für den Nierenzugang**  
Gebrauchsanweisung
- EL**  
5 **Καθετήρας νεφρικής πρόσβασης σχήματος κόμπρας**  
Οδηγίες χρήσης
- ES**  
6 **Catéter Cobra de acceso renal**  
Instrucciones de uso
- FR**  
7 **Cathéter Cobra pour accès rénal**  
Mode d'emploi
- HU**  
8 **Veséhez való hozzáférést biztosító Cobra katéter**  
Használati utasítás
- IT**  
9 **Catetere per accesso renale Cobra**  
Istruzioni per l'uso
- NL**  
10 **Cobra-niertoeegangskatheter**  
Gebruiksaanwijzing
- NO**  
11 **Cobra-kateter for nyretilgang**  
Bruksanvisning
- PL**  
12 **Cewnik Cobra do dostępu nerkowego**  
Instrukcja użycia
- PT**  
13 **Cateter Cobra de acesso renal**  
Instruções de utilização
- SV**  
14 **Renal Cobra-kateter för åtkomst**  
Bruksanvisning

T \_ R A C C \_ R E V 2

## **RENAL ACCESS COBRA CATHETER**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.**

### **DEVICE DESCRIPTION:**

The Renal Access Cobra Catheter is composed of a polymer reinforced with stainless steel braided wire. The proximal end has a Luer lock hub, and the distal end is formed into a hook shape. The catheter is 65 cm long and is 6.0 French in diameter.

### **INTENDED USE**

This device is intended for delivery of contrast media as well as access, advancement, or exchange of wire guides in the renal area.

### **CONTRAINDICATIONS**

This device is contraindicated for use in patients unsuitable for catheterization.

### **WARNINGS**

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

### **PRECAUTIONS**

- This product is intended for use by physicians trained and experienced in renal access techniques. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- If you encounter resistance while advancing or withdrawing the catheter, stop. Determine the cause of the resistance before proceeding.
- Do not withdraw the catheter while it is deflected in the cystoscope/endoscope.

### **POTENTIAL ADVERSE EVENTS**

Adverse events that may result from catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

### **INSTRUCTIONS FOR USE**

#### **Percutaneous Placement**

1. Flush the catheter.
2. With a wire guide already in place within the renal pelvis, advance the catheter over the wire guide and into the renal pelvis.
3. Retract the wire guide to allow the tip of the catheter to form.
4. Advance the catheter to the desired location within the renal pelvis. Slightly advance and retract the wire guide to help navigate the catheter.
5. Once the catheter is in the desired location, deliver the working wire guide of choice through the catheter.
6. Maintain access with the wire guide while removing and discarding the catheter.

### **HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd