COOK

MEDICAL

CE
0123

EN
3

Ring-McLean Sump Drainage Set

Instructions for Use

CS
4

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DA
5

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Brugsanvisning

DE
6

Ring-McLean Abszessdrainage-Set

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EL
7

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T _ R M S _ R E V S

RING-MCLEAN SUMP DRAINAGE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Ring-McLean Sump Drainage Sets consist of:

- Double-lumen Ultrathane® sump catheter
- Stiffening cannula with stainless steel trocar stylet
- TFE-coated curved Newton LLT wire guide, .038 inch, 100 cm long, 15 mm J
- Connecting tube, 14 Fr, 30 cm long
- StatLock® fixation device*

INTENDED USE

The Ring-McLean Sump Drainage Set is intended for abscess drainage.

CONTRAINDICATIONS

None known

WARNINGS

- **Any Ring-McLean Sump Drainage Set with the suffix “ACL” contains natural rubber latex, which may cause allergic reactions.**

PRECAUTIONS

- During catheter placement or positioning, leave the blunt cannula in place while aspirating to occlude the sump.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous catheters should be employed.

INSTRUCTIONS FOR USE

1. Localize the collection to be drained by palpation, fluoroscopy, ultrasound or CT. Anesthetize the skin and subcutaneous tissues overlying the collection and make a skin incision about 1 cm in length. Separate the tissue with a hemostat.

For Direct Trocar Access

1. For large, superficially located fluid collections or abscesses, the sump is introduced directly.
2. With the trocar stylet in place, advance the sump directly into the collection. Remove the sharp stylet and aspirate the inner cannula to confirm that the tip is in the collection. Leave the blunt cannula in place during initial aspiration to occlude communication with the sump lumen of the catheter.
3. Advance the wire guide through the cannula for several cm into the collection; it will eventually buckle and loop around in the cavity. The blunt cannula is detached from the sump catheter fitting and held with one hand as the sump is advanced over the wire guide with the other hand. Once the sump is in position, remove the blunt cannula.

For Seldinger-Type Access

1. For deep or small collections, it may be advisable to access initially with a needle or access set (e.g., NPAS-104-RH-NT) and then exchange for the sump catheter over a wire guide once the collection has been entered.
2. Access the target collection until fluid return is obtained. Inject contrast medium to confirm an intracavity location.
3. Once position is confirmed, advance .038 inch wire guide through the access device into the collection space.
4. Remove the access device (needle or access set).
5. Remove the trocar stylet of the sump and introduce the sump catheter with the inner cannula in place over the wire guide. The wire is held taut until the sump has passed through the skin and subcutaneous tissue.
6. Once the tip of the sump is in the collection, the cannula is detached from the sump fitting and held as the sump is advanced over the wire guide into the collection. Secure the sump catheter to the patient.

NOTE: Keep the second lumen open to allow air to pass to the distal tip of the catheter, resulting in free flow of the fluid and reduction of negative pressure.

A syringe may be attached to the side-arm fitting to flush the second lumen.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use

the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

* StatLock® is a registered trademark of Venetec International, Inc., San Diego, California.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd