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Bruksanvisning

ROADRUNNER® PC WIRE GUIDES WITH AQ® HYDROPHILIC COATING

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Roadrunner PC Wire Guide is coated with AQ (a biocompatible hydrophilic coating), which when activated becomes lubricious. The wire guide is made of a nitinol core with a distal platinum inner spring coil tip, which is completely covered by a bismuth-impregnated polymer jacket.

	Taper Length (cm)
Firm	10.5
Firm LT	18.5
Nimble	11.5
Nimble Floppy	9

INTENDED USE

The Roadrunner PC Wire Guide is used for catheter positioning and exchange in diagnostic and interventional procedures, exclusive of the coronary arteries.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the coronary arteries.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Use medical imaging when manipulating the wire guide. Do not advance or torque the wire guide without visual evidence of the corresponding movement of the distal tip.
- When using the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- Always keep the surface of the wire guide wet during use for optimal results.

INSTRUCTIONS FOR USE

Activating the AQ® Hydrophilic Coating


Fig. 1

1. Flush the wire guide holder by attaching or pressing a syringe with saline, heparinized saline solution, or sterile water to the fitting or opening of the wire guide holder. Inject enough solution to wet the wire surface entirely. This will activate the AQ® coating.
2. Carefully remove the wire guide from the holder.
3. Insert a wire guide insertion tool through the valve assembly of the guiding sheath or guiding catheter. Insert the distal tip of the wire guide through the insertion tool.

4. Attach a torque device to the wire guide (if provided).
5. Standard wire guide endovascular techniques may now be employed.
NOTE: An initial flush of heparinized saline solution through the catheter lumen will enhance initial movement of the wire guide through the catheter.
6. When the wire guide is removed from a catheter, wipe it with a sterile gauze pad moistened with saline, heparinized saline solution, or sterile water and place it back into the original holder. Fill the holder with heparinized saline solution to maintain the lubricous state of the surface of the product for subsequent use.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd