en  
5

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


## ROADRUNNER® PC WIRE GUIDES

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

### DEVICE DESCRIPTION

Roadrunner® PC Hydrophilic Wire Guides and Roadrunner® PC Double Flexible Tip Wire Guides are designed with a flexible or stiff shaft and are manufactured in 0.035 and 0.038 inch diameters and a length of 145 cm. These wire guides are manufactured as fixed core wire guides with polymer coated centerless ground nitinol core secured to a platinum spring coil at the distal tip. These wire guides are designed with either an angled or straight distal tip. The devices identified by a part number containing “-I” have incremental ink marks and devices identified by a part number containing “-DF” have double flexible tips. The spacing of the ink marks is intended as reference for visual aid only. These devices are manufactured with hydrophilic coating on the surface of the polymer jacket. These devices are packaged with an Olcott Torque Device, which is a polymer ancillary device that can assist in wire guide manipulation, and a wire guide holder.

### Performance Characteristics

- The devices identified by a part number containing “-DF” are double flexible wire guides, which allow for increased flexibility.
- Hydrophilic coating to enhance lubricity of device.
- The Olcott Torque device provides additional control when rotating the wire guide.
- The devices identified by a part number containing “-I” have ink marks for visual reference.
- Platinum coiling material for resistance to kinking.

### Device Compatibility

These wire guides are compatible with sterile water and saline, as well as appropriately sized needles, sheaths, dilators, endoscopes stents and catheters. The wire guide holder is compatible with a syringe. The Olcott Torque Device accepts 0.015 to 0.045 inch diameter wire guides.

### Patient Population

These devices are used for adult patients with a urological condition who require a procedure using a wire guide.

### Intended Users

These devices are intended for use by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

### Contact with Body Tissue

These devices interact with the structures and tissues of the upper and lower urinary tract.

### Operating Principle

These devices are inserted into the body through the lumen of a percutaneous entry needle, guiding sheath, or other interventional device. Only the flexible tip is intended as the initially inserted end. After insertion, medical imaging is used to visualize the movement of the wire guide to the desired location within the body. Catheters or other devices can then be inserted over the wire guide and directed to the intended location.

### INTENDED USE

These devices are platinum tipped, polymer and hydrophilically coated, nitinol fixed core wire guides that assist the insertion, removal, and positioning of other devices within the urinary tract. These devices act as a guide to access a portion of the urinary tract during a urological procedure to allow placement, replacement, or exchange of medical devices.

## INDICATIONS FOR USE

Multiple medical conditions may result in the need for a urological procedure. Urology wire guides are used as a guide to access a portion of the urinary tract to allow placement, replacement, or exchange of medical devices during urological procedures such as percutaneous nephrostomy (PCN), percutaneous nephrolithotomy (PNL), ureteroscopy (URS), and stent placement or exchange.

## DEVICE AND CLINICAL BENEFITS

To assist with urinary tract access for minimally invasive placement, replacement, or exchange of medical devices (without the wire guide itself having a direct diagnostic or therapeutic function).

## CONTRAINDICATIONS

None

## WARNINGS

- Do not pull the wire guide through a needle with a sharp edge when the wire is not aligned straight with the end of the needle. Extreme caution should be observed any time the wire guide is used within a needle. Manipulation or removal of the wire guide through a needle may result in destruction/separation of the polymer jacket of the wire requiring retrieval.
- Do not pass the wire through a damaged scope or ancillary device. A compromised lumen may cause damage or separation to the polymer jacket of the wire guide.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Do not force the device during the procedure; proceed with caution if excessive resistance is encountered. Excessive force can lead to the wire guide fracturing, bending, kinking, unraveling or uncoiling. Additionally, excessive force can cause the wire guide to become entrapped.
- If any resistance is felt or if the wire behavior and/or location seem improper, stop manipulating the wire and determine the cause. Failure to exercise proper caution may result in bending, kinking, separation of the wire guide jacket, damage to ancillary devices or damage to the urinary system. If necessary, remove the wire guide and ancillary device or scope as a complete unit to avoid possible damage and/or complications.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

## PRECAUTIONS

- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement or migration.
- Manipulation of wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access. Over manipulation can cause the wire guide to become entrapped.
- Avoid retracting the wire guide through a metal cannula or scope when the devices are not aligned straight, as the metal device can damage the polymer jacket of the wire guide.
- The surface of the wire guide is not lubricious unless it is wet. Before taking it out of its holder, activate the wire guide hydrophilic coating with sterile water or sterile saline. Removing the wire guide before activating the coating may result in damage to the polymer jacket of the wire guide.
- End hole size and length of device must be taken into consideration to ensure a proper fit between the wire guide and the procedural device.

## POTENTIAL ADVERSE EVENTS

- Acute Bleeding / hematoma
- Acute pyelonephritis
- Hematuria
- Infection / UTI / urosepsis
- Injury to adjacent organs
- Mucosal tear / laceration / abrasion

- Pain
- Ureteral laceration / avulsion / tear / lesion
- Ureteral perforation / false passage / urinoma
- Ureteral stenosis
- Ureteral stricture
- Urinary tract injury
- Vascular injury

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

## INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

## INSTRUCTIONS FOR USE

### Instructions for activating hydrophilic coating:

The hydrophilic coating on the wire guide is activated by immersion in sterile water or sterile saline solution.

1. Prior to using the wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder.
2. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. (Fig. 1)

#### NOTE:

- For Roadrunner PC Double Flexible Tip Wire Guide, the pink end of the wire guide is not hydrophilically coated or intended for insertion into body.
- Hydrophilic coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- For optimal performance, rehydrate the hydrophilic coated wire guide after exposure to ambient environment or after extended use; replace it with a new hydrophilic coated wire guide.

### Instructions for using olcott torque device:

1. With thumb-slide in a forward position, place Olcott Torque Device over proximal end of wire guide to desired location.
2. Pull thumb-slide (Fig. 2) back until wire guide is held securely in place. (Fig. 3)

### Instructions for using a wire guide:

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. Remove the wire guide from the holder.
3. Advance the flexible end of the wire guide into the device and position to the desired anatomical location.
4. After wire guide use is complete, remove from patient and discard.

## DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

## PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

## SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
**Cook Medical Europe Ltd**
**William A. Cook Australia Pty Ltd**